

# Soil Extraction Cell: An Alternative Method of Soil Extraction for Organics

Joe Boyd

Environmental Express

Charleston, SC

# Various Extraction Techniques



- Sonication by SW846 Method 3550C
- Soxhlet by SW846 Method 3540C
- Accelerated Solvent Extraction by SW846 Method 3545A
- Microwave Extraction by SW846 Method 3546

# Advantages to Method 3546



- Significant decrease in solvent consumption
- Multiple extractions at one time
- Short extraction time
- Closed environment

# Method 3546 Specifics



- Written for the microwave
- Uses elevated temperature and pressure to achieve analyte recoveries
- Applicable to semi-volatile compounds including organophosphorous pesticides, organochlorine pesticides, chlorinated herbicides, phenoxyacid herbicides, substituted phenols, PCBs, and PCDDs/PCDFs

# Recommended Extraction Conditions



- Temperature: 100 – 115°C
- Pressure: 50 – 150 psi
- Time at Temperature: 10 – 20 minutes
- Cooling: To room temperature
- Filtering / Rinsing: With same solvent system

# Equipment Specifications



- Section 11.9.1 of Method 3546 states “In general, the pressure is not a critical parameter, since it is a result of the solvent system vapor pressure at the elevated temperature.”
- Section 6.1.1 requires that the system “be capable of sensing the temperature to within +/- 2.5°C and automatically adjusting the .... output power within 2 sec of sensing. Temperature sensors should be accurate to +/- 2°C.”
- The HotBlock™ meets the requirements of Section 6.1.1

# Performance Based Measurement System (PBMS) Approach



Goal of PBMS is to break down barriers

- Conveys “what” needs to be accomplished, but not prescriptively “how” to do it
- Allows for flexibility in method selection
- Expedites approval of new and emerging technologies to meet mandated monitoring requirements
- Develops and uses cost-effective methods that meet program requirements and their associated performance criteria

# Method 3546 Statements



- Section 1.1 states “Other systems and other types of vessels may be used, provided that the analyst demonstrates appropriate performance for the specific application.”
- Section 1.5 states “This method has been validated using a solvent mixture of hexane and acetone (1:1) .... This solvent system or other solvent systems may be employed, provided that adequate performance is demonstrated for the analytes of interest.”



# Soil-XCell System



# Soil Extraction Cell Procedure




- Weigh 10 – 30g of sample into the Soil Extraction Cell.
- Add approximately 30 mls of solvent.
- Add any surrogates or spikes required by the determinative method.
- Hand-tighten the outer cap onto cell.
- Heat HotBlock to 130°C. This will yield an internal solvent temperature of 100 - 115°C.
- Heat cells in HotBlock for 10 – 20 minutes once the temperature has stabilized.
- Remove cells and allow to cool to room temperature.
- Filter and rinse with same solvent system.



# Soil-XCell and Method 3546

- The Soil-XCell System creates the same condition inside the soil extraction cells as the microwave creates in it's vessels
- Soil-XCell methods validated using Methylene Chloride as solvent system



# Temperature and Pressure Studies



- Temperature and pressure studies were performed using 30 g of soil and 30 mls of Methylene Chloride
- Solvent temp = 110°C
- Temp at top of cell = 79°C
- Pressure calculated to 91 psig using Antoine's Coefficient



# Results of a CRM for PAHs



- Wibby Lot # 7037-13
- 14 replicates analyzed using the Soil-XCell System for prep and Method 8270 for the analysis
- All 47 compounds present in the CRM were detected within the certified acceptance range
- 9 additional compounds that were not in the CRM were analyzed for and resulted in non-detects
- Data available at [www.envexp.com](http://www.envexp.com)

# Results of a CRM for TPH



- Wibby Lot # 6038-10
- Certified value of 914 mg/kg
- 15 replicates analyzed using the Soil-XCell System
- Average result of 698 mg/kg with all 15 replicates falling within the acceptance range of the standard

# TPH CRM data



Replicate	Result	Certified Value	% Recovery	Acceptance Range
Blank	ND	N/A	N/A	N/A
Rep 1	690 mg/kg	914 mg/kg	75.5%	236 – 1650 mg/kg
Rep 2	698 mg/kg	914 mg/kg	76.4%	236 – 1650 mg/kg
Rep 3	680 mg/kg	914 mg/kg	74.4%	236 – 1650 mg/kg
Rep 4	663 mg/kg	914 mg/kg	72.5%	236 – 1650 mg/kg
Rep 5	695 mg/kg	914 mg/kg	76.0%	236 – 1650 mg/kg
Rep 6	644 mg/kg	914 mg/kg	70.5%	236 – 1650 mg/kg
Rep 7	605 mg/kg	914 mg/kg	66.2%	236 – 1650 mg/kg
Rep 8	662mg/kg	914 mg/kg	72.4%	236 – 1650 mg/kg

# TPH CRM data



Replicate	Result	Certified Value	% Recovery	Acceptance Range
Rep 9	747 mg/kg	914 mg/kg	81.7%	236 – 1650 mg/kg
Rep 10	662mg/kg	914 mg/kg	72.4%	236 – 1650 mg/kg
Rep 11	658 mg/kg	914 mg/kg	72.0%	236 – 1650 mg/kg
Rep 12	718 mg/kg	914 mg/kg	78.6%	236 – 1650 mg/kg
Rep 13	750 mg/kg	914 mg/kg	82.1%	236 – 1650 mg/kg
Rep 14	804 mg/kg	914 mg/kg	88.0%	236 – 1650 mg/kg
Rep 15	796 mg/kg	914 mg/kg	87.1%	236 – 1650 mg/kg
Average	698 mg/kg	914 mg/kg	76.4%	236 – 1650 mg/kg



# Results of Real-World Samples for PAHs Compared to Method 3545A



- “Real-World” samples analyzed side-by-side using the Dionex™ ASE™ System and Soil-XCell System
- Samples extracted using Soil-XCell method with an internal solvent temperature of 110°C
- Methylene Chloride used as the solvent system

# Results of Real-World Samples for PAHs Compared to Method 3545A



Compound	ASE Replicate	Soil-XCell Rep 1	Soil-XCell Rep 2	Soil-XCell Rep 3
2-Fluorobiphenyl	66 % Recovery	70 % Recovery	66 % Recovery	71 % Recovery
Nitrobenzene-d5	66 % Recovery	63 % Recovery	68 % Recovery	64 % Recovery
Terphenyl -d14	39 % Recovery	66 % Recovery	66 % Recovery	66 % Recovery
2-Methylnaphthalene	< 260 mg/kg	< 260 mg/kg	900 mg/kg	710 mg/kg
Acenaphthene	840 mg/kg	951 mg/kg	3000 mg/kg	3100 mg/kg
Acenaphthylene	< 260 mg/kg	< 260 mg/kg	330 mg/kg	470 mg/kg
Anthracene	2600 mg/kg	2700 mg/kg	6900 mg/kg	9400 mg/kg
Benz(a)anthracene	5200 mg/kg	5600 mg/kg	19000 mg/kg	20000 mg/kg
Benzo(a)pyrene	4600 mg/kg	4800 mg/kg	16270 mg/kg	16000 mg/kg
Benzo(b)fluoranthene	5700 mg/kg	6800 mg/kg	22000 mg/kg	21000 mg/kg
Benzo(ghi)perylene	3200 mg/kg	3400 mg/kg	9100 mg/kg	8800 mg/kg
Benzo(k)fluoranthene	2100 mg/kg	2100 mg/kg	5800 mg/kg	5400 mg/kg

# Results of Real-World Samples for PAHs Compared to Method 3545A



Compound	ASE Replicate	Soil-XCell Rep 1	Soil-XCell Rep 2	Soil-XCell Rep 3
Chrysene	5100 mg/kg	5200 mg/kg	17000 mg/kg	17000 mg/kg
Dibenz(a,h)anthracene	860 mg/kg	730 mg/kg	2800 mg/kg	2400 mg/kg
Fluoranthene	18000 mg/kg	14000 mg/kg	28000 mg/kg	27000 mg/kg
Fluorene	1200 mg/kg	1200 mg/kg	3300 mg/kg	3800 mg/kg
Indeno(1,2,3-cd)pyrene	2800 mg/kg	2700 mg/kg	8100 mg/kg	8000 mg/kg
Naphthalene	430 mg/kg	570 mg/kg	2800 mg/kg	2000 mg/kg
Phenanthrene	10000 mg/kg	9200 mg/kg	24000 mg/kg	24000 mg/kg
Pyrene	14000 mg/kg	11000 mg/kg	23000 mg/kg	24000 mg/kg

- More data available at [www.envexp.com](http://www.envexp.com)

# Results of Real-World Samples for PCBs Compared to Method 3545A



- “Real-World” samples analyzed side-by-side using the Dionex ASE System and Soil-XCell System
- Samples extracted using Soil-XCell method with an internal solvent temperature of 110°C
- Methylene Chloride used as the solvent system

# Results of Real-World Samples for PCB's Compared to Method 3545A



Compound	ASE Replicate	Soil-XCell Rep 1	Soil-XCell Rep 2	Soil-XCell Rep 3
<b>DCBP</b>	60% Recovery	104% Recovery	132% Recovery	68% Recovery
<b>TCMX</b>	55% Recovery	56% Recovery	55% Recovery	58% Recovery
<b>PCB - 1016</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1221</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1232</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1242</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1248</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1254</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1260</b>	2500 mg/kg	1900 mg/kg	2000 mg/kg	1300 mg/kg
<b>PCB - 1262</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg
<b>PCB - 1268</b>	<440 mg/kg	<440 mg/kg	<440 mg/kg	<440 mg/kg

- More data available at [www.envexp.com](http://www.envexp.com)

# Soil-XCell vs. Microwave



- Monetary investment is much less for the Soil-XCell
- Stainless Steel cells last longer than Teflon vessels
- Solvent systems used in the microwave systems must have a percentage (>10%) of a polar component as pure hydrocarbon solvents do not absorb microwave energy
- Soil-XCell System utilizes an external heat source so any solvent system may be employed creating more flexibility

# Conclusions



- Soil-XCell System follows SW846 Method 3546
- Utilizes the PBMS approach to method flexibility
- Creates exact same conditions in soil cells for extraction as the microwave does in it's vessels
- Is consistent with Method 3446 which allows for other types of heating devices and vessels to be used
- Shows acceptable recoveries for CRM materials
- Achieves similar analytical results for real-world samples when compared to Method 3545A



# Questions?



**ENVIRONMENTAL EXPRESS, INC.**

1-800-343-5319

[info@envexp.com](mailto:info@envexp.com)

[www.envexp.com](http://www.envexp.com)