

# Emerging Interest in GC-Triple Quadrupole Mass Spectrometry for Environmental Analysis

*National Environmental Monitoring Conference  
August 18, 2011  
Jason Cole, Product Manager GC/MS*

# Presentation Overview

---


- What are the analytical strengths of GC/QqQ?
- Example applications demonstrating strengths
  - Organotins in water
  - PAH in rubber
  - Dioxins in food and feed



## **Analytical Strengths of GC Triple Quadrupole**

# Triple Quadrupole MS – Principle of Operation

## Quantitation of target compounds in matrix samples



# Analytical Strengths of Triple Quadrupole in SRM

Analytical Strength	Analytical Reason	Application Benefit
High sensitivity	SRM has a high duty cycle	Lowered LOQ's/LOD's
High selectivity	MS/MS eliminates most isobaric interferants	Lowered LOQ's/LOD's Less sample prep High speed of analysis
High specificity	MS/MS eliminates most isobaric interferants	Less false positives

*SRM process offers opportunity for lowered LOQs, reduced sample prep, higher speed of analysis, and less risk of false positives*



● **Organotins in Drinking Water**



● **The world leader in serving science**

# Application Challenge and Approach

---

- Lowered limits set by EU Water Directive
  - Maximum allowable concentration – 1.5 ng/L (ppt)
  - Maximum annual average – 0.2 ng/L (ppt)
- Approach
  - Extraction
    - 400 ml of water at pH5
    - Extraction with pentane
    - Evaporated to 400  $\mu$ l
    - Derivatized with sodium tetra-ethyl borate
  - 3  $\mu$ L injected into TSQ Quantum XLS, GC triple quadrupole

# Customer Limits of Detection Using SIM


---

Compound	LOD Using SIM
Monobutyltin	10 ng/L
Dibutyltin	5 ng/L
Tributyltin	5 ng/L
Triphenyltin	1 ng/L

*Greater than order of magnitude improvement in detection limits needed*



# Organotin at 0.5 ng/L



	Compound	Retention time	Transition
1	Monobutyltin	7.93	235.08>178.95
2	Dibutyltin	9.42	261.03>205.03
3	Tributyltin	10.67	289.09>176.95
4	Monooctyltin	11.35	289.09>176.95
5	Tetrabutyltin	11.74	289.09>176.95
6	Diocyltin	14.94	261.03>148.98
7	Triphenyltin	16.23	351.15>196.94
8	Tricyclohexyltin	16.29	233.08>150.98

**Injection of 0.2 pg absolute on column**

**Corresponding with 0.05ng/L in water sample**

**New EU directive: 0.2ng/L for tributyltin and here 0.05ng/l is shown**




● **PAHs in Rubber**



● **The world leader in serving science**


# The Case of the Problem Sample

- Sample matrix: rubber
- Looking for PAH's that come from extender oils
- Full sample prep applied to sample




# Comparison of Fullscan and SRM


Pyrene in rubber  
by fullscan



Pyrene in rubber  
by SRM



Pyrene in Standard  
by SRM



## Compounds Detected – SRM vs Fullscan


---

Compound	SRM*	Fullscan
Naphthalene	X	X
Acenaphthylene	X	
Acenaphthene	X	
Anthracene	X	
Phenanthrene	X	
Pyrene	X	

\* All compounds confirmed with retention times and ion ratios (2)

# Benefits of Triple Quad for the Analysis

- Less sample prep required
- Ability to measure “worst case” samples
- Greater confidence in confirmation





● **Dioxins in Food and Feed**



● **The world leader in serving science**

# Why GC-Triple Quad for Dioxins Analysis?


- Screening

- Higher quality of screening results
  - TEQs can be calculated from results of screen [1]
  - Congener specific - Results depend on PCDD/F and PCB congener profile
- Improvement of screening sensitivity for false negative results below 1%
  - Could result in a improved specificity, thus affecting the overall cost. [2]

- Confirmation

- Lower upfront cost
- More versatile approach

- References:
  - [1] Kotz, Malisch 2010, Dioxon Conference
  - [2] Fochi, Brambilla 2008





# Maximum Allowed Levels

---

**Commission Regulation (EC) No. 1881/2006**  
**Commission Directive 2006/13/EC)**


	<b>Maximum levels</b> WHO-PCDD/F-TEQ
<b>Food</b> (Meat, liver, eggs, milk products, fats and oils, fish)	<b>0.75 – 6</b> pg/g fat <b>4</b> pg/g wet weight (fish)
<b>Feed</b> (Compound feed, premixtures, additives ...)	<b>0.75 – 6</b> ng/kg product (12 % moisture content)

# TSQ Quantum XLS - TCDD in Buffalo Milk Samples

Blank (GC-MS/MS)

Buffalo Milk (GC-MS/MS)

Buffalo Milk (GC-HRMS)




# Buffalo Milk Sample, HxCDD

GC-MS/MS (5 µl injection PTV solv. split) → GC-HRMS 1/10 concentration


## GC-MS/MS

RT: 30.33 - 32.83 SM: 5G




## GC-HRMS

RT: 25.25 - 27.74 SM: 5G




# TSQ Quantum XLS vs. GC-HRMS of Certified Material



# TSQ Quantum XLS Ultra – Food Dioxin Screening

- Courtesy Alexander Kotz, EURL Freiburg, Germany

- Presented at the 6th POPs Users Meeting, Niagara-on-the-Lake, Canada, Apr 28/29, 2011



# Conclusion From Webinar by Dr. Kotz

---

- GC-MS/MS systems in principle applicable for PCDD/F and DL-PCB analysis in food and feed samples
- Sensitive enough to check for dioxins at maximum allowed levels
- GC/HRMS still preferred at lower concentration

# Presenation Summary

---

- Laboratories are showing increasing interest in GC triple quad technology because of inherent advantages in the technology:
  - Low LOQ's from combination of high selectivity with high sensitivity
  - High specificity for more confident results
- Examples of application of technology
  - Organotins in drinking water for reduced LOQ's
  - PAH in rubber for confident detection in worst case matrix scenario
  - Dioxins in food and feed for higher sensitivity/specificity screening

---

**Thank you very much for  
your attention!**