



## **Sensitive Determination of Hexavalent Chromium in Drinking Water**

**Brian De Borba, Lipika Basumallick, Jeffrey Rohrer**

# Outline

---

- Why do we need a sensitive method for hexavalent chromium analysis?
- U.S. EPA Method 218.6 and optimizations proposed in 2003
- Modifications proposed in 2011 to permit a detection limit of 1 ppt

# Health Effects of Chromium

---

## **Chromium-3**

**A nutritionally essential element  
often added to dietary vitamin supplements**

## **Chromium-6 (Chromate)**

**Strong oxidizers  
Genotoxic carcinogens**


# Hexavalent Chromium in the Media

---

- Original case filed in 1993 (Hinkley, CA, 0.58 ppm)
- Chemical compound in Erin Brockovich case (2000)
- Environmental Working Group Report, December 2010




# Widespread Detection of Chromate in U.S. Tap Water



Environmental Working Group Report, December 2010

**31/35 cities had chromate in their water.**

# Chromium Regulations at a Glance



1977	Maximum Contaminant Level (MCL) for total chromium was established at 50 µg/L in California, adopted from what was then a National Interim Drinking Water Standard for chromium.
1991	Federal MCL for total chromium was raised to 100 µg/L, but California stayed at 50 µg/L.
1999	California EPA Office of Environmental Health Hazard Assessment (OEHHA) established a PHG at 2.5-µg/L for total chromium. California Department of Public Health identified chromium as a contaminant for possible MCL revision, and included hexavalent chromium among the unregulated chemicals requiring monitoring.
2001	National Toxicity Program announced it would conduct long-term studies to evaluate the potential carcinogenicity of ingested hexavalent chromium.
2008	U.S. EPA conducted a review of the health effects of hexavalent chromium based on toxicity studies performed by National Toxicology Program.
2009	California OEHHA proposed a PHG of 0.06 µg/L for hexavalent chromium.
2010	U.S. EPA released Toxicological Review of hexavalent chromium. California OEHHA issued new PHG for hexavalent chromium at 0.02 µg/L.
2011	U.S. EPA will carefully review the conclusions and consider all relevant information to determine if a new standard needs to be set. In the interim period, U.S. EPA provided guidelines for monitoring hexavalent chromium (and continuing to monitor total chromium).

# Available Methods for Hexavalent Chromium Detection

---

## U.S. EPA Method 218.4

- Chelation extraction, atomic absorption
- Positive interference from some metals
- Cumbersome, not automated
- Modest detection limits (~5 ppb)

## U.S. EPA Method 218.6 (ASTM Method D5257-03)

- Ion chromatography separation of chromate coupled with postcolumn reaction (diphenylcarbazide) and UV-Vis detection (530 nm)

# Summary of U.S. EPA Method 218.6

---

An aqueous sample is filtered through a 0.45  $\mu\text{m}$  filter and the filtrate is adjusted to a pH of 9 to 9.5 with a buffer solution.  
Buffer Solution: 2500 mM Ammonium Sulfate and 1000 mM Ammonium Hydroxide

A measured volume of the sample (50-250  $\mu\text{L}$ ) is introduced into the ion chromatography system.

A guard column removes organics from the sample before the Cr(VI) as  $\text{CrO}_4^{2-}$  is separated on an anion-exchange-separator column.

Postcolumn derivatization of the Cr(VI) with diphenylcarbazide is followed by detection of the colored complex a 530 nm.

**\*\*Samples must be stored at 4 °C and analyzed within 24 hours of collection.**




# Chromatography Conditions for U.S. EPA 218.6

---

<b>Columns:</b>	
Guard Column	Thermo Scientific Dionex IonPac™ NG1
Separator Column	Thermo Scientific Dionex IonPac AS7
Eluent	250 mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> 100 mM NH <sub>4</sub> OH
Flow Rate	1.5 mL/min
Postcolumn Reagent	2 mM Diphenylcarbohydrazide 10% v/v CH <sub>3</sub> OH 1 N H <sub>2</sub> SO <sub>4</sub>
Flow Rate	0.5 mL/min
Detector	Visible 530 nm
Retention Time	3.8 min

# System Configuration for Hexavalent Chromium by U.S. EPA Method 218.6



# Thermo Scientific Dionex IonPac AS7 Anion-Exchange Column

---

## The Dionex IonPac™ AS7:

- Separates a wide variety of polyvalent anions, including polyphosphates, polyphosphonates, and other polyvalent complexing agents such as EDTA and NTA using acidic eluent (eliminating metal interferences) with postcolumn derivatization and UV-Vis detection.
- Has a unique polymer packing that provides superior performance for separating ionic and polar compounds and offers high speed, efficiency, and loading capacity at moderate backpressures.
- Is useful for Cr(VI) in environmental matrices, such as ground water, wastewater, and soil extracts.

# Method Detection Limit (MDL) for Cr(VI) with U.S. EPA Method 218.6

---

Matrix Type	Concn. Used to Compute MDL* (µg/L)	MDL* (µg/L)
Reagent Water	1	0.4
Drinking Water	2	0.3
Ground Water	2	0.3
Primary Sewage Wastewater	2	0.3
Electroplating Wastewater	2	0.8


\* MDL = (Standard Deviation) × ( $t_{s, 99}$ ), where  $t_{s, 99} = 3.14$  for  $n = 7$ .

# Modified Version of U.S. EPA 218.6


Method	Adjustment Buffer	Postcolumn Reaction (PCR) Coil Volume $\mu\text{L}$	Eluent Flow Rate mL/min	PCR Flow Rate mL/min	Injection Volume $\mu\text{L}$
U.S. EPA 218.6	2500 mM $(\text{NH}_4)_2\text{SO}_4$ 1000 mM $\text{NH}_4\text{OH}$	Not specified	1.5	0.5	50-250
Modified Version of U.S. EPA 218.6	250 mM $(\text{NH}_4)_2\text{SO}_4$ 1000 mM $\text{NH}_4\text{OH}$	750	1.0	0.3	1000

- Use a lower-sulfate buffer to adjust sample pH.
- Increase PCR coil to 750  $\mu\text{L}$ .
- Reduce eluent flow rate to 1 mL/min.
- Reduce PCR flow rate to 0.33 mL/min.
- Increase sample size to 1000  $\mu\text{L}$ .

# Effects of Flow Rate, Reaction Coil Volume, and Injection Volume on Hexavalent Chromium Response




**Reaction Coil Volume and Chromate Response**



**Injection Volume and Chromate Response**

# Effect of Sulfate and Chloride on Chromate Peak Response Using Modified Version of U.S. EPA Method 218.6



# Hexavalent Chromium Determination Using Modified Version of U.S. EPA Method 218.6

Columns: Thermo Scientific Dionex IonPac™  
NG1, AS7

Eluent: 250 mM (NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>  
100 mM NH<sub>4</sub>OH

Flow Rate: 1.0 mL/min


Inj. Volume: 1000 µL

Postcolumn Reagent: 2 mM diphenylcarbazide  
10% CH<sub>3</sub>OH  
1N H<sub>2</sub>SO<sub>4</sub>

Reaction Coil: 750 µL

Detector: UV-Vis (530 nm)

Sample: 1.0 µg/L Cr(VI)



17327

**Table 1. Method Detection Limit (MDL) for Chromate  
Based on a 1000 µL Injection**

Chromate Conc. (µg/L)	Std. Dev. (µg/L)	RSD (%)	MDL* (µg/L)
0.1	0.0060	6.996	0.018
0.2	0.0056	3.193	0.018

\* MDL = (Standard Deviation) × (t<sub>s, 99</sub>), where t<sub>s, 99</sub> = 3.14 for n = 7.

**New MDL in reagent water is 0.018 µg/L  
15 × lower than 218.6.**



---


# How to Achieve a Lower Method Detection Limit?

# Achieving a Lower Method Detection Limit (MDL)

---

<b>Method</b>	<b>Column Set</b>	<b>Reaction Coil Volume (μL)</b>	<b>Eluent Flow Rate (mL/min)</b>	<b>Postcolumn Flow Rate (mL/min)</b>
<b>Modified Version of 218.6</b>	Thermo Scientific Dionex IonPac™ NG1 Guard 4 × 50 mm Dionex IonPac AS7 Analytical 4 × 250 mm	750	1.0	0.3
<b>Current Study</b>	Dionex IonPac AG7 Guard 2 × 50 mm Dionex IonPac AS7 Analytical 2 × 250 mm	125	0.36	0.12

# Chromate Determination Using a 2 mm Column Format




Columns: Thermo Scientific Dionex IonPac™ AG7 (2 × 50 mm), AS7 (2 × 250 mm)  
Eluent: 250 mM (NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>  
100 mM NH<sub>4</sub>OH  
Flow Rate: 0.36 mL/min  
Inj. Volume: 1000 µL  
Postcolumn Reagent: 2 mM diphenylcarbazide  
10% methanol  
1N sulfuric acid  
Reaction Coil: 125 µL  
Flow Cell: Semi-micro (PEEK) 2.5 µL


A. 0.1 µg/L Cr(VI) in DI water  
B. 0.1 µg/L Cr(VI) in HIW

0.1 µg/L Cr(VI) in (A) DI water and (B) high ionic-strength water (HIW)

# Chromate Detection




Blank and low-level chromate




Comparison of chromate in high ionic-strength water (HIW) and Sunnyvale tap water

# Chromate in Tap Water



Sunnyvale, CA	0.051 μg/L	San Jose, CA	1.3 μg/L
---------------	------------	--------------	----------

# LCMRL and Method Detection Limit (MDL)



**Table 1. MDL for Chromate in High Ionic-Strength Water Based on a 1000 µL Injection**

Chromate Conc. (µg/L)	Std. Dev. (µg/L)	RSD (%)	MDL* (µg/L)
0.001	0.0003	10.03	0.0009
0.005	0.0004	6.62	0.0013

\* MDL = (Standard Deviation) × ( $t_{s, 99}$ ), where  $t_{s, 99} = 3.14$  for  $n = 7$ .

**Lowest Concentration  
Minimum Reporting Limit = 0.019 µg/L**


**MDL = 1 ppt**

# Postcolumn Delivery Options

---

- Pneumatic delivery—a pressurized chamber that uses pressure to deliver reagent
  - Advantages
    - Lower baseline noise
  - Disadvantages
    - Not software controlled
    - Requires monitoring to maintain accurate flow rate
- Single Piston Pump (e.g., AXP)
  - Reliable flow rate
  - Software controlled
- Dual Piston Pump
  - Flexibility for additional applications
  - Software controlled

# Alternative Mode of Postcolumn Reagent Delivery



Column: Thermo Scientific Dionex IonPac™  
AS7 (2 × 50 mm), AS7 (2 × 250 mm)  
Eluent: 250 mM (NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>,  
100 mM NH<sub>4</sub>OH  
Flow: 0.36 mL/min  
Inj. Vol: 1000 µL  
Postcolumn Reagent: 2 mM diphenylcarbazide  
10% methanol  
1N sulfuric acid  
Reaction Coil: 125 µL  
Flow Cell: Standard (PEEK), 2.5 µL

A. 0.1 µg/L Cr(VI) in DI water  
B. 0.1 µg/L Cr(VI) in HIW

A) DI water and B) high ionic-strength water (HIW) on a Thermo Scientific Dionex ICS-2100 system

**Postcolumn reagent delivered by an AXP pump**



# Detection Limit Using Different Postcolumn Delivery Mechanisms

Method Detection Limit (MDL) Comparison Between Postcolumn Pumps				
Column Format	Chromate Conc. (µg/L)	Pneumatic Pump MDL* (µg/L)	Single Piston MDL* (µg/L)	Dual Piston MDL* (µg/L)
4 mm	0.1	0.018	NA	NA
	0.2	0.018		
2 mm	0.005	0.0013		

\* MDL = (Standard Deviation) × (t<sub>s, 99</sub>), where t<sub>s, 99</sub> = 3.14 for n = 7.

# Column Dimension Comparison

---

Method Detection Limit (MDL) for Chromate Based on a 1000 µL Injection				
Column Format	Chromate Conc. (µg/L)	Std. Dev. (µg/L)	RSD (%)	MDL* (µg/L)
4 mm	0.1	0.0060	6.986	0.018
	0.2	0.0056	3.193	0.018
2 mm	0.005	0.0004	6.62	0.0013

\* MDL = (Standard Deviation) × (t<sub>s, 99</sub>), where t<sub>s, 99</sub> = 3.14 for n = 7.

1 ppt Detection Limit with 2 mm Column Format

# Summary

---

- Modifications proposed: use 2 mm column format and proportional reduction of the flow rates and reaction coil volume.
- Postcolumn reagent delivery can be configured three ways.
- Method detection limit for chromate at 1 ppt will allow a minimum quantitation limit of 3 ppt.
- Modifications allow sufficient sensitivity for determining hexavalent chromium at the proposed California PHG level of 20 ppt.
- We are working with the U.S. EPA Office of Ground Water and Drinking Water to incorporate these configurations into a new U.S. EPA Method 218.7.