

Productivity Enhancements in Microwave Assisted Extractions (EPA Method 3546) Of Semi Volatile Organic Compounds From Environmental Samples

Sumedh Phatak
Milestone Inc



MILESTONE
productivity tools



Agenda

- Introduction to EPA 3546
- Conventional Extraction Methods
- Microwave Assisted Extraction
- Data Comparison
- Improvements To Increase Throughput
- Summary
- Q&A

What is EPA Method 3546?

- Method for extracting water insoluble or slightly soluble organic compounds from soils, clays, sediments, sludge and solid wastes.
- Employs closed vessel microwaves to achieve the extraction of analytes at high temperatures and pressures
- Recovery equivalent to Soxhlet extraction (EPA method 3540) with – less solvent, shorter time and better precision in most cases

Typical Parameters – EPA 3546

Application: Soils, Clays, Sediments,
Sludge, Solid Wastes

Sample Size: 2 – 20 g

Solvent Volume: 25 ml

Temperature: 100 – 115 °C

Pressure: 50 – 150 psi

The need for EPA 3546...

- Extraction demands
- Data quality
- Representative sampling
- Productivity enhancement
- Cost effectiveness

Conventional Techniques

Soxhlet:

- Regarded as the standard method
- Inexpensive equipment
- Able to process large amounts of sample
- Separate filtration is unnecessary



Conventional Techniques

Limitations of Soxhlet Extractions:

- Long operation times - ~4 – 48 hours
- Large solvent consumption – 200 – 500 ml
- Large operation cost
- Only extraction capabilities
- High exposure to toxic materials



Conventional Techniques

Sonication:

- Fast method – 30 – 60 min
- Inexpensive equipment
- Able to process large amounts of sample
- Matrix Independent



Conventional Techniques

Limitations of Sonication Extractions:


- Large solvent consumption – 200 – 500 ml
- Labor intensive
- Filtration required
- High exposure to toxic materials



Microwave Theory

How Microwaves Work:

- Polar molecules rapidly oscillate with the oscillating electromagnetic field.
- Rotating molecules push, pull, collide other molecules and generate heat



Microwave Extraction

- Specially designed microwaves
- Rotors with multiple vessels
- Samples + mixture of solvents added
- Closed vessels ensure no loss of volatiles
- Concentration
- Analysis




Advantages of Microwave Assisted Extraction

- High throughput
- Reduction of operation time
- Minimum consumption of solvents
- Save of space
- Reduction of operation cost
- No expensive glassware
- Broad range of sample sizes
- Flexible platform (multiple chemistries)
- Conforms with EPA's pollution prevention goals

Milestone's Ethos EX

- Chassis
- Safety features
- Color touch-screen controller
- PID Control



MILESTONE
productivity tools

Improvements to Enhance Productivity – Rotors & Vessels

Milestone Ethos EX Rotors - Specifications




| SK-10 | SK-12 | Q-20 | PRO-16/24 | LVR | MUP-41 |
|---|--|--|--|--|--|
| Max. temp & pressure Max. Press. 100 bar Max. Temp. 300°C | Max. temp & pressure Max. Press. 35 bar Max. Temp. 300°C | Max. temp & pressure Max. Press. 40 bar Max. Temp. 250°C | Max. temp & pressure Max. Press. 30 bar Max. Temp. 260°C | Max. temp & pressure Max. Press. 10 bar Max. Temp. 200°C | Max. temp & pressure Max. Press. 30 bar Max. Temp. 260°C |
| Number of Samples 10 | Number of Samples 12 | Number of Samples 20 | Number of Samples 24 | Number of Samples 6 | Number of Samples 41 |
| Vessel Volume 100 mL | Vessel Volume 100 mL | Vessel Volume 45 mL | Vessel Volume 75 mL | Vessel Volume 270 mL | Vessel Volume 70 mL |
| Minimum vol. 10 mL | Minimum vol. 10 mL | Minimum vol. 3 mL | Minimum vol. 10 mL | Minimum vol. 10 mL | Minimum vol. 10 mL |
| Max. sample wt. Up to 20 g | Max. sample wt. Up to 15 g | Max. sample wt. Up to 2 g | Max. sample wt. Up to 15 g | Max. sample wt. Up to 30 g | Max. sample wt. 10 g |
| Vessel Material TFM | Vessel Material TFM | Vessel Material TFM | Vessel Material TFM | Vessel Material TFM | Vessel Material TFM or PFA |
| Vessel Technology Vent-and-Re seal | Vessel Technology Vent-and-Re seal | Vessel Technology Vent-and-Re seal | Vessel Technology Vent-and-Re seal | Vessel Technology Vent-and-Re seal | Vessel Technology Self-regulating |



Productivity Enhancement – High Sample Size & Solvent Volume

- Patented by Milestone
- Cap is held in place by a PEEK spring which deforms as pressure pushes the cap upwards
- Overpressure is gently and controllably released
- Pressure is maintained in the vessel even while venting
- No ejection of contents or loss of volatiles





Vent and Reseal Technology



Productivity Enhancement – Temperature Control

- Fiber Optic Probe
- Direct temperature measurement
- PID Control



Technology Comparison

Soxhlet vs Microwaves

| | Method 3540C | | Microwave Extraction | |
|------------------------|--------------|--------------|----------------------|--------------|
| | 12 Samples | 1000 Samples | 12 Samples | 1000 Samples |
| Labor Cost | \$140 | \$11,666 | \$41 | \$3,383 |
| Solvent Cost | \$108 | \$9,000 | \$8.40 | \$400 |
| Total Cost | \$248 | \$20,666 | \$49 | \$3,783 |
| Total Time Required | 316 hours | 13,167 hours | 3.08 hours | 257 hours |
| Total Solvent Consumed | 5.4 L | 450 L | 420 ml | 15.96 L |

Data Analysis

- Instrument setup
- 15g of Ottawa Sand
- Disposable Glass inserts were used


Run 1:

- 25ml of 1:1 hexane acetone
- 10 min to 100 °C + 20 min at 100 °C

Run 2:

- 25 ml of 1:1 DCM acetone
- 10 min to 100 °C + 20 min at 100 °C

Further Enhancement – Disposable Glass Vials



MILESTONE
productivity tools

Data Run - 1

- Sample: 15 g soil
- Solvent: 15 ml Hexane + 15 ml Acetone
- Dilution Factor: 1
- Final Volume: 1 ml
- Spikes analyzed in duplicates: Aniline, Phenol, 2-Chlorophenol, 1,4-dichlorobenzene, 1,2-dichlorobenzene, 2-methylphenol, 2-nitrophenol, 2,4-dimethylphenol, 4-chloroaniline, 2,4,6-trichlorophenol, 2,4,5-trichlorophenol, 4-nitrophenol, azobenzene, pentachlorophenol, anthracene, fluoroanthene
- Average % Recovery: 78.98%
- Surrogates analyzed: 2-fluorophenol, phenol d-6, nitrobenzene d-5, 2,4,6-tribromophenol, 4-terphenyl – d14
- Average % Recovery: 85.93%

Data Run - 2

- Sample: 15 g soil
- Solvent: 15 ml Hexane + 15 ml Acetone
- Dilution Factor: 1
- Final Volume: 1 ml
- Spikes analyzed in duplicates: Phenol, 2-Chlorophenol, 1,4-dichlorobenzene, 1,2-dichlorobenzene, 2-methylphenol, 2-nitrophenol, 2,4-dimethylphenol, 4-chloroanile, 2,4,6-trichlorophenol, 2,4,5-trichlorophenol, 4-nitrophenol, anthracene, fluoroanthene
- Average % Recovery: 75.45%
- Surrogates analyzed: 2-fluorophenol, phenol d-6, nitrobenzene d-5, 2,4,6-tribromophenol, 4-terphenyl – d14
- Average % Recovery: 75.15%

Advantages of Disposable Glass Vials

- No cleaning step required post extraction
- No need to purchase additional rotors
- Clean environment – low blanks
- Comparable or better recoveries
- Inexpensive

Summary

- Microwave extraction offers significant benefits over traditional extraction techniques.
- Modifications and upgrades are made to enhance productivity – Rotor design, safety features, temperature measurement, software
- The use of disposable glass vials further increase the productivity of a lab – no downtime, no cleaning, great results

Thank You!

Q&A