


The QuEChERS Extraction Approach and Comprehensive Two-Dimensional Gas Chromatography of Persistent Organic Pollutants in Human Breast Milk

Chris Rattray, Michelle Misselwitz, Jack Cochran
and Julie Kowalski




Polychlorinated Biphenyls




- ✓ Lipophilic
- ✓ Bioaccumulate
- ✓ Persistent

Brominated Flame Retardants



Pollutants in Human Milk



**NRDC**
NATURAL RESOURCES DEFENSE COUNCIL
THE EARTH'S BEST DEFENSE

Enter Search

NRDC gets top ratings from

ABOUT USNEWSISSUESPOLICYSMARTER LIVINGSMARTER BUSINESSNRDC.TVBLOGS

ISSUES: HEALTH
Healthy Milk, Healthy Baby
Chemical Pollution and Mother's Milk
Chemicals in Mother's Milk

NEWSREFERENCEEDUCATION

Medscape Today
NEWS

News Search Medscape

DrConnect™ Stay connected to your patients through secure, online information.
Click here to learn more.

Environmental Contaminants in Breast Milk: Persistent Organic Pollutants (POPs)
Authors and Disclosures
Print This

Death Caps
An Ed Rosenberg San Francisco Mystery

Persistent Organic Pollutants (POPs)
POPs are a family of man-made chemicals that are persistent in the environment and can accumulate in the fatty tissues of animals, including humans.

TOXIC POLLUTANTS IN BREAST MILK
For nutrition and infant health, authorities agree that the best food is breast milk.(1) Breast milk provides disease-fighting immune factors that infants' immature immune systems cannot provide for themselves. Compared with bottle-fed babies, those who are breast-fed have fewer infections of the intestine, brain, middle ear, respiratory tract, and urinary tract. They gain more benefit from vaccines,(2) and enjoy greater protection against allergies,(3) type 1 diabetes,(4) and ulcerative colitis and Crohn's disease.(5) Breast-feeding also enhances brain development.(6)

Current Methodology



- Liquid-Liquid Extraction, Silica SPE Cleanup
 - Multiple extractions, time consuming and solvent intense
 - Emulsifications, acid washing



- Focused Microwave Extraction, Silica and Alumina SPE Cleanup
 - Capital Equipment
 - Minimized solvent usage



- Semi-Automated Solid Phase Extraction with Multi-Column automated cleanup
 - Large capital equipment
 - Fast sample throughput

QuEChERS

Quick, Easy, Cheap, Effective, Rugged, and Safe

**PAHs and PCBs
in Seafood**

**Pesticides in Fruits and
Vegetables**

**Pesticides in
Dietary
Supplements**

➤ Simplified shake extraction
for multi- pesticide residue
analysis

➤ Extracted and partitioned
using acetonitrile and
salt/buffer solution

➤ Dispersive solid phase
extraction cleanup

**Veterinary
Drugs in Milk**

**BFRs and PCBs
in Milk??**

1. M. Anastassiades, S.J. Lehotay, D. Stajnbaher, F.J. Schenck, J. AOAC International 86, p. 412-431 (2003).

2. EN 15662, Foods of Plant Origin—Determination of Pesticide Residues Using GC-MS and/or LC-MS/MS Following Acetonitrile Extraction/Partitioning and Clean-up by Dispersive SPE—QuEChERS method.

3. QuEChERS-A Mini-Multiresidue Method for the Analysis of Pesticide Residues in Low-Fat Products. [http:// www.quechers.com](http://www.quechers.com)

4. AOAC Official Method 2007.01, Pesticide Residues in Foods by Acetonitrile Extraction and Partitioning with Magnesium Sulfate.

QuEChERS for POP Monitoring in Human Milk

- ✓ Low initial start-up costs
 - Centrifuge
 - 50 mL PTFE Tubes
- ✓ Low solvent consumption
 - Approx. 50 mL /sample
- ✓ Fast sample preparation
 - Easily 6 samples/person/day
- ✓ Multi-residue method
 - Simultaneous extraction of multiple-compound classes (ie PCB, BFR, OCP)



QuEChERS Extraction Approach

5 mL Milk + 5 mL LC/MS Water



Add 100 μ L IS (PCB 130)

10 mL Hexane:Acetone (1:1)

Shake 1 min



Add Original
Unbuffered
QuEChERS
Packets
4g MgSO₄
1g NaCl

Shake 1 min


Centrifuge 5 min @ 3000g



2.5 mL aliquot for cleanup

Cartridge SPE Cleanup

Evaporate to
“dry”




Add 0.5 mL
Hexane

Condition Silica
SPE Cartridge
3mL/500mg Silica

Add ~ 1cm
of Na_2SO_4
to tube

Elute with 5 mL
Hexane

Concentrate to
50 μL in isooctane



Instrument Conditions

GC X GC Conditions

Injection

Temp: 250°C

Injection: 1 µL fast, splitless (1.0 min)

Liner: Sky Single Taper w/ wool

Column

Rxi-XLB 30m x 0.25mmID x 0.25µm

Rxi-17Sil MS 1.0m x 0.25mm x 0.25µm

Conditions

Oven: 80°C (1 min) to 120°C @

10°C/min to 300°C @ 3°C/min

Secondary Oven Offset: 5°C

Carrier: He, constant flow 1.4 mL/min

Modulator Temperature Offset: 20°C

Modulation Time: 3 sec

LECO Pegasus® TOFMS Conditions

Transfer Line Temperature: 300°C

Source Temperature: 225°C


Electron Ionization: 70 eV

Stored Mass Range: 45 to 750 u


Acquisition Rate: 50 spectra/sec




Basics of GCxGC




GC x GC Analysis




GC x GC Analysis




GC x GC Analysis




GC x GC Analysis




GC x GC Analysis



GC x GC Analysis



GC x GC Analysis



Cow Milk vs Human Milk

Contents	Human	Cow
Protein	1.2	3.3
Lactose	7.0	4.8
Fat	3.8	3.8
Ash (minerals)	0.21	0.71
Total solids	12.4	12.8

Table Written by Harold Eddleman, Ph. D., President, Indiana Biolab, 14045
Huff St., Palmyra IN 47164

Gravimetric Analysis Non Volatile Extractables

NIST (Human Milk): 3.21%
QuEChERS (Cow Milk): 3.39%

Contents	Human	Cow
fatty acids 8C (%)	trace	6
polyunsaturated fatty acids (%)	14	3

Constituents of Human Milk United Nations University Centre <http://unu.edu>

Cow Milk - Percent Recoveries (n=3)

	0.1 ppb	0.5 ppb	1 ppb	10 ppb
PCB 18	138	103	152	90
PCB 52	131	122	163	97
PCB 44	111	107	157	98
PCB 66	108	112	165	99
PCB 101	114	111	162	101
PCB 110	110	112	167	102
PCB 151	109	111	161	103
PCB 153	114	116	154	104
PCB 138	149	118	155	106
PCB 187	118	116	155	106
PCB 183	110	118	154	106
PCB 180	111	124	144	107
BDE 47	186	140	83	98
PCB 170	114	125	138	107
BDE 100	131	167	63	96
BDE 99	130	154	53	94
PCB 206	114	136	101	104
PBB 153	107	175	38	86
BDE 153	119	180	32	82
Average	122	129	126	99
%RSD	16	18	38	7

Blanks

	LC/MS Water	
	ng/kg	pg on-column
PCB 52	3.88	0.97
PCB 99	3.28	0.82
PCB 110	5.2	1.3
PCB 118	4.52	1.13
PCB 153	2.32	0.58
PCB 105	1.48	0.37
BDE 47	67.4	16.9
BDE 100	4.48	1.12
BDE 99	29.4	7.34

Fortified Human Milk – SRM 1954







National Institute of Standards & Technology

Certificate of Analysis


Standard Reference Material® 1954

Organic Contaminants in Fortified Human Milk




GC x GC TOFMS (m/z 67)

Cow Milk




SRM 1954



PSA Cartridge SPE Cleanup

PSA = Primary, Secondary Amine




- ✓ Fatty Acids
- ✓ Sugars
- ✓ Weak Anion Exchange

Condition PSA Cartridge
6mL/ 500mg PSA
with Acetone


Load 2.5 mL sample extract
Pass through cartridge

Elute with 5mL Hexane




GC x GC TOFMS of SRM 1954

No Clean



PSA + Silica SPE



	QuEChERS (n=3) (ng/kg)	NIST Certified (ng/kg)	Percent Error
PCB 28	562	519	8
PCB 74	643	563	14
PCB 66	384	428	10
PCB 99	655	558	17
PCB 110	399	458	13
PCB 118	530	597	11
PCB 146	503	492	2
PCB 153	1080	977	11
PCB 105	360	482	25
PCB 138	982	639	54
PCB 178	514	449	14
PCB 167	478	467	2
PCB 177	470	447	5
PCB 156	516	511	1
PCB 172	578	443	31
PCB 157	450	467	4
PCB 180	896	696	29
BDE 47	2970	2570	16
BDE 100	1220	1280	5
BDE 99	637	739	14
BDE 154	414	464	11
PBB 153	297	474	37
BDE 153	1340	1440	7



National Institute of Standards & Technology

Certificate of Analysis

Standard Reference Material[®] 1953


Organic Contaminants in Non-Fortified Human Milk

SRM 1953


- ✓ Non-fortified milk
- ✓ Target concentration 10x less than fortified sample
 - Final concentration in low pg/ μ L amt

GC x GC ECD


Masses: s1



Masses: s1



Concurrent Solvent Recondensation Large Volume Splitless Injection




No modifications necessary to inject up to 30 μL on Agilent split/splitless injection port


P. Magni, T. Porzano, Concurrent Solvent Recondensation Large Sample Volume Splitless Injection, J. Sep. Sci. 26 (2003) 1491.

Patent No: US 6,955,709 B2.

- ✓ Fast autosampler injection
- ✓ Single taper/gooseneck inlet liner with wool
- ✓ 5 meter guard column / retention gap
- ✓ Initial oven temperature below boiling point of solvent



CSR-LVSI



Instrument Conditions

GC X GC Conditions

Injection

Temp: 250°C

Injection: 10 µL fast, splitless (1.0 min)

Liner: Sky Single Taper w/ wool

Column

IP Deact Guard 5m x 0.53mmID

Rxi-XLB 30m x 0.25mmID x 0.25µm

Rxi-17Sil MS 1.0m x 0.25mm x 0.25µm

Conditions

Oven: 80°C (1 min) to 120°C @

10°C/min to 300°C @ 3°C/min

Secondary Oven Offset: 5°C

Carrier: He, constant flow 1.4 mL/min

Modulator Temperature Offset: 20°C

Modulation Time: 3 sec

LECO Pegasus® TOFMS Conditions

Transfer Line Temperature: 300°C

Source Temperature: 225°C

Electron Ionization: 70 eV

Stored Mass Range: 45 to 750 u

Acquisition Rate: 50 spectra/sec



	QuEChERS (n=3) (ng/kg)	NIST Certified (ng/kg)	Percent Error
PCB 28	46.7	63.1	26
PCB 74	117	149	21
PCB 66	31.6	32.8	4
PCB 99	109	138	21
PCB 110	19.9	12.4	61
PCB 118	198	213	7
PCB 146	42.4	55.6	24
PCB 153	459	480	4
PCB 105	47.8	45.7	5
PCB 138	412	317	30
PCB 178	18.5	22.2	16
PCB 167	15.4	15.8	3
PCB 177	18.1	24.5	26
PCB 156	51.9	59.7	13
PCB 172	11.3	14.8	24
PCB 157	14.5	14.6	0
PCB 180	251	236	7
BDE 47	2700	2200	23
BDE 100	985	839	17
BDE 99	392	342	15
BDE 154	40.6	42.8	5
PBB 153	30.9	36.6	16
BDE 153	968	984	2

QuEChERS Extraction with PSA and Silica SPE Cleanup

○ Overview

- 5 mL Milk + 5 mL Water
- 10 mL Hexane:Acetone (1:1)
- Unbuffered QuEChERS Extraction
- PSA cSPE pass through
- Silica cSPE cleanup
- GC x GC TOFMS analysis

Summary

- ✓ Good agreement with SRM certified values
 - Avg Error = 15%
- ✓ Low solvent usage
 - ~ 50 mL total
- ✓ Cost Effective

Conclusion

- ❖ The QuEChERS extraction approach using Hexane:Acetone is a viable approach for multi-analyte screening in human milk
- ❖ Using a combination of PSA then Silica cSPE provides the needed removal of fatty acid interferences in human milk
- ❖ Using GC x GC allows simultaneous determination of PCBs and BFRs
- ❖ Concurrent solvent recondensation – large volume splitless injection can greatly increase detectability without the need for a specialized injection port. However extract cleanup efficiency may need to be revisited to increase system robustness.

Thank You!

