

A New GC/MS System Designed for Helium Carrier Gas Conservation



Approaches for the Smarter
Use of Helium

Dale R. Walker
GCMS Applications Specialist
August, 7, 2013

Helium Facts

- How it's produced
- Where did it all go?
- What can you do conserve?
- Is there an alternative gas?

Helium Facts

A brief history of Helium

Created by radioactive decay of heavy radioactive elements

- Primarily Thorium and Uranium
- Produces Alpha particles which form He nuclei
- Trapped in Natural Gas Fields – mainly in American Great Plains
- Up to 7% by volume
- Separated from natural gas by low temperature fractional distillation

Helium Facts continued.

National Helium Reserve - Amarillo, Texas (Dept. of Interior)

- Fed by 425 mile pipeline from Kansas to Texas
- Five Helium plants feed the pipeline

1996 Congress decided to phase out the Helium Reserve

- Helium Reserve was \$1.4 Billion in debt in 1996
- Plan is to empty the reserve to recover the funds
- Expected to be empty by 2018
- As of 2012 the Reserve supplied 30% of the worlds Helium

Helium Facts cont.

For many years the U.S. supplied 90% of the world's Helium

- New plant in Algeria can supply all of Europe's demand
- Other plants are now on line in Canada and Poland

U.S. maintains reserves in Kansas, Oklahoma, and Texas

- New plants coming on line – Wyoming, Qatar, and Russia
- New plants are not expected to ease the shortage

Helium Facts cont.

Most Helium is used for cryogenics (MRI's) and welding

- Also controlled atmospheres such as wind tunnels
- Gas Chromatography is a minor consumer

Prices are expected to continue to rise

- From 2002 to 2007 world prices doubled
- From 2007 to 2013 world prices have more than doubled again

Helium Facts continued

What can you do?

There are only two realistic choices:

Option 1: Implement Helium conservation measures. This approach extends the lifetime of existing Helium supplies, reduces operating costs, and ensures business continuity.

Option 2: Switch to a different carrier gas such as Hydrogen.

Option 1: Helium Conservation

Environmental Protection Agency, Region 6 – Houston, Texas

The Environmental Protection Agency (EPA) Region 6 Laboratory is working with manufacturers to find solutions to the ongoing Helium shortage that is affecting many laboratories.


Option 1: Helium Conservation

Approach 1: Manual switching valve

A manual switching valve is connected to both a Nitrogen and Helium cylinder and is manually switched at the end of each day, this requires someone to be physically present to switch the valve.

Option 1: Helium Conservation

An example installation



The EPA tested manual valve switching and found it suitable

Option 1: Helium Conservation

This simple cost saving alternative to switching carrier gases can be implemented immediately, with no time needed for method development, and without any noticeable effects on hardware or performance.

Depending on the total number of instruments and the time between sequences, this procedure has the potential to save laboratories thousands of dollars per year by conserving Helium.

Option 1: Helium Conservation


Approach 2: Automated valve

A better approach is to use an automated switching device that delivers Helium or Nitrogen to the GC flow control modules.

This type of automatic valve switching is preferred to the manual valve because it does not require user interaction. This switch may be programmed to engage or disengage at the beginning or the end of your analytical run or sequence.


Option 1: Helium Conservation

Automated switching valve: SW controlled




Option 1: Helium Conservation

Switching the Carrier with the *Carrier Saver* switch



Option 1: Helium Conservation

Switching the Carrier with the *Carrier Saver* switch



Option 1: Helium Conservation

Switching the Carrier with the *Carrier Saver switch*

Approach 3: Automated Helium conservation

Seamless integrated system includes GC hardware and software control which allows laboratories to integrate programmable Helium conservation into the daily operation of the GC or GC/MS system.



Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

Combines with 7890B Sleep-Wake feature


Automatically switches carrier gas supply to N₂ during GC idle time

TWO Critical Components

- **Carrier gas switching EPC module**
- **7890B Sleep/Wake Feature**


Plus ...

- **Chromatographic grade N₂ gas**



Carrier Gas Switch EPC Module


Ensures Business Continuity



- Built on 5th generation EPC
- Fully integrated and controlled by OpenLAB CDS
- Purge channel prevents cross contamination of gases
- Precise pressure control between tank and GC
- In AWAKE mode, switches between gases within 15-30 min for most detectors including MSD

Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module



Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

Greater reliability

Based on proven 5th generation AUX EPC

Provides Nitrogen substitution when the GC is idle

It also provides additional warnings if set points are not reached.

This safety feature is very important when using Hydrogen. If a leak is detected an auto shutdown occurs and it remains in that state until reset by the user.

Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

Greater performance

Purge channel prevents cross contamination of gases and acts as an intermediate pressure regulator from the tank to inlet EPC to ensure greater analytical precision.

Because there are no changes during the analytical sequence this eliminates the need for method revalidation.

Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

Seamless integration

Fully integrated within the 7890B GC
(Open LAB, Chem Station, Mass Hunter)

Method identifies both gases and uses the appropriate flows
for each gas.

Switching between each gas is easily implemented using new
OpenLAB Sleep/Awake functions.

Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

The new Agilent Open Lab CDS software provides many possible Sleep/Wake program configurations and operates much like a programmable home heating system.


Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

You program the 7890B GC to load the Sleep method when not in operation, this eliminates the use of Helium by using Nitrogen as a carrier gas.

How Does It Work ?

Helium Conservation Mode (Sleep Mode with Nitrogen Carrier)




Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module

You program the Wake method to load just before the start of the workday, this restores helium as the carrier gas so the GC system is completely ready and waiting to run samples.

How Does It Works ?

Normal Operation Mode (Wake Mode with Helium Carrier)



Option 1: Helium Conservation

Introducing the Programmable Helium Conservation Module


For GC and GC/MS systems, it only takes 15 minutes to completely remove all of the Nitrogen from the internal GC flow paths. The GCMS will pass the Nitrogen background test and be available to start the analytical sequence within 20 minutes of loading the wake method.




Programmable Helium Conservation Module

Performance Testing: Switching from Helium to Nitrogen

Running Helium only




Running Nitrogen in the injection port overnight




Running over a weekend with Nitrogen purge in the inlet shows no back-diffusion of Nitrogen into the MS

Performance: No Change in Chromatography After N₂ Carrier Sleep Method (GC/FID)


Day 1 - Original He carrier gas run




Day 2 – First He carrier gas run after overnight N₂ Sleep.M method



Day 3 – First He carrier gas run after overnight N₂ Sleep.M method



Performance: Pass MS Tune Within 15min After Switching From N₂ to He as Carrier (GC/MSD)



Time (min)	Counts of Nitrogen Ion			
	5 mL/min He	Relative to Saturation	2 mL/min He	Relative to Saturation
3	1735168	20.69%	8388096	100.00%
4	1033280	12.32%	4959232	59.12%
5	590080	7.03%	1618944	19.30%
6	354112	4.22%	722944	8.62%
7	228480	2.72%	333696	3.98%
10	56984	0.68%	102576	1.22%
15	9052	0.11%	17080	0.20%

Carrier Gas Savings Calculator

Save Helium by > 90% compared to a typical use

Helium supplies last **MUCH** longer

✓ Gas saver mode: 9 times longer

✓ N₂ standby: another 3.3 times more

✓ Total save: last **30 times** longer



Agilent Technologies

Carrier Gas Savings Calculator

Change values in gray boxes to calculate savings for your operating parameters

Method: ASTM D7504 - Analysis of Trace Impurities in Monocyclic

Column: HP-INNOWax, 30 m x 0.32mm ID x 0.25 um

Gas Flow Conditions

He Carrier Flow (mL/min):	1.2
He Split flow (mL/min):	190
Gas Saver Flow (mL/min):	20
Gas Saver On (min):	3
Run Time(min.):	30
Gas Volume in Cylinder (L):	8000
Runs per Day:	15
He Cylinder Cost (\$):	300
N2 Cylinder Cost (\$):	60

Parameter	No Gas Saver	Gas Saver On	Gas Saver + N ₂ Standby
Daily He Usage (L)	275	30	9
He Cylinder Life (days)	29	265	873
Daily N ₂ Usage (L)	0	0	21
N ₂ Cylinder Life (days)	0	0	381
Yearly He Cost (\$)	3,769	413	125
Yearly N2 Cost (\$)	0	0	57
Yearly Total Gas Cost (\$)	3,769	413	183

Option 2: Hydrogen Carrier

Hydrogen is great for many GC methods such as ECD and FID.

Many people have been successful switching their GCMS methods to Hydrogen carrier gas as well

- Vince Marti at Region 8 lab for 8270 and 8260
 - Vince uses cylinder Hydrogen
- Lisa Wool at Region 6 for Pesticide's, PCB's, TPH
 - Lisa uses a Hydrogen generator

Option 2: Hydrogen Carrier

However, Hydrogen is not the best gas for GCMS methods. Depending on the method requirements

- The main issue with Hydrogen carrier is that it is not inert
- There is some reactivity with specific compounds
- Reactivity with Dichloromethane forms HCL in the GC inlet
(Vince Marti in Region 8 discovered cooler inlet temps avoid this problem)
- BFB tuning anomalies can occur (ion 96)
- DFTPP tuning seems to be less effected

Pertinent Application Notes

Reducing Helium Costs more than 10-fold with the PCT Gas Saver mode

Agilent Technologies, Inc., April 1, 2010

5990-5444EN

Prevent Business Disruptions by Managing your Helium Usage

Agilent Technologies, Inc., April 30, 2013

5991-2256EN

Successfully converting from Helium to alternate gases when analyzing VOC's following EPA 524.2 methodology

Agilent Technologies, Inc., August 2013

Original Designers of Each System

Fred Feyerherm, PhD

Harry Prest, PhD

Jim McCurry, PhD

Bruce Quimby, PhD