

Analysis of Hexabromocyclododecane Diastereomers and Tetrabromobisphenol-A Using a Novel Method Combining Supercritical Fluid Chromatography with Tandem Mass Spectrometry

**Douglas Stevens¹, Lauren Mullin^{1,2}, Kenneth Rosnack¹, Ingrid
Ericson Jogsten², Dawei Geng², Andy Aubin¹, Bert van Bavel²**

¹Waters Corporation

²MTM Research Centre, Örebro University, Sweden


Overview



- Supercritical fluid chromatography (SFC) applicability to hexabromocyclododecane (HBCDD) and tetrabromobisphenol-A (TBBPA) analysis
- Chromatography and MS/MS method development
- Sample analyses results
- New/preliminary results: Chiral separations and additional flame retardants
- Conclusions


What is a Supercritical Fluid?

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®




Liquid/gas

> temp and pressure



Critical point


> temp and pressure



Supercritical fluid


	Diffusivity (cm ² /s)	Viscosity (g/cm x s)
Gas	10 ⁻¹	10 ⁻⁴
Supercritical Fluid	10 ⁻⁴ - 10 ⁻³ Liquid Like	10 ⁻⁴ - 10 ⁻³ Gas Like
Liquid	< 10 ⁻⁵	10 ⁻²

High diffusivity, low viscosity = fast, efficient chromatography



UPC² TQ-S



Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]



HBCDD Isomeric Separation

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]

- HBCDD added to Stockholm Convention list of Persistent Organic Pollutants (POPs)
- Predominant HBCDD diasteromers and enantiomers




- Proportion dependent on sample type

Chromatographic Efficiency: 3 min separation of HBCDD and TBBPA

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]


100 ppb HBCD and TBBPA tetra

05102013rcc3_48 Sm (Mn, 1x2)



Faster Analysis Time and Chromatographic Orthogonality


Waters
THE SCIENCE OF WHAT'S POSSIBLE.®



Method Development

UPC² TQ-S

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]




Chromatography Optimization


	α -HBCDD		β -HBCDD		γ -HBCDD		TBBPA	
	gradient	isocratic	gradient	isocratic	gradient	isocratic	gradient	isocratic
HSS C18 SB 1.8 μ m 3.0x100 mm	●	●	●	●	●	●	●	●
BEH 1.7 μ m 3.0x100mm	●	NA	●	NA	●	NA	●	NA
BEH 2-EP 1.7 μ m 3.0x100mm	●	●	●	●	●	●	●	●

● = peak detected,
baseline separation

● = peak present,
but not baseline
separation

● = peak not
detected






Ionization Type Comparison

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®


[M-H]⁻ Ionization Technique Results



Make-up Flow Composition

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®


Post Column Make-Up Flow Optimization




Method Limit-of-detection (LOD) and quantification (LOQ)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]

**0.100 pg/ μ l
LOD for α - $,$ β - and γ -
HBCDD**




**0.250 pg/ μ l
LOQ for β -HBCDD**



LOD
Peak-to-peak S/N >3

LOQ
Peak-to-peak S/N >10

**0.500 pg/ μ l
LOQ for α -, γ -
HBCDD**



ES⁻ MS/MS conditions

	Compound Name	Parent (m/z)	Daughter (m/z)	Au	Dwell (s)	Cone (V)	Collision (V)	PIC	Comments
1	HBCD	640.6000	78.9000	<input checked="" type="checkbox"/>	0.050	30	20	<input type="checkbox"/>	
2	HBCD	640.6000	80.9000	<input checked="" type="checkbox"/>	0.050	30	20	<input type="checkbox"/>	

	Compound Name	Parent (m/z)	Daughter (m/z)	Au	Dwell (s)	Cone (V)	Collision (V)	PIC	Comments
1	TBBPA	542.5000	78.9000	<input checked="" type="checkbox"/>	0.050	30	50	<input type="checkbox"/>	
2	TBBPA	542.5000	80.9000	<input checked="" type="checkbox"/>	0.050	30	50	<input type="checkbox"/>	
3	TBBPA	542.5000	419.5000	<input checked="" type="checkbox"/>	0.050	30	40	<input type="checkbox"/>	

- Capillary Voltage: 2.0
- Cone Voltage: 30
- Source Temperature: 150 °C
- Desolvation Temperature: 500 °C
- Desolvation Gas: 1000 L/hr
- Cone Gas: 150 L/hr
- Points per peak: 15

UPC² Gradient and Conditions

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]

Time (min.)	% A	% B	Flow Rate (mL/min.)
Initial	98	2	2
2	90	10	2
2.01	90	10	2.5
2.3	98	2	2.5
2.31	98	2	2

- Column: HSS C18 SB (1.8µm 3.0 x 100 mm)
- Column Temperature: 40 °C
- Sample Temperature: 25 °C
- Phases: CO₂ (A) and methanol (B)
- Injection Volume: 1 µL
- Make-up flow composition: 0.1% NH₄OH in 2-propanol
- Make-up flow rate: 0.2 mL/min

Application to Biological Matrices

Sample Analyzed

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®


- Can method separate and detect HBCDD isomers in complex matrices?
- Wide solvent compatibility with SFC so final extract used without solvent exchange
- Extracts of human plasma and whale blubber were donated from prior GC analyses




Human Plasma Detection: Low Level

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®


Quantifier MRM




Confirmatory MRM



Quantifier MRM



Confirmatory MRM



S
A
M
P
S
E
T
A
N
D

Sample	Calculated Concentrations ng/g lipid weight		
	α -HBCDD	β -HBCDD	γ -HBCDD
HS1	4.21	<LOD	<LOD


Calibration Curve Results

Concentration (pg/ μ l)	%RSDs ($n=5$)		
	α -HBCDD	β -HBCDD	γ -HBCDD
0.25	<LOQ	7.9	<LOQ
0.5	5.6	11.6	16.7
1	12.6	4.0	13.1
5*	6.8	2.3	5.3
10	3.0	2.2	3.7
25	5.5	2.0	1.6
50	2.3	2.2	1.8
100	3.7	1.3	1.0
Fit	Linear	Linear	Linear
Weighting	1/x	1/x	1/x
R2	0.998	0.999	0.999
*n=4			


Whale Blubber Detection: High Level

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]


Quantifier MRM




Confirmatory MRM




Quantifier MRM



Confirmatory MRM



S
A
M
P
S
E
A
N
D



Calibration Curve Results

Concentration (pg/ μ l)	%RSDs ($n=5$)		
	α -HBCD	β -HBCD	γ -HBCD
0.25	<LOQ	11.7	<LOQ
0.5	9.8	6.4	15.1
1	16.8	9.4	5.7
5	4.3	1.0	2.3
10	5.4	1.0	4.4
25	3.4	2.2	3.4
50	1.6	0.6	1.3
100	3.4	2.6	3.5
Fit	Linear	Linear	Linear
Weighting	1/x	1/x	1/x
R ²	0.998	0.999	0.999

Ion Ratio Confirmation

						Ion Ratios
	Calculated Concentrations pg/ μ l			α -HBCDD	β -HBCDD	γ -HBCDD
Human Serum Samples	α -HBCDD	β -HBCDD	γ -HBCDD	[1.123]	[1.097]	[1.11]
HS 1	0.575	ND	ND	1.342	ND	ND
Whale Blubber Samples				[1.112]	[1.163]	[1.127]
WB 1	12.88	ND	ND	1.108	ND	ND
WB 2	8.33	ND	ND	1.109	ND	ND
WB 3	47.7	ND	ND	1.117	ND	ND
WB 4	19.8	ND	ND	1.16	ND	ND
WB 5	45.5	ND	ND	1.088	ND	ND
WB 6	85.9	ND	ND	1.003	ND	ND


- Ion Ratio = Quantifier MRM area/Confirmatory MRM area
- Pass = +/- 20% of expected ion ratio
- Expected ratio determined by averaging all points' ion ratio in calibration series

Chiral Separations Preliminary Results

Chiral HBCDD Separation

- Prior LC method 4.0 x 200mm 5µm PM- β -CD column with >15 min run
- UPC² prototype cellulose 2.1 x 150mm 2.5µm column on standards
- (+/-) α - and γ - separated using same co-solvent, (+/-) β - with different co-solvent
- Further method optimization underway


Chiral Separation EtOH co-solvent



(+/-) α -HBCDD in Whale Blubber

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]

α -HBCDD Enantiomers in Whale Blubber




Additional Flame Retardants Preliminary Results

Investigation of Traditionally GC-Amenable Flame Retardants

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]



- Wide variety of FRs monitored in environment and biota
- Polybrominated diphenyl ethers (PBDEs) traditionally require GC-MS



Diverse Flame Retardants on Single Platform: Preliminary Detection

Waters
THE SCIENCE OF WHAT'S POSSIBLE.[®]

Flame Retardant	Current Chromatographic Technique	Optimum Ionization Technique (UPC ²)	Predominantly Observed Ion	Retention Time (min.)
triphenylphosphate	LC	ESI ⁺	[M+H] ⁺	1.01
α -, β -, γ -HBCDD	LC	ESI ⁻	[M-H] ⁻	0.88,1.01,1.33
TBBPA	LC	ESI ⁻	[M-H] ⁻	2.11
tris (1,3-dichloro-2-propyl) phosphate	GC and LC	ESI ⁻	[M-H] ⁻	0.62
tris (2,3-dibromopropyl) phosphate	GC and LC	ESI ⁻	[M-H] ⁻	1.11
BDE 28	GC	APCI ⁻	[M-H] ⁻	0.88
BDE 33	GC	APCI ⁻	[M-H] ⁻	0.67
BDE 47	GC	APCI ⁻	[M-Br+O] ⁻	1.12
BDE 99	GC	APCI ⁻	[M-Br+O] ⁻	1.4
BDE 100	GC	APCI ⁻	[M-Br+O] ⁻	1.43
BDE 209	GC	APCI ⁻	[M-C ₆ Br ₅] ⁻	4.19



Conclusions

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

- Rapid and efficient 3 min separation of α -, β - and γ -HBCDD demonstrated effective for complex samples
- Ability to inject wide range of solvents eliminates need for solvent exchange required for RP-LC
- Chiral separation method development underway for HBCDD enantiomers
- UPC² shows potential for single chromatographic platform for GC and LC amenable flame retardants

Acknowledgements

- Samira Salihovic (human serum extracts) and Anna Rotander (whale blubber extracts), MTM Research Centre, Örebro Universitet, Örebro Sweden
- Waters Corp - Baiba Cabovska, Mike Jones, Ken Fountain and Mark Baynham, Chris Hudalla
- Reference for APCI Ionization of PBDEs:
S. Zhou, E. Reiner, C. Marvin, P. Helm, N. Riddell, F. Dorman, M. Misselwitz, L. Shen, P. Crozier, K. MacPherson and I. Brindle, *Analytical and Bioanalytical Chemistry*, 2010, 396, 1311-1320.