

**Surfactant Analysis in Hydraulic Fracturing  
Flowback Waters Using SPE, 2-D UHPLC  
Followed by Corona CAD, MS/MS and Orbitrap  
HRAM Detection**

Richard Jack, Ph.D.


Director – Vertical Marketing, Environmental and  
Industrial

Lee Ferguson, Ph.D.


Associate Professor - Pratt School of Engineering &  
Nicholas School of the Environment  
Duke University

# Hydraulic Fracturing Technique

- Hydraulic fracturing or fracking is a technique that consists of pumping water into a wellbore to create small fractures into which fluids can migrate and leave openings to release trapped gas
- Fracking facilitates the extraction of natural gas from shale plays in which it is unreachable with conventional technologies



# What is Shale Gas?




is defined as a “dark brown calcareous siltstone, with a sn et al., 1993)



Typical shale in surface outcrop

- Muddy sediments rich in organic matter are deposited in seas and lakes
- As these rocks are buried, organic matter is converted to gas, and porosity is generated
- These rock layers are often the **source rocks** for conventional oil and gas fields

Microscope view



# Hydraulic Fracturing Technique

- With the spreading of fracking, concern is raising about its safety and sustainability for the environment, in particular for air and water quality



- Opponents say there are potential adverse effects on the environment – and perhaps surrounding communities – because of the use of chemicals and large amounts of water that are injected into the subsurface
- Growing concern among the scientific community is that the fracking technology, itself, may result in the leakage of methane into the atmosphere or into water reservoirs



# Hydraulic Fracturing Site in WV



# Produced Water

**High Total Dissolved Solids (60–250,000 mg/L)**

**Chloride, bromide, strontium, barium**

Table A. Typical Composition Data from Several Randomly Selected Flowback Water Streams.

Measurement	--- Randomly Selected Flowback Water Samples ---				
	A	B	C	D	E
pH	5.89	5.83	5.95	5.93	7.0
Sodium, mg/l	54,629	1,477	34,548	43,108	3,310
Calcium, mg/l	15,200	15,680	6,800	3,600	241
Magnesium, mg/l	4,730	1,707	899	6,062	49
Barium, mg/l	98	112	127	547	1
Iron, mg/l	92	60	105	1,274	4
Manganese, mg/l	1.8	1.4	1.7	99.6	na
Bicarbonate, mg/l	195	183	348	415	1,098
Sulfate, mg/l	60	10	20	10	48
Chloride, mg/l	125,000	35,000	68,000	93,000	5,000
Sulfide, mg/l	na	na	na	na	na
Total Dissolved Solids, mg/l	200,006	54,230	110,847	148,016	9,751

# Analysis of Flowback waters – the challenge

---

- Extremely dirty!
- 20x salt conc. Compared to sea water
- Change over time
- Vary by subsurface
- Additives vary due to cost, personal preference and other factors
- Organics – validated methods are not available






# Hydraulic Fracturing Fluid Composition

---

- Mostly Water and Sand
- <2% chemicals



- Hundreds of chemicals in fracking solutions – typically 10-15 used per Fracking operation.
- The EPA has narrowed the list of compounds down to less than 20.
- The EPA is in the process of developing analytical methods

Component	Percentage
Other	0.49%
Surfactant	0.085%
KCl	0.06%
Gelling agent	0.056%
Scale inhibitor	0.043%
pH adjusting agent	0.011%
Breaker	0.01%
Crosslinker	0.007%
Iron control	0.004%
Corrosion inhibitor	0.002%
Biocide	0.001%
Acid	0.123%
Friction reducer	0.088%

# Compounds found in Hydraulic Fracturing Fluids

Typical Chemical Additives Used in Frac Water

Compound	Purpose	Common application
<b>Acids</b>	Helps dissolve minerals and initiate fissure in rock (pre-fracture)	Swimming pool cleaner
<b>Sodium Chloride</b>	Allows a delayed breakdown of the gel polymer chains	Table salt
<b>Polyacrylamide</b>	Minimizes the friction between fluid and pipe	Water treatment, soil conditioner
<b>Ethylene Glycol</b>	Prevents scale deposits in the pipe	Automotive anti-freeze, deicing agent, household cleaners
<b>Borate Salts</b>	Maintains fluid viscosity as temperature increases	Laundry detergent, hand soap, cosmetics
<b>Sodium/Potassium Carbonate</b>	Maintains effectiveness of other components, such as crosslinkers	Washing soda, detergent, soap, water softener, glass, ceramics
<b>Glutaraldehyde</b>	Eliminates bacteria in the water	Disinfectant, sterilization of medical and dental equipment
<b>Guar Gum</b>	Thickens the water to suspend the sand	Thickener in cosmetics, baked goods, ice cream, toothpaste, sauces
<b>Citric Acid</b>	Prevents precipitation of metal oxides	Food additive; food and beverages; lemon juice
<b>Isopropanol</b>	Used to increase the viscosity of the fracture fluid	Glass cleaner, antiperspirant, hair coloring



Source: DOE, GWPC: Modern Gas Shale Development in the United States: A Primer (2009).



## 5.4. Analytical Method Development

### *US EPA: Study of the Potential Impacts of Hydraulic Fracturing on Drinking Water Resources: Progress Report*

**Table 45.** Chemicals identified for analytical method testing activities. Selection criteria for the chemicals included, but were not limited to, frequency of occurrence in fracturing fluids and wastewater, toxicity, environmental mobility, and availability of detection systems for the chemical.

Chemical Class	Chemical Name(s)	CASRN	Purpose in Hydraulic Fracturing	Reason Selected
Alcohols	Propargyl alcohol	107-19-7	Corrosion inhibitor	Toxicity, frequency of use
	Methanol	67-56-1		
	Isopropanol	67-63-0		
	t-Butyl alcohol	75-65-0	Byproduct of t-butyl hydroperoxide	
Aldehydes	Glutaraldehyde	111-30-8	Biocide	Toxicity, frequency of use
	Formaldehyde	50-00-0	Biocide	
Alkylphenols	Octylphenol	27193-28-8	Surfactant	Toxicity, frequency of use
	Nonylphenol	84852-15-3		
Alkylphenol ethoxylates	Octylphenol ethoxylate	9036-19-5	Surfactant	Frequency of use
	Nonylphenol ethoxylate	26027-38-3		
Amides	Thiourea	62-56-6	Corrosion inhibitor	Toxicity
	Acrylamide	79-06-1	Friction reducer	Toxicity, frequency of use, requested by EPA researchers
	2,2-Dibromo-3-nitropropionamide	10222-01-2	Biocide	
Amines (alcohol)	Diethanolamine	111-42-2	Foaming agent	Frequency of use
Aromatic hydrocarbons	BTEX, naphthalene, benzyl chloride, light petroleum hydrocarbons		Gelling agents, solvents	Toxicity, frequency of use, requested by EPA researchers
Carbohydrates	Polysaccharides		Byproduct	Requested by EPA researchers
Disinfection byproducts	Trihalomethanes, haloacetic acids, N-nitrosamines*		Byproduct	Toxicity
Ethoxylated alcohols	Ethoxylated alcohols, C8-10 and C12-18	68954-94-9	Surfactant	Frequency of use

*Table continued on next page*

## 5.4. Analytical Method Development - 2012

**Table 46.** Existing standard methods for analysis of selected hydraulic fracturing-related chemicals listed in Table 45. The EPA will analyze samples using existing methods to determine if the procedure meets the quality assurance criteria for the current study.

Chemical Class	Standard Method*
Alcohols	SW-846 Methods 5030 and 8260C
Aldehydes	SW-846 Method 8315
Alkylphenols	No standard method
Alkylphenol ethoxylates	No standard method
Amides	SW-846 Methods 8032A
Amines (alcohols)	No standard method
Aromatic hydrocarbons	SW-846 Methods 5030 and 8260C
Carbohydrates	No standard method
Disinfection byproducts	DWA Methods 521, 551, and 552
Ethoxylated alcohols	ASTM D7485-09
Glycols	Region 3 Draft Standard Operating Procedure
Halogens	SW-846 Method 9056A
Inorganic elements	SW-846 Methods 3015A and 6020A
Radionuclides	SW-846 Method 9310


\* DWA methods can be found at <http://water.epa.gov/scitech/methods/cwa/index.cfm>. SW-846 Methods can be found at <http://www.epa.gov/epawaste/hazard/testmethods/sw846/online/index.htm>.

<http://www2.epa.gov/hfstudy>

# Alkylphenol Ethoxylate analysis

---

- DeArmond P.D., DiGoregorio A.L. Rapid liquid chromatography-tandem mass spectrometry-based method for the analysis of alcohol ethoxylates and alkylphenol ethoxylates in environmental samples. *Journal of Chromatography A*, Aug 30; 1305: 154-63 (2013).
- Assement ... 2015 <http://www2.epa.gov/hfstudy> discusses the prevalence and use of alkylphenols, nonylphenols in the US



# Ethoxylates reported in FracFocus

---

- **Ethoxylation** is an industrial process in which [ethylene oxide](#) is added to alcohols and phenols to turn it into a [surfactant](#)
- 4-Nonylphenol polyethoxylate
- 2-(2-Butoxyethoxy)ethanol
- 2-(2-Ethoxyethoxy)ethanol
- 2-(2-Ethoxyethoxy)ethyl acetate
- Alcohols, C11-14-iso-, C13-rich, butoxylated ethoxylated
- Alcohols, C11-14-iso-, C13-rich, ethoxylated
- Alcohols, C12-13, ethoxylated
- Alcohols, C12-14, ethoxylated
- Alcohols, C12-14, ethoxylated propoxylated
- Alcohols, C12-14-secondary
- Alcohols, C12-14-secondary, ethoxylated
- Alcohols, C12-15, ethoxylated .....
- ETC

## Prevalence in HF – voluntary disclosure in fracfocus.org

---


**Table C-2. Chemicals reported to FracFocus in 10% or more of disclosures for oil-producing wells, with the number of disclosures where chemical is reported, percentage of disclosures, and the median maximum concentration (% by mass) of that chemical in hydraulic fracturing fluid.**

Name	Number of disclosures	% of disclosures	Median Maximum concentration in HF fluid
Nonyl phenol ethoxylate	2,829	16.3 %	20%
Ethoxylated propoxylated C12-14 alcohols	2,342	13.5%	2%
Poly(oxy-1,2-ethanediyl)-nonylphenyl-hydroxy (mixture)	1668	9.6%	5%

<http://www2.epa.gov/hfstudy>

# Separation



---

- *Prototype*: SEC 80 Size Exclusion Stationary Phase
  - Polymethacrylate-based polymer
  - 5- $\mu\text{m}$  diameter
  - 80Å pore size
- 
- Thermo Scientific Acclaim™ Surfactant Plus RP
    - anionic surfactants (alkylbenzenesulfonates, alkyl sulfates, alkylether sulfates) and
    - cationic surfactants (alkyl quaternary ammonium salts, benzylalkylammonium salts, pyridinium salts, and quaternary imidazolium compounds)
    - nonionic surfactants, and polyethylene glycols (PEGs).



# Detection - Charged Aerosol Detection


- Detection of non-volatile and many semi-volatile analytes
- Consistent response independent of chemical structure
- Applicable to both HPLC and UHPLC
- Compatible with gradient conditions Sub-nanogram sensitivity



# Detection - MS/MS and HRAM


- Quadrupole-equipped HR/AM instruments
- HR/AM analyzer permits parallel detection of all target product ions in one concerted high resolution mass analysis

## A SRM



*Serial* monitoring


## B PRM



*Parallel* monitoring

# Identification of Emerging Contaminants

**High-resolution mass spectrometry  
powered by Orbitrap technology helps contaminants “emerge”:**



## **Target screening**

known contaminants



## **Target screening**

Know contaminants  
with greater sensitivity



## **Non-target screening for unknowns**

Assign molecular and  
structural formula to  
chromatographic  
features


## However:


- Environmental samples are highly complex
- Limits on LC-HRMS of complex samples:
  - Poor chromatographic resolution
  - Matrix effects/ion suppression

## Solution:


- Combine two separation mechanisms into a single analysis

**=> 2D HPLC**






Separation based  
on **Molecular  
Size**



Separation based  
on  
**Hydrophobicity**

Hydrophobicity



Size


# Advantages and Objectives

---

- Higher resolution in chromatography -> takes full advantage of MS technology
- 2D retention time offers another, complementary source of analyte chemical information
- Matrix reduction due to enhance separation: Gain in signal intensity
- Higher confidence in non target analysis due to reduced co-elution
- Develop a comprehensive 2D HPLC-HRMS approach
- Optimize this technique for analysis of nonionic surfactants
- Test the quantitative ability of this method




# How does a 2D HPLC system work?



For effective separation:

- Separation mechanisms have to be **orthogonal**.
- Example: Size and Hydrophobicity or Hydrophilic interaction and Hydrophobicity.
- While eluting from the first column – requires strong retention on the second column


# How does a 2D HPLC system work?




For effective separation:

- Separation mechanisms have to be **orthogonal**.
- Example: Size and Hydrophobicity or Hydrophilic interaction and Hydrophobicity.
- While eluting from the first column – requires strong retention on the second column

## 2D HPLC-HRMS Setup



## 2D HPLC-HRMS Setup




# Two candidates for the first separation

## HILIC

Hydrophilic Interaction Chromatography


- Dionex Acclaim HILIC, 1.0 x 15 cm  
Flow rate: 50  $\mu\text{l}/\text{min}$   
99.5 % ACN in  $\text{H}_2\text{O}$
- Separation based on polarity – apolar elute first
  - High organic mobile phase




## SEC

Size Exclusion

- *prototype* Acclaim SEC80, 4.6 x 300mm  
Flow rate: 350  $\mu\text{l}/\text{min}$   
20% ACN in  $\text{H}_2\text{O}$
- Separation based on size – large molecules elute first



# SEC of OPEO



**OPEO**


for n=1

$C_{16}H_{26}O_2$

Mass: 250.38


m/z: 268.2271 ( $NH_4^+$ -Adduct)

**Chromatograms of OPEO 7-8**



## Advantages for SEC

- Longer run => more RP time
- Compatible mobile phase: No make-up flow required for RP 2<sup>nd</sup> dimension






# Qualitative and qualitative analysis of nonionic surfactants

---


- Nonionic surfactants are complex, polymeric series consisting of alcohol- or alkylphenol hydrophobes that are polyethoylated
- They are applied ubiquitously in consumer, industrial, and commercial applications
- Some nonionic surfactants (alkylphenol ethoxylates) are precursors of known endocrine disruptors
- Analysis of these species by HPLC-MS is very challenging due to difficulty in separations, problems with matrix effects

# Separating Surfactants by 1D-H


- 4 Surfactant standards for optimization:




Octylphenol  
ethoxylate  
(OPEO)



Nonyl phenol  
ethoxylate  
(OPEO)




Lauryl alcohol  
ethoxylate  
(LaOE)




Cetyl alcohol  
ethoxylate  
(CetOE)

Hydrophobicity




## Separation on SEC:




## Separation on RP:



# Separating Surfactants by 2D-HPLC



# Resolving and quantifying mixtures of nonionic surfactants using the Charged Aerosol Detector




# Application to environmental samples

---


## Hydrofracking flowback wastewater sample

- Fracking waters often contain nonionic surfactants as emulsifiers
- To date these surfactants have not been well-characterized in wastewaters discharged from fracking operations






# Hydrofracking wastewater- Results


## Fracking Fluid Sample - Basepeak




## XIC



## Surfactant Standards



- Environmental sample gives comparable surfactant bands
- Retention times



# Conclusions

---

- Two dimensional HPLC combined with HRMS is a powerful tool for resolving micropollutants in complex mixtures
- Size exclusion chromatography combined with reverse phase HPLC provides excellent orthogonality for nonionic surfactant analysis
- Nonionic surfactants can be resolved and identified in hydrofracking wastewaters by 2D-HPLC-HRMS.
- 2D-HPLC-HRMS spectra show lower noise, fewer matrix effects, and less co-elution than 1D HPLC-HRMS analysis
- The developed method can be applied for both target analysis as well as suspect screening of micropollutants in complex environmental samples.

# Questions?

---

**ThermoFisher**  
S C I E N T I F I C

The world leader in serving science