



# Increase Analyte and Solvent Recoveries for US EPA 8082



NEMC July13-17 2015





- US EPA 8082 background
- Overview of PCBs
- Evaporation step after extraction
- Importance of solvent recovery
- US EPA 8082, analytical results
- Analyte recovery and reproducibility
- Summary
- Questions

Jason Wagner





#### Scope of application

- US EPA method 8082 is used to determine the concentration of polychlorinated biphenyls (PCBs) in extracts prepared from many types of solid, tissue, and aqueous matrices.
- A solvent exchange from dichloromethane (DCM) to *n*-hexane is required so the sample can be efficiently analyzed through a gas chromatography / electron capture detector (GC / ECD).



PCBs





### What are PCBs

- Polychlorinated biphenyls are a group of chemicals that contain 209 individual compounds with varying harmful effects.
- From 1930 to 1977 Monsanto marketed Aroclors touting chemical stability, low flammability and electrical insulating properties.
- In this time over 600,000 tons of PCBs were produced in the US<sup>1</sup>
- Many capacitors, transformers and other electrical products, as well as their waste from production, have contributed to massive amounts of PCBs released into the environment

<sup>1.</sup> Breivik, K; Sweetman, A; Pacyna, J; Jones, K (2002). "Towards a global historical emission inventory for selected PCB congeners — a mass balance approach1. Global production and consumption". *The Science of the Total Environment* **290** (1–3): 181–98.

### **BUCHI** Toxicity associated with PCBs



- Mixtures of PCBs have shown toxic, carcinogenic and mutagenic effects in testing and research.
- Main routes of exposure are by contact, consuming contaminated food and by breathing contaminated air.
- PCBs have not been produced or used in the US since the late 1970s but redistribution of PCBs already present in soil and water continues.
- PCBs released to waste water find their way to soil and sediment which can lead to slow leaching and bioaccumulation in fish.

Jason Wagner

EPA Technology Transfer Network. (2000). *Polychlroinated biphenyls (PCBs)(Aroclors)*. Retrieved from http://www.epa.gov/ttnatw01/hlthef/polychlo.html





### Workflow

Sample Preparation	<ul><li>Dependent on sample matrix</li><li>Aqueous, solid or tissue</li></ul>
Extraction	<ul><li>Dependent on sample matrix</li><li>Liquid-liquid, Soxhlet, PSE, MAE, etc.</li></ul>
Concentration	<ul> <li>Parallel evaporation</li> </ul>
Analysis	• GC/ECD





#### Concentration Step

- The extraction step for US EPA 8082 uses dichloromethane (DCM) as a solvent as well as variants with acetone.
- Use of electron capture detection (ECD) requires a solvent exchange because even a small amount of DCM will overwhelm the ECD.
- DCM from the extraction step is concentrated to a residual volume followed by the addition of n-hexane and further concentration to a residual volume to remove any traces of DCM.



### Overview of parallel evaporation





### Utilization of a vacuum source



### Why evaporate under vacuum?

- Lowering the pressure lowers the boiling point
- Not as much heat or cooling is needed
- Nearly a closed system



## Utilization of a vacuum source



### Advantages:

- Protect temperature-sensitive samples
- Speed up evaporation
- Eliminate N<sub>2</sub> usage
- Minimized exposure



# **BUCHI** Overview of parallel evaporation



Syncore<sup>®</sup> vortex evaporator overview



# **BUCHI** Overview of parallel evaporation



Syncore<sup>®</sup> vortex evaporator overview



NEMC July13-17 2015

# **BUCHI** Overview of parallel evaporation



Syncore<sup>®</sup> vortex evaporator overview



# **BUCH** Overview of parallel evaporation



Syncore<sup>®</sup> vortex evaporator overview



# Importance of solvent recovery



"Emissions Reduction Program"

### US EPA Method 3500C and other concentration methods

 Summary of Methods; section 2.2 → "Solvent recovery apparatus is recommended for use in evaporative concentrators. EPA recommends that incorporation type of reclamation system as a method to implement an emissions reduction program."

#### Source:

http://www.epa.gov/osw/hazard/testmethods/sw846/pdfs/3500c.pdf

## **BUCHI** Importance of solvent recovery



Environmental aspects

No air pollution by organic or chlorinated solvents

- Chlorinated solvents are (potential) carcinogens
- Safe laboratory
- Clean atmosphere







Analytical results: Accutest<sup>®</sup> Laboratories

Spiked (0.05/0.5 ng) and concentrated laboratory control samples were analyzed for all compounds mentioned in US EPA 8082

Analytical Method	CAS Number Equivalent	Parameter Name	Spike Level (ng)	% Recovery LCS #1	% Recovery LCS #2	% Recovery LCS #3/ 2:1	% Recovery LCS #4/ 2:1
Surrogate Compounds							
SW8082	877-09-8	% Tetrachloro-m-xylene	0.05	86	88	80	82
SW8082	2051-24-3	% Decachlorobiphenyl	0.05	134	120	114	116
Target Analytes							
SW8082	12674-11-2	Aroclor 1016	0.5	94	96	86	86
SW8082	11096-82-5	Aroclor 1260	0.5	124	113	107	109

\* Recoveries over 100% are due to analytical uncertainty





Which one would you choose to sample for analysis?







6

(5)

4

(3)

(2)

C

### Flushback module



Jason Wagner





### Chilled appendage technology

- Less supervision
- High reproducibility
- Low risk of losing sample
- Optical sensor not required









- Gradient function
- Reproducible results
- Environmentally friendly

- Vacuum pump
- Vacuum controller
- Chiller







- > 95 % solvent recovery
- High analyte recovery
- Elimination of N<sub>2</sub> cost
- Essential reproducibility
- Advanced automation

\*





Questions and answers

Find out more: http://www.buchi.com/us-en/syncore-epa

Additional "Thank You" to: Lionel Thomas Laboratory Manager

Accutest<sup>®</sup> Laboratories