

Combining SPME and Derivatization for Analysis of THC and Metabolite in Surface Waters

Katherine K. Stenerson
NEMC, 2015



sigma-aldrich.com/analytical

1

Agenda

1. What is Solid Phase Microextraction?
2. On-fiber derivatization
3. Quantitative analysis of THC and metabolite in surface waters.

SPME Overview

- Solvent-free extraction technique for nearly any sample or matrix
- Alternative to head-space GC and solid phase extraction (SPE) techniques
- Directly interfaced with GC analysis
- Non-destructive to sample
- Reusable (100+ times)
- Inexpensive
- Fast



Assembled SPME fiber and holder with fiber immersed in a liquid sample.



Manual SPME holder and inlet guide.

The SPME Concept

SPME



Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

The SPME Concept

SPME



Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

The SPME Concept

SPME



Sample Adsorption


Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

The SPME Concept

SPME



Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

The SPME Concept

SPME




Sample Adsorption

Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device

©2003, Sigma-Aldrich, All Rights Reserved

The SPME Concept



Sample Adsorption


Please click on the numbered steps below for an animated sequence of the instruction.

- 1 Drill down septum piercing needle to avoid breakage
- 2 Insert needle into container
- 3 Adjust needle depth for aqueous sampling or headspace sampling
- 4 Extend plunger to expose fiber
- 5 Retract fiber before removing to avoid damaging the fiber.
- 6 Remove SPME Device


©2003, Sigma-Aldrich, All Rights Reserved

SPME Fiber Coating: The Business End

- Not an exhaustive extraction technique
- An equilibrium is set up between analytes dissolved in the sample (solution or gas phase) and in the liquid coating on the fiber.
- The fiber coatings consist of:
 - Polymer films (e.g. PDMS)
 - Particles + binder (e.g. carbons or DVB in PDMS)



Adsorption Mechanism for SPME



Number of Moles of Analyte Extracted by Fiber (n)

$$n = (K_{fs} V_f V_s C_0) / (K_{fs} V_f + V_s)$$

K_{fs} = Distribution constant between fiber and sample

C_f^∞ = Equilibrium concentration on fiber

V_f = Volume of fiber coating

C_0 = Initial concentration of sample

V_s = Sample volume

$$K_{fs} = C_f^\infty V_f / C_s^\infty V_s$$

K \longrightarrow affinity of analyte for stationary phase on fiber

Absorbent vs. Adsorbent Fibers


Absorbent-type fibers (**F**ilm-type fibers)

Analytes are extracted by partitioning

- Liquid phase
- Retains by thickness of coating

Analytes do not compete for sites

Fibers can have high capacity




Adsorbent-type fibers (**P**article-type fibers)

Physically traps or interacts with analytes

- Porous particles
- High surface area

Analytes may compete for sites

Fibers have limited capacity



Derivatization & SPME

- Necessary for some compounds
 - Enhance thermal stability
 - Enhance detection
- Three approaches using SPME & derivatization
 1. Derivatize in-matrix and extract derivatives
 2. Simultaneous extraction & derivatization
 3. On-fiber derivatization after extraction

On-fiber derivatization techniques




Or, put another way..... approaches to derivatization with SPME:


1. Derivatize in-matrix and extract derivatives



2. Simultaneous extraction & derivatization





3. On-fiber derivatization after extraction



SPME & Derivatization;

Example 1: *in-matrix* derivatization followed by SPME

- Amphetamines in plasma, 50 ng/L
- Pre-extraction derivatization with ethyl chloroformate
- SPME headspace extraction & GC-MS analysis





SPME & Derivatization;

Example 2: on-fiber derivatization, reagent loaded pre-extraction

- Analysis of aldehydes in beer
- Load derivatization reagent (PFBOA*) on fiber prior to extraction
- Headspace extraction of aldehydes and derivatization on PDMS-DVB fiber


*pentafluorobenzyl)hydroxylamine hydrochloride



SPME & Derivatization;

Example 3: on-fiber derivatization, reagent loaded post-extraction

- Analysis of resveratrol in red wine
- Extracted resveratrol by immersion, 85 um polyacrylate SPME fiber
- Expose fiber to BSTFA:TMCS after extraction



Fiber Selection & On-Fiber Derivatization

- Some fiber coatings may swell when exposed to certain derivatization reagents
 - Example: PDMS can swell when exposed to silylating reagents
- PEG fiber could be damaged by derivatization reagents that react with hydroxyl groups
- Best coating choices for on-fiber derivatization:
 - Polyacrylate : resists swelling
 - Adsorptive coatings (DVB, Carboxen)

Cannabis



- Schedule 1 substance under federal law
 - High potential for abuse
 - No currently accepted medical treatment in the US
- Legal for recreational and medical use in CO, WA, AK, OR.
- Of interest in effluent as part of sewage epidemiology to study drug use.
- Principal psychoactive constituent is **tetrahydrocannabinol (THC)**.
 - Identified by chemist Raphael Mechoulam in 1963
- Major metabolite is **11-nor-9-Carboxy-THC (THCCOOH)**


Why SPME instead of SPE?

- No organic solvents required
- More easily automated
- Less “hands-on” sample preparation time
- Highly sensitive
- Compatible with existing equipment
(i.e. GC/MS)




SPME Method Development for THC & THCCOOH

Approach: On-fiber derivatization after extraction




Method Optimization Studies; THCCOOH response



Fiber selection


- Best response using adsorbent fibers
- 2 cm DVB/CAR/PDMS was the final selection




Extraction time

- High K value = long extraction time
- 60 minutes chosen for method

Method Optimization Studies; Derivatization





- Silylation using MSTFA; forms TMS derivatives
- Derivatives are stable
- Water can affect reaction



- Derivatization time limited to 30 min.

EI Spectra of TMS Derivatives



Abundance



MW 386

m/z-->
Abundance

Abundance



MW 488

Final SPME method

Fiber: DVB/CAR/PDMS, 2 cm

Sample: 8 mL, pH=2 in 10 mL vial

Extraction: immersion, 60 min w/agitation

Derivatization: 30 min at 60°C, MSTFA (500 µL in 10 mL vial)

Desorption: 260 °C, 3 min

Fiber post-bake: 260°C, 10 min

Analysis: GC-MS/SIM on SLB-5ms (20 x 0.18mm x 0.18µm)


Excess reagent to compensate for water

To reduced carryover:


THC < 0.1%

THCCOOH < 0.01% (usually not detected)

Method Linearity



10-100 ng/L (Feb. 2015)



5-100 ng/L
(April 2015)

Method Detection Level Study

- Deionized water spiked at 10 ng/L
- $n=8$
- Calculated LOD using students t value (99% confidence level) x std. dev. Calculated LOQ using 10 x std. dev.

	Avg. amount measured (ng/L)	Std. Dev. (ng/L)	RSD	Accuracy	LOD (ng/L)	LOQ (ng/L)
THC	9.84	1.12	11%	98%	3.4	11.2
THCCOOH	9.52	1.35	14%	95%	4.0	13.5

RSD = Relative standard deviation


LOD = Limit of detection

LOQ = Limit of quantitation

- Accuracy of > 95% for both compounds
- LODs < 5 ng/L, LOQs < 15 ng/L


10 ppT THC and THCCOOH in deionized water

Abundance



THC - TMS derivative

Abundance




THCCOOH - TMS derivative

Analysis of Surface Water Samples

- Analyzed 2 sample sets
 1. **Creek** in Colorado; adjacent to discharge point for treated effluent from wastewater treatment plant
 2. **Drainage basin** in State College/University Park, PA
- Samples collected at two sampling spots
- Analyses done in duplicate
- MS/MSD run with each set; prepared from one sample for each.
 - 50 ng/L spiking level for creek water set
 - 10 ng/L spiking level for drainage basin set

Sampling sites

Drainage area



Source: opp.psu.edu



Source: bouldercolorado.gov

Results: Drainage Basin Samples

	Site 1 avg. n=2 (ng/L)	RPD	Site 2 avg. n=2 (ng/L)	RPD	Site 1 avg. MS/MSD (ng/L)	Accuracy	RPD
THC	4.1	85%	0	--	9.1*	91%	1%
THCCOOH	0	--	0	--	10.3*	103%	3%

RPD= reproducibility of 2 measurements

MS/MSD = matrix spike and matrix spike duplicate, spiked at 10 ng/L;

**amount reported after sample subtraction*

- THC detected at site 1; just above LOD
- THCCOOH not detected at either site
- MS/MSD good accuracy and reproducibility


Results: Creek Samples

	Site 1 avg. n=2 (ng/L)	RPD	Site 2 avg. n=2 (ng/L)	RPD	Site 1 avg. MS/MSD (ng/L)	Accuracy	RPD
THC	8.8	2%	7.9	1%	42.0*	84%	5%
THCCOOH	4.5	50%	---		41.5*	83%	5%

**less unspiked*

- THC detected at both sites > LOD but < LOQ; good reproducibility of measurement
- THCCOOH detected at site 1
- MS/MSD good accuracy and reproducibility

Results: Creek Samples



Conclusions

- SPME with on-fiber derivatization can be used for simultaneous analysis of THC and the metabolite THCCOOH.
- The detection limit from water was 3-4 ng/L for both compounds with a quantitation limit of 10-15 ng/L. When applied to surface water samples, the method was able to detect both THC and THCCOOH at low ng/L levels.
- The method requires minimal hands-on time, and is cost effective with minimal waste in that it does not require the use of organic solvents or single-use extraction cartridges.
- In addition to THC and THCCOOH, the SPME method could possibly be extended to include additional cannabinoids.

Acknowledgements

Bob Shirey
Yong Chen
Gary Oishi



Thank you!

Questions?

