



# Concentration of Samples Complying to US EPA 8082

## Parallel evaporation technology for high solvent and analyte recoveries


Jason Wagner<sup>1</sup>, Jeff Reid<sup>1</sup>, Lionel Thomas<sup>2</sup>, Natascia Turrà<sup>3</sup>  
<sup>1</sup>BUCHI Corporation, 19 Lukens Drive, New Castle, DE 19720, United States  
<sup>2</sup>Accutest® Laboratories, 500 Ambassador Caffery Parkway, Scott, LA 70583  
<sup>3</sup>BUCHI Labortechnik AG, Meierseggrasse 40, Postfach, CH-9230 Flawil 1, Switzerland

**Introduction**  
US EPA method 8082 is used to determine the concentration of polychlorinated biphenyls (PCBs) in extracts prepared from many types of solid, tissue, and aqueous matrices. A solvent exchange from dichloromethane (DCM) to *n*-hexane is required so the sample can be efficiently analyzed through a gas chromatography / electron capture detector (GC / ECD).

A reproducible concentration of surrogate compounds 2,4,5,6-tetrachloro-*m*-xylene (TCMX) and decachlorobiphenyl (DCB), is necessary to obtain accurate data of the testing compounds, Aroclor 1016 and Aroclor 1260. When performing concentrations in compliance with US EPA 8082 it is important to comply with the EPA's "Emissions Reduction Program."

A collaboration with BUCHI Corporation and Accutest® Laboratories was developed to optimize analyte and solvent recoveries for US EPA 8082. Accutest® Laboratories is one of the nation's largest environmental testing laboratories, combining advanced technology and experienced personnel to deliver "Total Performance You Can Count On." The BUCHI Syncore® Analyst, a parallel vortex evaporator equipped with a Solvent Vapor Recovery system with integrated control (SVR-N), was used to evaporate up to 12 samples simultaneously. A programmed pressure and time profile optimized for a DCM and *n*-hexane mixture eliminates the otherwise required solvent exchange, which drastically saves time and yields high analyte recoveries.

### Experimental set-up




**Method 3500C summary of methods; section 2.2**  
"The resultant extract is dried and concentrated in a Kuderna-Danish (K-D) apparatus. Other concentration devices or techniques may be used in place of the Kuderna-Danish concentrator if they give acceptable results for the intended application."


"NOTE: Solvent recovery apparatus is recommended for use in methods that use Kuderna-Danish or other evaporative concentrators. EPA recommends the incorporation of this type of reclamation system as a method to implement an emissions reduction program."

Source: <http://www.epa.gov/osw/hazard/testmethods/sw846/pdfs/3500c.pdf>

**Cooled appendix technology**  
Due to the cooled appendix technology the sample is automatically concentrated to the predefined volume of either 0.3, 1 or 3 mL (A→B). The appendix is actively cooled by the BUCHI Recirculating Chiller. After concentration the sample is transferred to the sample vial (C→D).



**Flushback effect**  
The Flushback module partially condenses the solvent vapor at the top of the sample vessel generating a continuous rinsing along the glass wall. With this, adsorption of analytes at the glass wall is avoided and high analyte recovery rates are obtained.




**Parallel concentration process**  
Application specific on-site support is offered to effectively prepare your environmental sample.

An environmental benign sample concentration process is performed by applying the following parameters and a programmed vacuum gradient.

Solvent	DCM/ <i>n</i> -hexane
Starting volume	110 mL
Final volume	10 mL
Heating plate	60 °C
Cover temp.	60 °C
Cooling temp.	10 °C
Orbital movement	250 rpm

Twelve samples are concentrated in parallel using the Syncore® Analyst R-12 in only 28 minutes, *i.e.*, 2:20 minutes per sample. The solvents, dichloromethane and *n*-hexanes, are recovered at 95% efficiency without the use of nitrogen.



**Analytical results**  
Spiked (0.5/0.05 ng) and concentrated laboratory control samples were analyzed by GC/MS for several volatile surrogate compounds and target analytes found in US EPA 8082.

Analytical Method	CAS Number Equivalent	Parameter Name	Spike Level (ng)	Percent Recovery LCS #1	Percent Recovery LCS #2	Percent Recovery LCS #3/ 2:1	Percent Recovery LCS #4/ 2:1
Surrogate Compounds							
SW8082	877-09-8	% Tetrachloro- <i>m</i> -xylene	0.05	86	88	80	82
SW8082	2051-24-3	% Decachlorobiphenyl	0.05	134	120	114	116

Target Analytes							
SW8082	12674-11-2	Aroclor 1016	0.5	94	96	86	86
SW8082	11096-82-5	Aroclor 1260	0.5	124	113	107	109

\* Elevated % recovery is due analysis uncertainty

**Conclusion**  
• High solvent recovery > 95%  
• High analyte recovery for PCBs  
• Excellent reproducibility  
• Compliance with waste minimization and pollution prevention  
• Automation and intensification → cost reduction  
• Also applicable for US EPA 8270/625 SVOC, 8081A,B/608 organochlorine pesticides, 8082/608 polychlorinated biphenyls, and 8015 nonhalogenated organics

