

# Passive Monitoring:

A Guide to Sorbent Tube Sampling for EPA  
Method 325

Nicola Watson

July 2015



## EPA 325 - Refinery perimeter monitoring

Revised Federal regulation (CFR 40, part 60 and 63) to be implemented September 2015, compliance within 3 years

- Requires continuous monitoring of vapour-phase organics (specifically Benzene) around the boundary of oil refineries
- US EPA Methods 325 A (Sampling protocol) and 325 B (Laboratory analysis)
- 2-week passive sampling using industry standard sorbent tubes.
- Subsequent analysis is by TD–GC(MS) analysis (MS recommended)

### Target Compounds:

- **Benzene**
- Hazardous air pollutants (HAPs) VOCs
- Vapour-phase organics present in refinery air (light/Middle fuel distillates)






# Monitoring method requirements

- 12–24 monitoring stations round each refinery

Further samplers:

- Replicates
- Blanks
- Calibration
- Alternative sorbents
- Additional shorter-term monitoring for pinpointing fugitive emissions





*Example of monitoring stations on a rectangular site of 750–1500 acres.<sup>2</sup> Monitoring sites (●) are placed just beyond the boundary at 20° intervals. Sources between two monitoring stations and within 50 m of the boundary (★) require that additional monitoring stations (●) are installed.*

# Setting up the field stations

## Monitoring industrial air using passive sorbent tubes

### 325 Field Station™

- Passive (diffusive) samplers deployed around the perimeter under weather proof hoods
- Housing up to five tubes: samples, duplicates, blanks and differing sorbents
- Robust, weather proof shelter
- Sample must sit 1.5-3 meters above the ground





# Diffusive Sampler Theory

## Ficks' 1<sup>st</sup> Law

- Diffusion is a molecular transport property. It is the process by which matter progresses along a concentration gradient until there is an equalisation of concentration within a single phase.
- Thus, as with any chemical system, a diffusing gas is spontaneously adopting its most probable energy distribution in its quest for an even, maximum dispersion and thus maximum entropy.
- The rate of this migration property is measured by its flux (J), which is the quantity of matter passing through a reference surface area per unit time.
- As  $J_x$  is the component of a vector and as matter flows down the concentration gradient away from its source, the coefficient of proportionality, D (the diffusion coefficient) in the matter flux expression must be negative.

$$J_x \propto \frac{dC}{dx}$$

therefore,

$$J_x = -D \frac{dC}{dx}$$

where:

$J_x$  - rate of diffusion ( $\text{molcm}^{-2}\text{s}^{-1}$ );

$dC/dx$  - concentration gradient ( $\text{molcm}^{-2}$ ).

D - diffusion coefficient ( $\text{cm}^2\text{s}^{-1}$ ).




# Application of Fick's First Law to passive samplers

## Axial type samplers

- A passive sampler is essentially a collection medium, either
  - solid sorbent,
  - liquid sorbent,
  - or chemically impregnated inert support,which is separated from the atmosphere of interest by a zone of still air.

The driving force for matter flux across the still-air gap is the induced concentration gradient formed between the sampler opening and the air-adsorbent interface where vapour phase matter is scavenged.






# Application of Fick's First Law to passive samplers

## Axial type samplers

For application of Fick's First Law to a diffusive sampler several simplifying assumptions are necessary:


- Ambient concentration of the analyte at the surface of the monitor ( $C_{amb}$ ) i.e. does not take matter from its surrounding environment faster than it can be replaced
- Zero concentration of the analyte at the surface of the sorbent, i.e. the adsorbent is a zero sink and therefore there is no saturation of the adsorbent ( $C_{ads} = 0$ )
- A linear concentration gradient between the two. steady state conditions always exist



# Passive sampling of volatiles onto sorbent tubes

Passive (diffusive) sampling for workplace and environmental air monitoring

- Vapours migrate across the air gap at a constant “uptake rate”
- Diffusive sampling is a slow process, typically sample for :
  - Occupational Hygiene
  - Workplace exposure
  - Personal exposure
  - Environmental air monitoring (days/ weeks)
- One sorbent ONLY



**Note: Diffusive (passive) sampling is gaining momentum and can drive TD-GC/MS system sales. It won't work with glass or SafeLok tubes. Use standard stainless steel or Silcosteel® tubes**

# Method 325 sampler

## Sample geometry

The uptake rate is directly proportional to adsorbent bed surface area and inversely proportional to diffusion path length. Commercially available passive samplers fall into two main categories

- low uptake rate tube-type devices,
- high uptake rate badge-type devices.



## Passive sampling of volatiles onto sorbent tubes

Uptake rate - a particular analyte being sorbed onto a particular sorbent under a set of monitoring conditions

**Well validated** for ambient air – 100's published uptake rates ISO 16017-2, ASTM D6196, EN 14662-4

**Robust** – Variable ambient conditions (temperature, wind speed, humidity, interferences) have minimal impact on uptake rate

**Low cost** – Samplers are re-usable more than 50 times and are inexpensive to buy and transport

**Versatile** – Sorbent tubes can be used for pumped or passive sampling and offer quantitative sampling & release of compounds over a wide volatility range.



During passive sampling, a diffusion cap is fitted to the sampling end of the tube, while the other end is kept sealed (*Note penclip is optional*)

Application Note 001: Uptake rates for tube-type axial diffusive samplers




## Analysis by TD-GC/GCMS

- Passive sampling sorbent tubes are analysed using a thermal desorption device.
- This works by heating the sample tube and releasing the compounds into a flow of inert gas.
- This process is slow, so a focusing step is needed to retain the compounds.
- Once the sample is focused, it is heated again and rapidly injected into the GC column in a narrow band of vapour.




# Sensitivity Enhancement

- VOCs from 100 L of air or gas can be introduced to the GC column in as little as 100  $\mu$ L of carrier gas
- $10^6$  concentration enhancement means ppt and sub-ppt detection limits depending on detector sensitivity






During stage 1, trapped analytes are desorbed from the heated sample tube and transferred to the electrically-cooled focusing trap



- The TD-100 heated valve is inert and low volume allowing quantitative recovery of high & low volatility compounds plus reactive species
- It also isolates the TD system allowing compliance with standard methods: leak testing, backflush trap desorption, purge to vent, overlap mode, etc.


During stage 2, the trap heats rapidly to transfer/inject analytes into the GC column. Split effluent is quantitatively re-collected on a clean tube ready for repeat analysis



- Repeat analysis of re-collected samples makes it easy to validate analyte recovery through the TD flow path
- A change to the overall VOC profile indicates any bias

# Long-term passive sampling of ambient and industrial air

2-week passive sampling of light hydrocarbons monitored around the perimeter of a major petrochemical installation



## Data interpretation

Two-week diffusive sampling uptake rates for benzene on a variety of sorbents

Sorbent	Uptake rate (mL/min)	Uptake rate (ng/ppm/min)
Carbograph™ 1TD or Carbopack™ B	0.64	2.02
Carbopack™ X	0.61	1.99

To determine the concentration of benzene, five-point calibrations are used to calculate the mass on tube from the peak abundance. The following equation is then used to determine the airborne concentration.

$$\text{Concentration (ppm)} = \frac{\text{Mass of sample on tube (ng)}}{\text{Uptake Rate (ng/ppm/min)} \times \text{Sampling time (min)}}$$

# Summary

- Method 325 is due to be released September 2015, with compliance needed within three years
- Sampling protocol is dependant on the size of the site, with each size bracket requiring a different number of sample sites
- Passive sampling is a widely used technique for occupational monitoring, but is gaining ground in ambient air as an alternative to TO-15 longer term sampling
- Many uptake rates are published for method 325, others can be determined through experimentation.

# Questions?

[www.markes.com](http://www.markes.com)

[enquiries@markes.com](mailto:enquiries@markes.com)

UK: +44 (0)1443 230935

USA: 866-483-5684 (toll-free)

