

Chemometrics for Improved Limits of Detection and Dynamic Range in Flow Injection Analysis with Spectrometric Detection

JONATHAN WOOD

FIALAB INSTRUMENTS INC

SEATTLE, WA

PRESENTED AT NEMC 2016

***FLA**lab[®]*



www.flowinjection.com

Special Thanks

National Environmental Monitoring Conference

Dr. Ilkka Lähdesmäki, Dr. Jamin Hoggard, Adam Bigelow, and
Phillip Janecek

Table of Contents

- **Introduction**
- Current Methodology
- New Technique
- Results
- Conclusion



Introduction – Purpose

- Demonstrate how multivariate chemometric techniques can improve detection limits in common environmental assays
 - FIAlyzer 1000 flow injection analyzer
 - UV-VIS spectrophotometer (CCD array or PDA based)
- Technique leverages the wide spectral range of spectrophotometers and real time data processing algorithms
- 3-5 fold noise reduction in common assays


Introduction – Flow Injection Analysis


- Flow Injection Analysis (FIA):
 - Automated technique for sample injection, reagent mixing, and signal detection by moving streams
 - Sample injected into a continuously moving stream
 - Reagent merged with sample, generating colored product
 - Color intensity measured in detector

Introduction – Flow Injection Analysis



Introduction – Flow Injection Analysis





Introduction – Chemometrics

- Application of data analysis techniques from statistics, mathematics, and computer science to analytical chemistry
- Examples: Multivariate analysis, classification & clustering, curve resolution


Table of Contents

- Introduction
- **Current Methodology**
- New Technique
- Results
- Conclusion

Current Methodology

- As each sample passes through a flowcell, a spectrophotometer measures absorbance:
 - $A = \log_{10} \frac{I_0}{I}$
 - Measured at a primary wavelength – typically peak wavelength of analyte's absorbance spectrum
 - Corrected using a reference wavelength – a wavelength where no absorbance is expected. Fluctuations due to measurement noise.

Reference Correction Example



Current Methodology - Continued

- Absorbance and concentration are related by Beer's Law:
 - $A = \log_{10} \frac{I_0}{I} = \epsilon lc$
 - ϵ = absorptivity, l = optical pathlength, c = concentration
- Construct calibration curve from absorbance and concentration of known standards
- Use calibration curve and absorbance of unknown samples to interpolate concentration

- Plots
 - View by run
 - View by peak
 - Processing profile
- Calibration
 - Main view
- Results
 - Detailed
 - Summary
- Report
- Method Log



Devices



Sample Table



Method Operation




Real Time Data View



Data Analysis

Channel 1 -540nm Channel 2 -580nm Channel 3 -595nm Channel 4 -600nm



Fit model: 1st order polynomial

Fit Parameters

Coeff A: 0.00073

Coeff B: 0.13690

Coeff C: 0

R2: 0.99966

☐ Apply Drift Correction

	Name	Known Concentration	Peak Response	Enabled
▶	Std 1	0	0.0007	<input checked="" type="checkbox"/>
	0.002 ppm N...	0.002	0.0009	<input checked="" type="checkbox"/>
	0.005 ppm N...	0.005	0.0010	<input checked="" type="checkbox"/>
	0.01 ppm N...	0.01	0.0019	<input checked="" type="checkbox"/>
	0.02 ppm N...	0.02	0.0028	<input checked="" type="checkbox"/>
	0.05 ppm N...	0.05	0.0070	<input checked="" type="checkbox"/>
	0.1 ppm N-N...	0.1	0.0142	<input checked="" type="checkbox"/>
	0.2 ppm N-N...	0.2	0.0271	<input checked="" type="checkbox"/>
	0.3 ppm N-N...	0.3	0.0420	<input checked="" type="checkbox"/>
	0.5 ppm N-N...	0.5	0.0699	<input checked="" type="checkbox"/>

Load Data

Save Data

Current Methodology - Improvement

- Current methodology is widely used and simple
- Our application uses a spectrophotometer for a detector
 - Sensitive across a wide spectral range (typical 200nm – 800nm UV-VIS range)
 - Spectral peaks are usually broad (Nitrate FWHM ~80nm)
 - However, we were only monitoring a single wavelength
- Can we improve detection limits by using multivariate data we are already collecting? Yes!

Table of Contents

- Introduction
- Current Methodology
- **New Technique**
- Results
- Conclusion

New Technique

- Goal: Use full absorbance spectrum data rather than single wavelength to reduce noise and increase detection limit
 - Most analytes have a well-defined absorbance spectrum
 - If we know what the absorbance spectrum looks like, we can use this as a model to fit experimental data against
 - Derive model spectrum from a matrix-matched standard
 - Fit all other experimental data against this spectrum
 - From the fitted spectrum, interpolate the absorbance at the wavelength of interest and continue as in the previous methodology

New Technique - Steps


1. Derive model spectrum from known, matrix-matched standard
2. For each spectrum, fit experimental spectrum using model spectrum
3. Interpolate the absorbance at the primary wavelength from fitted spectrum
4. We can now use the interpolated data, rather than raw data, to construct calibration curve and derive maximum absorbance for each sample
5. Continue as in the previous methodology

Table of Contents


- Introduction
- Current Methodology
- New Technique
- **Results**
- Conclusion


Results – Model and Fitted Data, Nitrate

MODEL SPECTRUM




LOW CONCENTRATION FITTED DATA






Results – Single vs Multivariate Data, Nitrate

SINGLE WAVELENGTH DATA, 2ppb



MULTIVARIATE FITTED DATA, 2ppb



Results – Concentration Estimates, Nitrate

True Concentration (ppb)	Single Wavelength Conc (ppb)	Multivariate Concentration (ppb)
0	1.679	0.661
2	3.087	2.489
5	5.960	5.631
10	10.66	10.32
20	19.34	20.21

Table of Contents

- Introduction
- Current Methodology
- New Technique
- Results
- **Conclusion**

Conclusion

- Standard analysis techniques rely on measuring absorbance at a single wavelength of interest
- Our analytical instruments employ spectrophotometers that collect absorbance data across a wide spectral range
- Use the entire absorbance spectrum and multivariate analysis
 - Derive a “model spectrum” from a known matrix-matched standard
 - Fit experimental spectral data against this model spectrum
 - Interpolate the primary wavelength from the fitted spectrum
 - Using interpolated data, construct calibration curve and interpolate sample concentrations as usual
- 3-5 fold reduction in noise and lower detection limits

Thank You!

For More Information:

FIALab Instruments Inc

2151 N Northlake Way

Seattle, WA 98103 USA

sales@flowinjection.com

(206) 258-2290

<http://www.flowinjection.com>