

Adventures in Citizen Science

Reagent Grade Enzymes and Open Source Lab Equipment



Ellen R. Campbell, CEO
Wilbur H (Bill) Campbell, PhD





CITIZEN SCIENTISTS NEED ACCURATE DATA

- Why gather data that is not accurate, reliable, and precise?
- Why use legacy methods that require hazardous reagents and large volumes?
- It's the 21st Century! We need data in digital, exportable formats.
- Accurate optical data requires portable photometers that accept standard 1 mL flat-sided cuvettes. Affordable device needed.

We've combined biotech with an inexpensive photometer to enable *sustainable science*.



NECi Superior Enzymes





REAGENT GRADE ENZYMES FOR CHEMISTRY

- **Selectivity**

“Find” target in complex mixtures

- **Sensitivity**

Low detection limits in complex mixtures

- **Specificity**

False negatives and false positives are rare

- **Safety**

For shipping, storage, handling, and disposal


Enzyme-based analytical chemistry is accurate, reliable, and environmentally benign.



FROM BASIC RESEARCH TO CITIZEN SCIENCE

- How do plants use nitrogen fertilizers?
- How does Nitrate Reductase (NaR) help plants make proteins and DNA?
- Led to investigation into enzyme structure and function
- Development of MAbs to enable purification
- Early studies in molecular biology (1980s)
- Development of nitrate assay using purified NaR (1988)
- Enzyme-based nitrate removal system published in *Nature* (1992)

NITRATE REDUCTASE




ENZYME BASED NITRATE ANALYSIS

- Nitrate reductase (NaR) catalyzes the reduction of **nitrate to nitrite**. The B vitamin **NADH** is cofactor.
- The resulting **nitrite** reacts with the Griess color reagents to produce a **highly colored product**. Absorbance measured at 540nm.
- **NaR and NADH replace the toxic cadmium** or zinc reagents used in conventional nitrate analysis methods.

Minimal change is required from standard methods






Samples are measured against a standard nitrate curve to determine results in ppm Nitrate-N

From Basic Research to Products

- NECi founded in 1993
- R&D funded by the Small Business Innovation program
- NIH funded cloning and expression of NaR
- USDA funded product development: lab reagents and test kits
- NSF funded development of the Photometer
- USDA funded expression and test kits for Phosphate
- NECi's nitrate detection reagents on the market since 1993. Recombinant NaR since 2001.






BENEFITS OF RECOMBINANT NITRATE REDUCTASE

- **Recombinant enzyme design & production**
 - Stringent quality control, specified traits & batch consistency
- **Nitrate reductase simply replaces cadmium**
 - Griess reaction follows nitrate reduction = minimal training to switch
- **Active site only binds nitrate**
 - High specificity, low sensitivity, rare interferences in complex matrices
- **Nitrate reductase naturally occurring enzyme in plants**
 - Non-toxic for lab personnel + no hazardous shipping or waste fees



USGS VALIDATION: ENZYME VS. CADMIUM




Nitrate by NaR vs
nitrate by
Cadmium.

Data courtesy of
USGS in Denver
CO

Groundwater
samples
analyzed by
Discrete
Analyzer

USGS VALIDATION: ENZYME VS. CADMIUM



Source: Colorimetric Determination of Nitrate Plus Nitrite in Water by Enzymatic Reduction, Automated Discrete Analyzer Methods. USGS 2011. <http://pubs.usgs.gov/tm/05b08/contents/TM5-B8.pdf>

Enzymes for simple yet accurate Test Kits

- Start with common validated lab methods
- Develop appropriate enzyme and assay
- Simplify lab methods for ease of use
- Simple extraction with water as solvent
 - *Measure available nitrate and phosphate*
- Sample size is small in relation to assay volume
 - *Generally 20-fold dilution to minimize interferences*
- Less handling = fewer chances to introduce error
 - *In protocol or calculations*



WHICH ONE WILL YIELD BETTER DATA?



ENZYMATIC TEST KITS FOR CITIZEN SCIENCE

- **Fast** Results in 20 minutes
- **Easy** Prepackaged, premeasured, pre-calculated
- **Versatile** Nitrate (and phosphate) in many sample types
- **Accurate** Strong correlation to laboratory data
- **Economical** ~ \$5 per sample (varies with test kit size)
- **Safe** Non-toxic reagents for handling & disposal




VARIOUS SAMPLE TYPES



- Well + Ag Water
 - Standard Range Levels
- Environmental Water
 - Low Range Levels
- Soil
- Plant Petioles
- Green or Dry Animal Forage

ENZYME BASED PHOSPHATE TEST KITS



- Methods X Publication April, 2015
- Project funded by U.S. Dept. of Agriculture
- Test kits now available for soil, low and standard range water, and laboratory analysis

21st CENTURY CITIZEN SCIENCE TECHNOLOGY


Legacy Technology + Test Kits	NECi Photometer + Test Kits
Add Reagents to 5-10 mL Sample	50 µL Sample into Buffer
Cadmium or Zinc Powder	Freeze-dried Enzyme
Round Glass Vials	Flat-Sided Cuvettes
2+ Points of Manual Data Entry/Transfer	Automated Data Input, Storage, Transfer
Acidification or Complex Preparation	Water Extraction Sample Preparation
Manual Sample Detail Entry	Sample Collection Time + GPS Stamp
Bulky, Expensive Equipment	Low-Cost, Handheld Instrument

NSF FUNDED OPEN-SOURCE PHOTOMETER

- Collaborative Effort with Michigan Technological University
- Open Source 3-D Printed Case Design by Dr. Joshua Pierce's Graduate Team
- SBIR/STTR program award #IIP-1417061
- \$65 in materials to build your own
- PLOS One publication summer 2015



SQUARE OPTICS ARE BETTER THAN ROUND!



NECi photometer optical bench

Uses Standard and micro 1 cm cuvettes

Double beam for accurate blanking and background

Dual wavelength LEDs for versatility and low power requirements

Design features lead to increased accuracy






COMMERCIAL HANDHELD PHOTOMETER

- Dual Wavelength for reliable nitrate and phosphate data anywhere
- Data is stored and tracked on Android devices for export and further analysis
- User-friendly, non-toxic test kits with improved mobile app interface
- Applications for alcohol & glycerol in biodiesel, fuel, and beverages in progress



Nitrate Analysis Equipment: Cost vs. Data Quality



METHOD VALIDATION DATA

Interagency round-robin studies showing accuracy + reliability of
nitrate reductase method for nitrate-n analysis



EPA CLEAN WATER ACT ATP SAMPLE MATRICES

Sample Type	NaR Reduction (mg N/L)	Cd Reduction (mg N/L)
Denver area treatment plant Influent wastewater	0.03	0.03
Denver area treatment plant Wastewater effluent #1	7.8	7.6
Denver Area treatment plant Wastewater effluent #2	0.23	0.26
Michigan paper mill waste stream effluent	0.04	0.03
Denver area metal finisher waste stream effluent	270.8	272.6
Denver area Commercial laundry waste stream effluent	4.8	4.8
Environmental Resources Associates #507 Hardness WasteWatR reference material	0.05	0.06
Michigan Confined Animal Feeding Operation (CAFO) effluent from tiled field	13.77	14.1
Low-nutrient seawater (collected offshore Hawaii)	0.027	0.030
ERA #608 Reference Standard	6.8	7.02
USGS PE N-116 (low nutrient-fortified river water)	0.45	0.48
USGS PE N-115 (high nutrient-fortified river water)	2.28	2.36

Source: U.S. EPA Alternate Test Procedure Validation Report for Enzymatic Reduction Method of Nitrate-N Analysis. Data Merged from Tables 1 + 11.



CLEAN WATER ACT EPA VALIDATION STUDY

Table 5 Initial Performance and Recovery (IPR) Summary

Lab	# Analyses	Spike Conc. (mg/L)	Mean Recovery (%)	RSD (%)	Minimum Recovery (%)	Maximum Recovery (%)
2	4	2.5	100.67	0.71	99.67	101.30
3	4	2.5	101.82	0.54	101.29	102.40
4	4	2.0	96.16	2.23	93.50	98.70
5	4	2.5	106.38	0.73	105.82	107.52
7	4	2.0	101.79	0.78	100.63	102.37
8	4	2.5	102.76	0.92	101.77	103.93
9	4	2.0	100.50	1.13	98.97	101.66
10	4	2.0	98.45	0.66	97.77	99.21
11	4	2.7	101.11	0.95	100.17	102.43
12	4	2.0	99.34	2.45	95.85	101.50

ADDITIONAL ANALYTES

- **We've adapted the photometer for alternative wavelengths.**
 - Soil health: microbial respiration test
 - Turbidity
- **We've developed recombinant forms of reagent grade enzymes for new analytes and sample types**
 - Glycerol in biodiesel for fuel QA/QC
 - Ethanol
 - Galactose and other sugars for food QA/QC
 - Peroxide



NECi Superior Enzymes Thanks the Following for Their Work and Support:

Small Business Innovation Research programs of
The US Dept of Agriculture (nitrate and phosphate reagents and test kits)
The National Science Foundation (photometer)

- Michigan's ETF program (travel to conferences)
- Validation partners at ASTM International, US EPA, & US Geological Survey

The people who produced this work:

- Dr JM Pearce and students at MI Technological University
- Dr CJ Patton, USGS in Denver, CO
- NECi Scientists and Technicians

Find more information at www.nitrate.com

