


Canister Cleaning Practices and Blank VOC Concentrations

Jason S. Herrington (presenter), Gary Stidsen, Joe Konschnik, and Steve Kozel

U.S. EPA Method TO-15

- A “*guidance document*” for a Performance Based Method (PBM)
- Means a laboratory can meet performance criteria by hook or crook



TO-15 Blank Criteria

- **8.4.1.6**

- At the end of the evacuation/pressurization cycle, the canister is pressurized to 206 kPa **(30 psig) with humid zero air.** The canister is then analyzed by a GC/MS analytical system. Any canister that has not tested clean (compared to direct analysis of humidified zero air of **less than 0.2 ppbv of targeted VOCs**) should not be used.

“Guidance” for Canister Cleaning

- **8.4.1** Evacuate down to 500 mTorr
- Hold under vacuum for 60 minutes
- Fill with humidified “zero air” to 30 psig
- Repeat cycle two additional times for a total of three cycles

- **8.4.1.8** As an **option** to the humid zero air cleaning procedures, the canisters are **heated** in an isothermal oven not to exceed 100 °C during evacuation of the canister...

Informal Customer Survey

- We do not use a humidified gas...
- We do not use heat...
- We use 14 cycles to clean our canisters...
- We use 6 cycles to clean...
- We use nitrogen, not air...
- Etc...
- The various iterations of cleaning regimens goes on at infinitum...

So why today's talk?

- Assumptions about environmental conditions are

- No risk

- BUT

Compound	1-in-1-million cancer risk (pptv)	Noncancer effects (pptv)
Acrolein	NA	9
Naphthalene	5.6	570
Propylene Dichloride	11	870
1,3-Butadiene	15	900
Acrylonitrile	6.8	920
Ethylene Dibromide	0.22	1200
Methyl Bromide	NA	1300
Benzene	39	9200
Carbon Tetrachloride	26	16000
Vinyl Chloride	44	38000
1,4-Dichlorobenzene	15	133000
Ethylbenzene	92	230000
1,1,2,2-Tetrachloroethane	2.5	NA
Ethylene Dichloride	9.5	NA
Tetrachloroethylene	25	NA
Trichloroethylene	93	NA

Current Study

- Last Summer @ NEMC - Wayne Whipple (U.S. EPA Region 5) and I coincidentally presented on canister cleaning
- Both of us had limited sample sets and/or test parameters
- So for today I evaluated 30 brand new electropolished stainless steel canisters under various cleaning regimens

Objectives

- Time Dependence
 - Only one obscure reference to blank holding times prior to analysis in Method TO-15. Specifically, in section 8.4 “Cleaning and Certification Program” it states **12 hours** of aging.
 - The only other time TO-15 mentions an ageing period is in section 6 “Interferences and Contamination,” where the method states that canisters should be aged for **24 hours**; however, this is for the qualification of “new” canisters.

Objectives (cont'd)

- Sweep Gas
 - Method TO-15 routinely states **"zero air."** However, most laboratories are using nitrogen from the headspace of their liquid nitrogen dewars, which are used for the air concentrators.
- Heat
 - Clearly listed as an **"option,"** which some laboratories are taking advantage of.

Experimental Design

Experiment	Sweep Gas	Humidification	Heat	# of Cycles
N ₂ /NO-H ₂ O/25	Nitrogen	0% RH	25 °C	1
N ₂ /H ₂ O/100	Nitrogen	50% RH	100 °C	12
Air/H ₂ O/100	Air	50% RH	100 °C	12
Air/H ₂ O/ 25	Air	50% RH	25 °C	12
Proprietary	Air	50% RH	100 °C	12

- All canisters filled to 30 psig with 50% RH zero air
- 800 mL sample volumes
- SIM calibrated from 10 to 300 pptv
- Analyzed on Day 0 and 14


Ethylene to Hexachlorobutadiene

Rxi-1 105 m x 0.25 mm x 1.0 μ m

-45 °C for 4.0 min

then 4.5 °C/min to 125 °C for 0 min


then 9.0 °C/min to 250 °C for 0 min



“Usual Suspects”

Compound	W. Whipple	Me
Acetone	X	X
Acrolein	X	X
Benzene	X	X
Benzyl Chloride	X	X
Dichlorodifluoromethane (F12)	X	X
Ethanol	X	X
Hexachlorobutadiene	X	X
Isopropyl Alcohol	X	X
Methylene Chloride	X	X
Naphthalene	X	X
Propylene	X	X
1,2,4-Trichlorobenzene	X	X
2-Butanone (MEK)	X	X
2-Hexanone (MBK)	X	X







Time




Time (cont'd)

- Observation:
 - Most VOCs “grow” significantly with time... not just acrolein
- Explanation:
 - KISS: The VOCs and/or precursors are still on the canister walls at low levels and take time to volatilize and/or react
- Recommendation:
 - Blank aging/holding times need to be more in line with sample turn-around time


VOC Sinks



Sweep Gas




Sweep Gas w/ Time




Sweep Gas (cont'd)

- Observation:
 - Using air as a sweep gas appears to work better relative to nitrogen, especially over time
- Explanation:
 - KISS: Air has oxygen, which carries an oxidative potential nitrogen lacks
- Recommendation:
 - Use air... as suggested years ago

Temperature



Temperature w/ Time



Temperature (cont'd)

- Observation:
 - No clear trend
- Explanation:
 - No KISS...
 - We know what 100 °C means for H₂O. This resulted in an absence of water vapor on the canister walls, thereby allowing polars to stick.
 - Final vacuum/heat stage

Temperature (cont'd)

- Recommendation:
 - Utilize a more moderate temperature (e.g., 70 – 80 °C)
 - Evaluate the potential of a final “cool” evacuation
 - Extend canister cleaning evaluation to silicon-lined canisters

Conclusions/Future Work

- Everything here is very preliminary...
- Time plays a major roll for most VOCs
- Use air as a sweep gas
- Heat needs further evaluation
- Evaluate silicon-lined canisters with new information at hand
- Evaluate more cleaning cycles