




# Improve Data Quality in Environmental Laboratories Using Triple Quadrupole ICP-MS Technology

Maura Rury, Regional Marketing Manager  
NEMC 2017


## As and Se in Environmental Samples

- The role of the As and Se in the environment
- Arsenic: A potential hazard in the food chain



- Many different chemical forms are known and they differ in toxicity and bioavailability
- Plants such as rice are well known for high accumulation of As from soils

- Selenium: An essential nutrient




- Knowledge of Se content in soil may prevent Se deficiency in both human and animal populations

# Handling Interferences with Collision Reaction Cell Technology

- Thermo Scientific™ iCAP™ QCell Technology
- He Kinetic Energy Discrimination (KED)

Analyte	LMCO	Interferences	Precursors
<sup>51</sup> V	35	<sup>35</sup> Cl <sup>16</sup> O, <sup>37</sup> Cl <sup>14</sup> N, <sup>34</sup> S <sup>16</sup> OH	H, N, O, S, Cl
<sup>56</sup> Fe	39	<sup>40</sup> Ar <sup>16</sup> O, <sup>40</sup> Ca <sup>16</sup> O	O, Ar, Ca
<sup>63</sup> Cu	45	<sup>40</sup> Ar <sup>23</sup> Na, <sup>12</sup> C <sup>16</sup> O <sup>35</sup> Cl, <sup>31</sup> P <sup>32</sup> S	C, N, O, Na, P, S, Cl, Ar
<sup>75</sup> As	47	<sup>40</sup> Ar <sup>35</sup> Cl, <sup>40</sup> Ca <sup>35</sup> Cl, <sup>40</sup> Ar <sup>34</sup> SH, <sup>37</sup> Cl <sub>2</sub> H	H, S, Cl, Ca, Ar

- Quadrupole filters out exact mass of target analyte
- Qcell (collision mode) uses pure He to alter kinetic energy




# Handling Interferences with Collision Reaction Cell Technology

- ...and Low Mass Cut Off

Analyte	LMCO	Interferences	Precursors
$^{51}\text{V}$	35	$^{35}\text{Cl}^{16}\text{O}$ , $^{37}\text{Cl}^{14}\text{N}$ , $^{34}\text{S}^{16}\text{OH}$	H, N, O, S, Cl
$^{56}\text{Fe}$	39	$^{40}\text{Ar}^{16}\text{O}$ , $^{40}\text{Ca}^{16}\text{O}$	O, Ar, Ca
$^{63}\text{Cu}$	45	$^{40}\text{Ar}^{23}\text{Na}$ , $^{12}\text{C}^{16}\text{O}^{35}\text{Cl}$ , $^{31}\text{P}^{32}\text{S}$	C, N, O, Na, P, S, Cl, Ar
$^{75}\text{As}$	47	$^{40}\text{Ar}^{35}\text{Cl}$ , $^{40}\text{Ca}^{35}\text{Cl}$ , $^{40}\text{Ar}^{34}\text{SH}$ , $^{37}\text{Cl}_2\text{H}$	H, S, Cl, Ca, Ar

- Quadrupole filters out exact mass of target analyte
- Qcell (collision mode) uses pure He to alter kinetic energy





# Analysis of As and Se by Single Quadrupole ICP-MS

- Determination of trace elements in drinking water, compliant with EPA 200.8 requirements
- Trace contaminants can be quantified in wastewaters and surface waters
- Robust interface to handle samples with high TDS matrices – determination of trace elements in seawater, 25% NaCl
- All analyses can be performed using single-mode analysis – He KED



***But what if the environmental samples have  
Rare Earth Elements present in the matrix?***

# Analysis of As and Se in Environmental Samples

- Ar<sub>2</sub> and ArCl interferences are easily removed using He KED mode analysis
- But if REE are present....

## Single Quad ICP-MS: KED

Typically enhances M<sup>2+</sup> Interferences

Doubly charged ions (e.g. <sup>156</sup>Gd<sup>++</sup>) appear at m/z 78 and hence interfere with the detection of <sup>78</sup>Se<sup>+</sup>

As and Se,  
m/z 75, 78, 80

Rare Earth Elements, m/z 140-176

→ <sup>150</sup>Sm, <sup>150</sup>Nd, <sup>156</sup>Gd, <sup>160</sup>Gd

m/z

# Analysis of As and Se in Environmental Samples

- Ar<sub>2</sub> and ArCl interferences are easily removed using He KED mode analysis
- But if REE are present....

## Single Quad ICP-MS: KED

Typically enhances M<sup>2+</sup> Interferences

### **Solution:**

Mass shift As and Se using O<sub>2</sub>

Doubly charged ions (e.g. <sup>156</sup>Gd<sup>++</sup>) appear at m/z 78 and hence interfere with the detection of <sup>78</sup>Se<sup>+</sup>

As and Se,  
m/z 75, 78, 80



<sup>75</sup>As<sup>16</sup>O<sup>+</sup>, <sup>78</sup>Se<sup>16</sup>O<sup>+</sup>, <sup>80</sup>Se<sup>16</sup>O<sup>+</sup>  
m/z 91, 94 and 96

Rare Earth Elements, m/z 140-176

→ <sup>150</sup>Sm, <sup>150</sup>Nd, <sup>156</sup>Gd, <sup>160</sup>Gd

m/z

# Analysis of As and Se in Environmental Samples

- Ar<sub>2</sub> and ArCl interferences are easily removed using He KED mode analysis
- But if REE are present....

## Single Quad ICP-MS: KED

Typically enhances M<sup>2+</sup> Interferences

### **Solution:**

Mass shift As and Se using O<sub>2</sub>

**Other interferences:** <sup>91</sup>Zr<sup>+</sup>, <sup>94,96</sup>Mo<sup>+</sup>, if present in the sample

Doubly charged ions (e.g. <sup>156</sup>Gd<sup>++</sup>) appear at m/z 78 and hence interfere with the detection of <sup>78</sup>Se<sup>+</sup>

As and Se,  
m/z 75, 78, 80



<sup>75</sup>As<sup>16</sup>O<sup>+</sup>, <sup>78</sup>Se<sup>16</sup>O<sup>+</sup>, <sup>80</sup>Se<sup>16</sup>O<sup>+</sup>  
m/z 91, 94 and 96

<sup>91</sup>Zr<sup>+</sup>, <sup>94,96</sup>Mo<sup>+</sup>

Rare Earth Elements, m/z 140-176


→ <sup>150</sup>Sm, <sup>150</sup>Nd, <sup>156</sup>Gd, <sup>160</sup>Gd

m/z




# Analysis of As and Se in Environmental Samples

- Q1
  - Control ions entering the collision/reaction cell and the remainder of the spectrometer
- Q2
  - Use  $O_2$  to efficiently convert As and Se to AsO and SeO
  - Reactions are selective – the REE<sup>++</sup> ions **do not react**
- Q3
  - Selectively detect AsO and SeO, free from REE<sup>++</sup> interferences



***Are you SURE a triple quadrupole is the best approach for this application?***




- Sample matrix contains 1 ppm Dy, Gd, Nd, Sm, Tb
- Increased BECs observed for all SQ-modes due to unresolved doubly charged REE interferences
- Hydrogen is suitable for

***Reaction with oxygen in TQ mode dramatically lowers BEC concentrations for both As and Se***

# Analysis of As and Se in Environmental Samples

***Are you SURE a triple quadrupole is the best approach for this application?***



- Sensitivity:
  - As: 3953 cps/ppb
  - $^{78}\text{Se}$ : 4443 cps/ppb
- BEC:
  - As = 30 ppt
  - $^{78}\text{Se}$  = 3 ppt
- IDL:
  - As = 0.17 ppt
  - $^{78}\text{Se}$  = 2.0 ppt

## *What about the accuracy?*


### Results for Certified Reference Sample

AGV-1	Content in original sample ( $\mu\text{g}\cdot\text{g}^{-1}$ )	Certified content ( $\mu\text{g}\cdot\text{g}^{-1}$ )
$^{75}\text{As}$	0.892	0.88
$^{78}\text{Se}$	< LOQ	-
Deep Sea Sediment		
$^{75}\text{As}$	1.303	-
$^{78}\text{Se}$	0.109	-

Spike recovery results in samples  
(1 ppb As and Se)

Analyte	AGV-1	Sediment
Arsenic	94.6 %	97.6 %
Selenium	93.4 %	97.6 %

# Analysis of Routine Environmental Samples



Analyte	Measurement Mode	Measured Conc $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
$^{23}\text{Na}$	SQ-KED	5085	5380	95
$^{24}\text{Mg}$	SQ-KED	2665	2540	105
$^{27}\text{Al}$	SQ-KED	55.3	49.5	112
$^{39}\text{K}$	SQ-KED	863	839	103
$^{56}\text{Fe}$	SQ-KED	93.2	91.2	102
$^{59}\text{Co}$	SQ-KED	0.05	0.05	107
$^{60}\text{Ni}$	SQ-KED	0.52	0.48	110
$^{63}\text{Cu}$	SQ-KED	18.2	17.4	105
$^{75}\text{As}$	TQ-O <sub>2</sub>	0.43	0.41	104
$^{78}\text{Se}$	TQ-O <sub>2</sub>	0.10	N/A	N/A
$^{208}\text{Pb}$	SQ-KED	0.08	0.081	99
$^{238}\text{U}$	SQ-KED	0.10	0.09	109

***Br River water reference material - JNR @ \$hR\$-Se!***



# Analysis of Routine Environmental Samples - CRM 1643F

Analyte	Measurement Mode	Measured Conc $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
$^{51}\text{V}$	TQ-O <sub>2</sub>	36.03	$36.07 \pm 0.28$	100
$^{52}\text{Cr}$	TQ-O <sub>2</sub>	18.06	$18.50 \pm 0.10$	98
$^{55}\text{Mn}$	SQ-KED	39.50	$37.14 \pm 0.60$	106
$^{57}\text{Fe}$	SQ-KED	106.11	$93.44 \pm 0.78$	114
$^{58}\text{Ni}$	SQ-KED	55.97	$59.8 \pm 1.4$	94
$^{59}\text{Co}$	TQ-O <sub>2</sub>	24.22	$25.30 \pm 0.17$	96
$^{63}\text{Cu}$	SQ-KED	21.32	$21.66 \pm 0.71$	98
$^{66}\text{Zn}$	TQ-O <sub>2</sub>	87.78	$74.4 \pm 1.7$	118
$^{75}\text{As}$	TQ-O <sub>2</sub>	57.10	$57.42 \pm 0.38$	99
$^{80}\text{Se}$	TQ-O <sub>2</sub>	11.97	$11.70 \pm 0.081$	102
$^{98}\text{Mo}$	TQ-O <sub>2</sub>	123.8	$115.3 \pm 1.7$	107
$^{107}\text{Ag}$	SQ-KED	0.94	$0.97 \pm 0.0055$	97
$^{111}\text{Cd}$	SQ-KED	5.91	$5.89 \pm 0.13$	100
$^{121}\text{Sb}$	SQ-KED	54.59	$55.45 \pm 0.40$	98
$^{125}\text{Te}$	SQ-KED	0.95	$0.977 \pm 0.0084$	97
$^{205}\text{Tl}$	SQ-KED	6.38	$6.892 \pm 0.035$	93
$^{208}\text{Pb}$	SQ-KED	18.19	$18.488 \pm 0.084$	98



# Analysis of Routine Environmental Samples – Seawater Reference NRC CASS-5


Analyte	Measurement Mode	Measured Conc $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
$^{51}\text{V}$	TQ-O <sub>2</sub>	1.37	1.32 $\pm$ 0.14	104
$^{52}\text{Cr}$	TQ-O <sub>2</sub>	0.069	0.106 $\pm$ 0.013	65
$^{55}\text{Mn}$	SQ-KED	2.56	2.62 $\pm$ 0.20	98
$^{57}\text{Fe}$	TQ-O <sub>2</sub>	1.49	1.44 $\pm$ 0.11	103
$^{63}\text{Cu}$	SQ-KED	0.378	0.380 $\pm$ 0.028	99
$^{75}\text{As}$	TQ-O <sub>2</sub>	1.17	1.24 $\pm$ 0.09	94
$^{98}\text{Mo}$	TQ-O <sub>2</sub>	10.80	9.82 $\pm$ 0.72	110
$^{111}\text{Cd}$	SQ-KED	0.0227	0.0215 $\pm$ 0.0018	106

# The Advantage of Triple Quadrupole Technology

- Single Quadrupole: KED
  - Removes polyatomic interferences
- Triple Quadrupole: TQ-O<sub>2</sub>
  - More specific interference removal


# The Advantage of Triple Quadrupole Technology

- Single Quadrupole: KED
  - Removes polyatomic interferences
- Triple Quadrupole: TQ-O<sub>2</sub>
  - More specific interference removal




# The Advantage of Triple Quadrupole Technology

- Single Quadrupole: KED
  - Removes polyatomic interferences
- Triple Quadrupole: TQ-O<sub>2</sub>
  - More specific interference removal



*Use of O<sub>2</sub> as a reactive gas improves detection sensitivity for <sup>75</sup>As in comparison to KED*

# The Advantage of Triple Quadrupole Technology – As Speciation



## iCAP TQ ICP-MS


Parameter	Value	
Mode	TQ-O <sub>2</sub>	
Gas Flow	100% O <sub>2</sub> @ 0.4 mL·min <sup>-1</sup>	
Cell Settings	QCell Bias	Quad Bias
	-7.5 V	-12 V

## Ultimate 3000 HPLC System

Parameter	Value
Column	Thermo Scientific Dionex AS7, 2x250mm
Mobile Phase	Ultra pure water, 20mM (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub> ; 200mM 20mM (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub>
Gradient	5-20mM in 7 minutes, 200mM for 3 minutes
Injection Volume	25μL

# The Advantage of Triple Quadrupole Technology – As Speciation

Species shown here are AsB and DMA at 20 ng·L<sup>-1</sup>




## **SQ Mode with He KED:**

Peak area precision 7-10% (N=5)

LOD: 3.6 ppt (AsB)

5.8 ppt (DMA)



## **TQ Mode with O<sub>2</sub>:**

Peak area precision < 5% (N=5)

LOD: 2.4 ppt (AsB)

1.9 ppt (DMA)



# The Power of Triple Quadrupole Technology

- Elemental impurities in Ni alloys
- Ti, Cr in high purity sulfuric acid
- As in Vitamin B12 (high Co matrix)
- Cd in the presence of high Mo concentrations
- As, Se in samples containing rare earth elements
- Measure nanoparticles at decreasing diameters
- P, Ti in high Si matrix
- Ti in human serum
- As, Cr, V in high purity hydrochloric acid
- S, P in steel and high concentrations of iron



Find out more:  
[thermofisher.com/iCAPTQ](https://thermofisher.com/iCAPTQ)

**The possibilities are endless!**