

ThermoFisher SCIENTIFIC

Kick Interference to the Curb: Tackle ICP-OES Applications with Ease

Maura Rury, Regional Marketing Manager NEMC 2017 What are spectral interferences?

- Produced when emission spectra from matrix components overlap with those of the analyte/s
- They result from interactions between components in the sample and the sample matrix
- The severity of the interferences is dependent on the analyte wavelengths being used





Spectral Interferences

Examples

- Background shift (flat, raised or lowered baseline)
- Background shift (sloped baseline)
- Partial peak overlap
- Full, complete peak overlap





Spectral Interferences

Do you need to correct for these?

- Spectral interferences alter the magnitude of the signal that reaches the detector
- The magnitude of the signal is related to the concentration present in the sample
- If the signal increases/decreases as a result of an interference, the instrument will produce an incorrect concentration



Correcting for interferences – probably a good idea!



Tackling Spectral Interferences

Use One or More Correction (IEC) Factors

- Just follow these simple steps:
 - Identify all elements that have interferences
 - Identify all elements that are causing interferences
 - Use carefully-prepared solutions to allow the instrument to measure the spectral overlap and calculate an accurate IEC factor
 - Repeat process for each interference that must be corrected
 - Re-calculate all IEC factors if any conditions change (plasma parameters, sample matrix, analyte mixture)



Isn't there another way to correct for interferences?



Use Careful Background Point Selection

- Background correction should mimic shape of background emission
- Avoid overcorrection or correction on a nearby emission peak





Tackling Spectral Interferences

Use Careful Wavelength Selection

- Careful wavelength selection ideally, interference free!
- Element Finder plug-in for Thermo Scientific[™] Qtegra[™] ISDS Software
 - Automatically selects interference-free wavelengths
 - Eliminates interferences before you know they exist





Spectral Interference Example

Analysis of Cd

- Where do we want it?
 - Rechargeable NiCd batteries
 - Paint pigments ("cadium yellow" or "cadmium red")
 - He-Cd lasers
 - Corrosion-resistant plating materials
- Where don't we want it?
 - Airborne particles
 - Drinking water
 - Fruits/vegetables grown in Cd-contaminated soil
 - Rice grown in Cd-contaminated fields
 - Food oils





Analysis of Cd – Theory









Analysis of Cd – What Happened?





Addressing Interference on Cd



Here's How Element Finder Helps:

- 1. Hey! You have an interference!
- 2. It's probably from Fe!
- 3. What do you want to do about it?
 - 1. Calculate an IEC? not recommended!
 - 2. Find an alternative wavelength? Let me help you!
- 4. Cd has a strong emission wavelength at 214 nm; Fe doesn't emit there try that one!



Another Spectral Interference Example

- Analysis of B in a high Fe sample matrix asymmetric peak for boron at 249.773 nm
- Clear interference on boron peak (visible at 249.780 nm and 249.790 nm)





Analyze three solutions:

- 1. The sample
- 2. A single element solution of the analyte at the concentration in the sample
- 3. A single element solution of the suspected interference at the concentration in the sample
 - If this is not the interfering analyte, further analysis might be needed to determine this

Then

• Measure the interference to calculate a correction factor

OR

• Select an alternative wavelength that is free from interferences





Automatic Interference Identification





Automatic Identification – Additional Interferences Found!





Automatic Interference Identification

| Wavelengths suggested for Analyte Elements that are suited to your analysis | | | | | | | | | | | |
|---|------------|---------|-------------|--------------|---------------------------|-----------------|--|--|--|--|--|
| 5 | Wavelength | Order T | Intensity 🔻 | Preference T | Available for analysis | Measure Mode | Remarks | | | | |
| В | 249.773 | 135 | 4,000,000 | Automatic | No | Radial | Interference with: Fe, 249.782; Interfer | | | | |
| Boron | 249.678 | 135 | 2,000,000 | Automatic | No | Radial | Interference with: Fe, 249.653; Interfer | | | | |
| Analyte | 208.959 | 461 | 1,500,000 | Automatic | Yes | Radial | | | | | |
| | 208.893 | 461 | 750,000 | Automatic | No | Radial | Interference with: Fe, 208.412; | | | | |
| | 182.641 | 484 | 660,000 | Automatic | No | Radial | | | | | |
| | 182.591 | 484 | 290,000 | Automatic | No | Radial | | | | | |
| | 181.837 | 485 | 73,000 | Automatic | No | Radial | | | | | |
| | 4 | | | | | | • | | | | |





17

Final Spectral Interference Example

Trace Elements in Pentanol

- Multiple interference challenges
 - 1. Physical
 - Differences in nebulization and/or transport efficiency between the standards and samples
 - 2. Chemical
 - Differences between the behavior of standards and samples when in the plasma
 - Easily ionized element (EIE) effects
 - Plasma loading
 - 3. Spectral
 - Atomic emission overlaps from other elements in the sample
 - Molecular emission from solvent





Final Spectral Interference Example

Trace Elements in Pentanol

- Multiple interference challenges
- 1. Physical
 - Differences in nebulization and/or transport efficiency between the standards and samples

2. Chemical

- Differences between the behavior of standards and samples when in the plasma
 - Easily ionized element (EIE) effects
 - Plasma loading
- 3. Spectral

Element Finder

- Atomic emission overlaps from other elements in the sample
- Molecular emission from solvent

Matrix match standards to samples

Use an internal standard

Carefully select sample introduction

Reduce sample uptake to reduce matrix loading



Full Frame of Pentanol Sample



- Emission from other analytes
- Emission from sample contaminants
- Molecular emission from solvent (CO, CO₂, CN)





Vanadium wavelength scan







Pentanol Results

- Successful tactics to address interferences
- Results: Good precision and accuracy, and spike recoveries ~100%

| Element and Wavelength (nm) | View | Spike Concentration (mg kg ⁻¹) | Measured Spike Concentration (mg kg ⁻¹) | Spike Recovery (%) | RSD on three Replicates of the Spike (%) | MDL (µg kg⁻¹) |
|-----------------------------------|--------|--|---|-----------------------|--|------------------|
| Ag 328.068 | Radial | 2.47 | 2.45 | 99.2 | 1.6 | 3.5 |
| AI 167.079 | Axial | 2.47 | 2.82 | 114.2 | 0.2 | 1.4 |
| As 189.042 | Axial | 2.49 | 2.58 | 103.6 | 0.2 | 7.6 |
| Ba 455.403 | Radial | 2.47 | 2.53 | 102.4 | 1.8 | 0.3 |
| Ca 393.366 | Radial | 2.47 | 2.46 | 99.6 | 1.9 | 0.2 |
| Cd 214.438 | Axial | 2.47 | 2.69 | 108.9 | 0.3 | 0.3 |
| Cr 267.716 | Radial | 2.47 | 2.48 | 100.4 | 1.6 | 4.8 |
| Cu 324.754 | Radial | 2.47 | 2.41 | 97.6 | 2.0 | 2.4 |
| Fe 238.204 | Radial | 2.47 | 2.50 | 101.2 | 1.6 | 4.8 |
| Hg 184.950 | Axial | 2.48 | 2.59 | 104.4 | 0.6 | 2.6 |
| K 766.490 | Radial | 2.47 | 2.41 | 97.6 | 2.4 | 55 |
| V 309.311 | Radial | 2.47 | 2.46 | 99.6 | 1.7 | 2.1 |



Kick Interference to the Curb for Successful ICP-OES Results

- Thoughtful method development will ensure that physical, chemical and spectral interferences are identified and corrected
- Wavelength scans and full frame images (captured in seconds) will help illustrate nearby emission from other analytes and sample matrix components and the solvent
- Background shifts and nearby spectral interferences can often be addressed with careful background correction
- When background correction is insufficient on its own:
 - Spend your valuable time and lose productivity calculating IEC factors
 OR
 - Let Element Finder identify possible interferences and automatically select alternative wavelengths

 Imagine Automated Method Development!

 New Element Finder plug-in for Qtegra Intelligent

 Scientific Data Solution sofware with ICP-OES analysis.

 Find your elements

 thermofisher.com/ICP-OES



