

The Coefficient of Determination (R^2) vs Relative Standard Error (RSE)

AUSTIN HABICH

FIALAB INSTRUMENTS INC

SEATTLE, WA

PRESENTED AT 2017 NEMC CONFERENCE

***FLA**lab®*



www.flowinjection.com

Introduction – Purpose


- To explain why trusting R^2 could result in inaccurate data
- To introduce a new calibration metric that avoids overlooking inaccurate data
- To illustrate the points above with an example
- Provide recommendations on each technique

Current Methodology – R^2

- Coefficient of Determination:

$$R^2 = 1 - \frac{\sum e_i^2}{\sum (y_i - \bar{y})^2}$$

- where y_i = known concentration at point i , \bar{y} = average concentration, e_i = error at point i
- Well defined acceptable limits;
For example: $R^2 \geq 0.999$



Disadvantage of R^2

- The magnitude of error e_i at point i matters

$$R^2 = 1 - \frac{\sum e_i^2}{\sum (y_i - \bar{y})^2}$$

- More concentrated standards have a larger affect R^2




Disadvantage of R^2

Table 1: 20% Absolute Error
Comparison

| Standard [=] ppm | Error: e_i | Error ² e_i^2 |
|---------------------|-----------------|-------------------------------|
| 1 | 0.2 | 0.04 |
| 5 | 1 | 1 |
| 20 | 4 | 16 |

Channel 2 -650nm Channel 1 -878nm



Fit model: 1st order polynomial

Fit Parameters

Coeff A: -0.00011

Coeff B: 0.21906

Coeff C: 0

R^2 : 0.99987

RSE (%): 219.58848

☐ Apply Drift Correction

| | Name | Peak Response | Known Concentration | Calculated Concentration | % Error | Enabled |
|---|---------|---------------|---------------------|--------------------------|---------|-------------------------------------|
| ▶ | 0 ppm | 0.0015 | 0 | 0.01 | NA | <input checked="" type="checkbox"/> |
| | 1 ppb | 0.0014 | 0.001 | 0.01 | 589.07 | <input checked="" type="checkbox"/> |
| | 2 ppb | 0.0007 | 0.002 | 0 | 80.09 | <input checked="" type="checkbox"/> |
| | 5 ppb | 0.0029 | 0.005 | 0.01 | 178.09 | <input checked="" type="checkbox"/> |
| | 10 ppb | 0.0024 | 0.01 | 0.01 | 16.39 | <input checked="" type="checkbox"/> |
| | 50 ppb | 0.0097 | 0.05 | 0.04 | 10.44 | <input checked="" type="checkbox"/> |
| | 100 ppb | 0.0251 | 0.1 | 0.12 | 15.05 | <input checked="" type="checkbox"/> |
| | 500 ppb | 0.1097 | 0.5 | 0.5 | 0.28 | <input checked="" type="checkbox"/> |
| | 1 ppm | 0.2132 | 1 | 0.97 | 2.63 | <input checked="" type="checkbox"/> |
| | 3 ppm | 0.6487 | 3 | 2.96 | 1.27 | <input checked="" type="checkbox"/> |
| | 5 ppm | 1.1013 | 5 | 5.03 | 0.56 | <input checked="" type="checkbox"/> |

Load Data

Save Data

Example

| | Name | Peak Response | Known Concentration | Calculated Concentration | % Error |
|---|-------------|---------------|---------------------|--------------------------|---------|
| ▶ | 0ppm N-NO3 | 0.0012 | 0 | 0.08 | NA |
| | 1ppm N-NO3 | 0.0112 | 1 | 0.96 | 4.33 |
| | 5ppm N-NO3 | 0.0568 | 5 | 4.95 | 1.08 |
| | 20ppm n-NO3 | 0.2291 | 20 | 20.02 | 0.08 |

- Vary the calculated concentration of the 1 ppm standard, observe R^2
- Vary the calculated concentration of the 20 ppm standard, observe R^2
- Compare the results




Figure 1: Standard Curve

At 15% error




Figure 2: 1 ppm error




Figure 3: 20 ppm error

At 25% error




Figure 4: 1 ppm error




Figure 5: 20 ppm error

At 30% error




Figure 5: 1 ppm error




Figure 7: 20 ppm error

Summary of Example

Table 2: 20 ppm vs. 1 ppm R² Comparison

| Absolute error | 1 ppm R ² | 20 ppm R ² |
|----------------|----------------------|-----------------------|
| 0% | 0.999995 | 0.999995 |
| 15% | 0.99989 | 0.99814 |
| 25% | 0.99978 | 0.99358 |
| 30% | 0.99971 | 0.98911 |

Summary of Example

- With the same relative error, more concentrated samples have a greater negative effect on the R^2
- R^2 may be deceptive in that, a high percent error of a low standard can still yield a high coefficient of determination




Figure 8: 20 ppm vs. 1 ppm R^2 Comparison

New Technique: Relative Standard Error

- Uses normalized error, so magnitude does not produce bias for concentrated standards

$$\%RSE = 100 * \sqrt{\frac{\sum_i^n \left(\frac{e_i}{y_i}\right)^2}{n - 2}}$$

Where

n = number of calibration points;

e_i = error at point i ;

y_i = known concentration at point i



New Technique: Relative Standard Error

Table 3: Absolute Error Comparison [15%]

| <i>Nominal Concentration [=]ppm</i> | <i>Absolute Error</i> | <i>Calculated Concentration [=]ppm</i> | $\left(\frac{e_i}{y_i}\right)^2$ |
|---|-----------------------|--|----------------------------------|
| 0.1 | 0.015 | 0.085 | 0.0225 |
| 5 | .75 | 4.25 | 0.0225 |
| 20 | 3 | 17 | 0.0225 |

New Technique: Relative Standard Error

- Normalized error reduces bias for RSE
- Concentrated standard's error is normalized by its concentration value

→ Its error has as much weight as less concentrated standards

- Rule of Thumb for good calibration
 - *Good Calibration: $RSE < 10\%$*




Figure 9: RSE for 20 ppm and 1 ppm variation

RSE / R² Comparison

| | Absolute error | R ² | RSE |
|--------------|----------------|----------------|--------|
| 20 ppm Error | 0% | 0.99995 | 2.92% |
| | 15% | 0.99814 | 11.00% |
| | 25% | 0.99358 | 17.92% |
| | 30% | 0.98911 | 21.41% |
| 1 ppm Error | 0% | 0.99995 | 0.71% |
| | 15% | 0.99989 | 10.63% |
| | 25% | 0.99978 | 17.69% |
| | 30% | 0.99971 | 21.23% |

Recommendations

- R^2 heavily weights concentrated standards relative to less concentrated standards of the same percent error
 → A good R^2 may not mean all standards are acceptable
- To avoid this, one of the following actions can be taken:
 - ✓ Reduce dynamic range
 - ✓ Implement two channels for low and high concentrated standards
 - ✓ Observe error of each standard compared to known value

| | Name | Peak Response | Known Concentration | Calculated Concentration | % Error |
|---|-------------|---------------|---------------------|--------------------------|---------|
| ▶ | 0ppm N-NO3 | 0.0012 | 0 | 0.08 | NA |
| | 1ppm N-NO3 | 0.0112 | 1 | 0.96 | 4.33 |
| | 5ppm N-NO3 | 0.0568 | 5 | 4.95 | 1.08 |
| | 20ppm n-NO3 | 0.2291 | 20 | 20.02 | 0.08 |

Conclusion

- RSE does normalize error terms to yield a less biased result
- RSE can catch errors in low standards since all error terms are relative to their standards
- Rule of Thumb for good calibration
 - ✓ *Good Calibration: $RSE < 10\%$*
 - ✓ Observe error of each standard compared to known value

Thank You!

For More Information:

FIAlab Instruments Inc

2151 N Northlake Way

Seattle, WA 98103 USA

sales@flowinjection.com

(206) 258-2290

<http://www.flowinjection.com>

