

EPA 625: Base, Neutral, Acid Semi-volatiles in Municipal and Industrial Waste Water by Solid Phase Extraction

Tom Hall, Fluid Management Systems, Watertown, MA



Overview

- Development of a fully automated extraction system for EPA 8270/625
- Establishment of an extraction procedure capable of implementation of various aqueous matrices
- Participation in ILI SPE Demonstration project for method suitability
- Implementation of platform across other EPA methodologies

Matrices

- DI and Tap
- ASTM 5909 synthetic Waste Water
- TCLP Fluid
- Pond, River and Reservoir surface water
- Real world industrial effluents and influents
- Groundwater



SPE Waste Water Study Participation

- Participated in both Phase I and II
- Submitted data for all matrices in each respective study
- 3rd Party participation using FMS equipment



Comparison of LLE/CLE vs. Automated SPE Methods

LLE/CLE

Open to laboratory background

Uses >360mls solvent

Shaking / Continuous process

Forms emulsions requiring centrifuging

Little Selectivity

Requires water removal

Automated SPE

Closed system

Uses <60mls solvent

Filtration process

No emulsions formed

Wide Selectivity (adsorbent)

In-line water removal

Fast and Unattended



Comparison of LLE/CLE vs. Automated SPE Methods

LLE/CLE

No Separation of waste

More volume to evaporate

Massive solvent emission

CLE uses a lot of power

Requires lots of solvent for cleaning

Automated SPE

Separates Aqueous and Organic Waste

<60mls solvent to evaporate

6 times less solvent emission

Easily Capture Solvent

Lower solvent costs

Lower Disposal Costs

Reduced Solvent Usage



FMS TurboTrace[®] ABN

Touch Screen Control

SD Card for Storing and Transferring
Methods

Six Elution and Wash Solvents

Positive Pressure Pump for precise,
consistent delivery of Sample and Solvent

Sample Size 2 ml to unlimited


Automatic Bottle Rinse

Real Time Graphical Display of
Extraction Steps


Automatic Liquid Level Sensor

ABN Cartridge

Carbon Cartridge



Dual SPE Cartridge Extraction



Direct to Vial
Concentration

SuperVap Features


- 6 (250ml) and 12 (50ml) position models for extractions.
- Dry bath heating element
- Independent secondary heater for extract nipple (Can be disabled).
- Sensor controlled
- Savable temperature log capability.

EPA 625/8270 Procedure

- Fully automates all processes, not manual steps
- Capable of delivering up to 3 independent fractions (or combine as 1)
- Can elute direct to FMS SuperVap evaporator
- Automated organic solvent rinse of sample bottle
- Can elute either cartridge independently or in parallel
- Modular and Expandable
- Direct to GC Vial extracts




Recoveries by Analyte Class




Recoveries for Individual Compounds with specific Analyte Classes




PAHs




Phthalates



Phenols




Ion Exchangers



Semi-Volatiles

◆ TCLP Fluid
 ■ River Water
 ▲ Influent
 ✕ Waste Water



Method Validation

Proficiency Data for Laboratory Validation



IDC/IPR Data

Entry Identifier	IDC #1		IDC #2		IDC #3		IDC #4		EPA Limit
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	
1,4-Dichlorobenzene-d4	20	100	20	100	20	100	20	100	
N-Nitrosodimethylamine-d6	20.75	83	21.11	84	20.22	81	20.33	81	
N-Nitrosodimethylamine	24.69	99	25.47	102	24.91	100	25.29	101	
Pyridine	24.44	98	26.04	104	25.45	102	25.75	103	
Phenol-d6	21.07	84	21.23	85	20.37	81	20.93	84	
Phenol	26.78	107	28.03	112	27.56	110	27.34	109	5-112
Aniline	25.37	101	27.59	110	26.23	105	26.26	105	
Bis(2-chloroethyl)ether-d8	19.3	77	19.77	79	19.67	79	18.69	75	
Bis(2-chloroethyl)ether	24.32	97	26.18	105	24.99	100	25.97	104	12-158
2-Chlorophenol-d4	19.6	78	20.05	80	19.52	78	19.48	78	
2-Chlorophenol	24.77	99	26.31	105	26.24	105	25.97	104	23-134
1,3-Dichlorobenzene	21.13	85	21.19	85	21.01	84	20.32	81	D-172
1,4-Dichlorobenzene	21.32	85	22.1	88	21.82	87	20.19	81	20-124
1,2-Dichlorobenzene	21.55	86	21.84	87	21.87	87	20.49	82	32-129
Benzyl Alcohol	28.6	114	29.15	117	29.16	117	28.71	115	
2,2'-oxybis(1-chloropropane)	27.13	109	27.65	111	26.77	107	26.68	107	36-166
2-Methylphenol	25.96	104	27.49	110	26.31	105	26.76	107	
4-Methylphenol-d8	20.95	84	21.18	85	21.74	87	21.63	87	

IDC/IPR Data

Entry Identifier	IDC #1		IDC #2		IDC #3		IDC #4		EPA Limit
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	
3/4-Methylphenol	26.59	106	27.33	109	26.7	107	27.02	108	
Hexachloroethane	22.29	89	21.84	87	22.07	88	20.52	82	40-113
N-Nitroso-di-n-propylamine	26.88	108	26.89	108	27.43	110	26.8	107	D-230
Naphthalene-d8	20	100	20	100	20	100	20	100	
Nitrobenzene	26.71	107	27.95	112	27.76	111	27.09	108	35-180
Nitrobenzene-d5	21.06	84	21.14	85	21.02	84	20.55	82	
2-Nitrophenol-d4	23.84	95	23.04	92	23.64	95	23.14	93	
Isophorone	27.63	111	28.14	113	29.13	117	27.89	112	21-196
2-Nitrophenol	28.05	112	29.36	117	29.97	120	28.35	113	29-182
2,4-Dimethylphenol	26.91	108	28.22	113	27.32	109	26.28	105	32-119
Bis(2-chloroethoxy)methane	27.1	108	29.2	117	29.36	117	28.2	113	33-184
2,4-Dichlorophenol-d3	20.56	82	21.03	84	20.97	84	21.37	85	
2,4-Dichlorophenol	27.14	109	28.35	113	27.39	110	27.27	109	39-135
1,2,4-Trichlorobenzene	21.7	87	21.57	86	21.78	87	20.81	83	44-142
Naphthalene	22.63	91	23.53	94	23.66	95	22.69	91	21-133
4-Chloroaniline-d4	20.76	83	21.9	88	20.05	80	21.24	85	
4-Chloroaniline	25.48	102	27.21	109	25.49	102	25.69	103	
Hexachlorobutadiene	20.55	82	21.81	87	21.73	87	20.53	82	24-116
Acenaphthene-d10	20	100	20	100	20	100	20	100	
2-Methylnaphthalene	23.73	95	23.87	95	24.17	97	23.29	93	

IDC/IPR Data

Entry Identifier	IDC #1		IDC #2		IDC #3		IDC #4		EPA Limit
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	
1-Methylnaphthalene	23.33	93	23.12	92	24.37	97	23.46	94	
4-Chloro-3-methylphenol	28.55	114	29.94	120	30.63	123	28.77	115	22-147
Hexachlorocyclopentadiene	21.87	87	22.14	89	22.03	88	20.55	82	
2,4,6-trichlorophenol	25.9	104	27.01	108	30	120	28.24	113	37-144
2,4,5-Trichlorophenol	27.47	110	26.86	107	29.38	118	28.89	116	
2-Chloronaphthalene	23.16	93	23.69	95	24.38	98	23.38	94	60-118
2-Nitroaniline	30.05	120	30.91	124	32.64	131	31.4	126	
1,4-Dinitrobenzene	32.72	131	33.68	135	35.76	143	34.63	139	
Dimethylphthalate-d6	19.42	78	19.77	79	21.57	86	20.63	83	
Dimethylphthalate	26.03	104	25.79	103	28.38	112	26.87	107	D-112
1,3-Dinitrobenzene	28.16	113	29.11	116	29.95	120	30.35	121	
Acenaphthylene-d8	17.81	71	18.31	73	19.72	79	18.74	75	
Acenaphthylene	23.51	94	23.71	95	25.11	100	24.24	97	33-145
2,6-Dinitrotoluene	28.46	114	28.78	115	30.64	123	30.38	122	50-158
1,2-Dinitrobenzene	26.02	104	27.77	111	29.21	117	28.86	115	
Acenaphthene	23.4	94	23.36	93	24.55	98	23.23	93	47-145
3-Nitroaniline	25.73	103	26.93	108	27.36	109	27	108	
2,4-Dinitrophenol	27.14	109	29.52	118	30.94	124	30.93	124	D-191
Dibenzofuran	23.43	94	23.7	95	24.81	99	24.16	97	
2,4-Dinitrotoluene	31.85	127	32.86	131	35.06	140	33.97	136	39-139
4-Nitrophenol-d4	25.73	103	25.43	102	28.2	113	25.06	100	
4-Nitrophenol	24.13	97	23.77	95	24.23	97	23.29	93	D-132

IDC/IPR Data

Entry Identifier	IDC #1		IDC #2		IDC #3		IDC #4		EPA Limit
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	
2,3,4,6-Tetrachlorophenol	28	112	29.36	117	31.24	125	29.17	117	
2,3,5,6-Tetrachlorophenol	29.6	118	29.26	117	32.17	129	31.82	127	
Fluorene-d10	18.96	76	18.06	72	18.19	73	17.55	70	
Fluorene	23.32	93	23.83	95	25.13	101	25.01	100	59-121
4-Chlorophenyl phenyl ether	22.1	88	22.11	88	24.34	97	23.41	94	25-158
Diethylphthalate	26.23	105	26.62	106	29.25	114	27.23	109	D-114
4-Nitroaniline	22.62	90	22.24	89	28.2	113	21.16	85	
Phenanthrene-d10	20	100	20	100	20	100	20	100	
4,6-Dinitro-2-methylphenol-d2	13.41	54	13.85	55	14.41	58	15.09	60	
4,6-Dinitro-2-methylphenol	27.15	109	28.13	113	28.54	114	29.74	119	D-181
Diphenylamine	22.39	90	22.22	89	24.11	96	24.56	98	
Azobenzene	25.33	101	26.85	107	26.32	105	27.48	110	
4-Bromophenyl phenyl ether	22.54	90	23.28	93	23.75	95	24.09	96	53-127
Hexachlorobenzene	20.46	82	21.09	84	21.96	88	22.21	89	D-152
Pentachlorophenol	27.19	109	27.83	111	29.47	118	28.81	115	14-176
Phenanthrene	23.17	93	24	96	24.6	98	24.9	100	54-120
Anthracene-d10	21.28	85	21.4	86	20.79	83	21.2	85	
Anthracene	24.76	99	25.31	101	25.87	103	25.82	103	27-133
Carbazole	28.52	114	28.93	116	30.79	123	30.7	123	
Di-n-butylphthalate	25.12	100	26.82	107	27.43	110	26.16	105	1-118
Fluoranthene	26.21	105	26.85	107	27.73	111	27.48	110	26-137
Pyrene	25.28	101	26.27	105	26.97	108	26.97	108	52-115
Benzyl butyl phthalate	28.78	115	29.41	118	30.4	122	30.71	123	D-152

IDC/IPR Data


Entry Identifier	IDC #1		IDC #2		IDC #3		IDC #4		EPA Limit
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	
Benzidine	17.25	69	17.25	69	13.16	53	6.18	25	
Chrysene-d12	20	100	20	100	20	100	20	100	
Benzo[a]anthracene	24.19	97	24.17	97	23.65	95	22.69	91	33-143
Chrysene	24.06	96	24.75	99	23.97	96	22.77	91	17-168
Bis(2-ethylhexyl)adipate	27.2	109	27.51	110	26.97	108	26.19	105	
3,3-Dichlorobenzidine	10.85	43	12.86	51	9.16	37	9.41	38	D-262
Bis(2-ethylhexyl)phthalate	24.16	97	25.11	100	23.56	94	23.25	93	8-158
Di-n-octyl phthalate	26.81	107	27.02	108	26.64	107	25.43	102	4-146
Perylene-d12	20	100	20	100	20	100	20	100	
Benzo[a]pyrene-d12	4.7	19	4.43	18	4.51	18	4.53	18	
Benzo[a]pyrene	25.13	101	24.45	98	24.41	98	23.83	95	17-163
Benzo[b]Fluoranthene	26.83	107	26.34	105	26.06	104	26.24	105	24-159
Benzo[k]fluoranthene	25.53	102	25.93	104	25.29	101	24.78	99	11-162
Indeno[1,2,3-cd]pyrene	7.28	29	25.76	103	24.12	96	26.37	105	D-171
Dibenz[a,h]anthracene	29.79	119	29.54	118	27.91	112	29.57	118	D-227
Benzo[g,h,i]perylene	26.84	107	26.77	107	26.25	105	26.44	106	D-219

Entry Identifier	MDL #1		MDL #2		MDL #3		MDL #4		MDL #5		MDL #6		MDL #7		MDL #8		EPA Limit	STD DEV	MDL	EPA
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec				
1,4-Dichlorobenzene-d4	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00	
N-Nitrosodimethylamine-d6	2.66	53	3.23	65	3.07	61	3.83	77	3.71	74	3.46	69	2.77	55	3.36	67		0.365455	1.10	
N-Nitrosodimethylamine	2.69	54	3.27	65	3.18	64	4.18	84	4.16	83	3.7	74	3.49	70	3.83	77		0.399404	1.20	
Pyridine	3.05	61	3.17	63	3.33	67	4.14	83	4.55	91	4.09	82	3.87	77	4.51	90		0.536732	1.61	
Phenol-d6	2.15	43	2.62	52	2.76	55	3.25	65	3.18	64	3.38	68	2.39	48	2.59	52		0.383642	1.15	
Phenol	2.72	54	3.18	64	3.32	66	4.02	80	3.93	79	4.14	83	3.52	70	3.76	75	5-112	0.365051	1.09	1.5
Aniline	2.52	50	2.72	54	2.81	56	3.66	73	3.58	72	3.54	71	3.48	70	3.67	73		0.406794	1.22	
Bis(2-chloroethyl)ether-d8	2.74	55	3.16	63	3.25	65	4	80	3.95	79	3.87	77	2.92	58	3.28	66		0.438102	1.31	
Bis(2-chloroethyl)ether	2.9	58	3.29	66	3.47	69	4.23	85	3.99	80	4.39	88	3.63	73	3.93	79	12-158	0.401734	1.20	5.7
2-Chlorophenol-d4	2.62	52	3.02	60	3.15	63	3.62	72	3.63	73	3.62	72	2.71	54	3.07	61		0.366151	1.10	
2-Chlorophenol	2.82	56	3.21	64	3.35	67	4.21	84	4.04	81	4.1	82	3.4	68	3.82	76	23-134	0.407256	1.22	3.3
1,3-Dichlorobenzene	2.33	47	2.87	57	2.82	56	3.51	70	3.36	67	3.47	69	2.8	56	3.33	67	D-172	0.320565	0.96	1.9
1,4-Dichlorobenzene	2.5	50	2.84	57	2.88	58	3.51	70	3.53	71	3.56	71	2.92	58	3.5	70	20-124	0.346046	1.04	4.4
1,2-Dichlorobenzene	2.46	49	2.95	59	3.02	60	3.73	75	4	80	3.66	73	3.09	62	3.58	72	32-129	0.409134	1.23	1.9
Benzyl Alcohol	2.25	45	2.28	46	2.88	58	3.54	71	4.28	86	3.4	68	2.72	54	4.07	81		0.727255	2.18	
2,2'-oxybis(1-chloropropane)	2.1	42	2.33	47	1.89	38	3.99	80	3.6	72	1.89	38	3.55	71	3.84	77	36-166	0.936318	2.81	5.7
2-Methylphenol	2.75	55	3.21	64	3.06	61	4.12	82	4.62	92	3.85	77	3.51	70	3.97	79		0.544111	1.63	
4-Methylphenol-d8	2.34	47	2.83	57	2.97	59	3.68	74	3.85	77	3.47	69	2.56	51	3.04	61		0.474482	1.42	
3/4-Methylphenol	2.56	51	3.02	60	3.18	64	3.64	73	4.32	86	3.6	72	3.49	70	3.84	77		0.428325	1.28	
Hexachloroethane	2.21	44	2.56	51	2.78	56	3.42	68	3.53	71	3.05	61	2.74	55	3.33	67	40-113	0.377334	1.13	1.6
N-Nitroso-di-n-propylamine	3.47	69	3.75	75	4.03	81	4.74	95	5.14	103	4.69	94	3.66	73	4.13	83	D-230	0.55752	1.67	
Naphthalene-d8	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00	
Nitrobenzene	2.66	53	3.17	63	3.29	66	4.05	81	4.3	86	4.31	86	3.7	74	3.8	76	35-180	0.454522	1.36	1.9
Nitrobenzene-d5	2.37	47	3	60	2.92	58	3.58	72	3.9	78	3.83	77	2.73	55	3.15	63		0.465669	1.40	
2-Nitrophenol-d4	1.79	36	2.28	46	2.47	49	2.99	60	3.08	62	3.05	61	2.68	54	2.98	60		0.31686	0.95	
Isophorone	2.66	53	3.14	63	3.21	64	4.05	81	4.2	84	4.12	82	3.62	72	3.92	78	21-196	0.435677	1.31	
2-Nitrophenol	2.08	42	2.54	51	2.9	58	3.19	64	3.45	69	3.44	69	3.4	68	3.73	75	29-182	0.399535	1.20	2.6
2,4-Dimethylphenol	1.92	38	2.42	48	1.76	35	3.89	78	3.38	68	2.22	44	3.6	72	3.73	75	32-119	0.847722	2.54	2.7
Bis(2-chloroethoxy)methane	2.79	56	3.3	66	3.42	68	4.29	86	4.33	87	4.47	89	3.76	75	4.01	80	33-184	0.460579	1.38	5.3
2,4-Dichlorophenol-d3	1.9	38	2.2	44	2.65	53	3.06	61	3.14	63	3.21	64	2.31	46	2.53	51		0.410505	1.23	
2,4-Dichlorophenol	2.38	48	3.14	63	2.86	57	3.63	73	3.57	71	3.44	69	3.15	63	3.66	73	39-135	0.304193	0.91	2.7
1,2,4-Trichlorobenzene	2.45	49	2.95	59	2.95	59	3.65	73	3.5	70	3.68	74	3.09	62	3.46	69	44-142	0.320669	0.96	1.9
Naphthalene	2.81	56	3.23	65	3.23	65	3.96	79	3.89	78	3.91	78	3.5	70	3.8	76	21-133	0.321188	0.96	1.6

Entry Identifier	MDL #1		MDL #2		MDL #3		MDL #4		MDL #5		MDL #6		MDL #7		MDL #8		EPA Limit	STD		MDL	EPA
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec		DEV			
4-Chloroaniline-d4	1.83	37	2.07	41	2.14	43	2.59	52	2.86	57	2.76	55	2.14	43	2.5	50		0.321907	0.97		
4-Chloroaniline	1.94	39	2.25	45	2.46	49	3.17	63	2.93	59	3.2	64	2.97	59	2.88	58		0.355327	1.07		
Hexachlorobutadiene	2.04	41	2.52	50	2.67	53	3.31	66	3.22	64	3.38	68	2.85	57	3.24	65	24-116	0.342331	1.03	0.9	
Acenaphthene-d10	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00		
2-Methylnaphthalene	2.98	60	3.46	69	3.34	67	4.14	83	4.21	84	4.42	88	3.06	61	3.62	72		0.509542	1.53		
1-Methylnaphthalene	2.95	59	3.43	69	3.39	68	4.18	84	4.12	82	4.31	86	3.09	62	3.6	72		0.469589	1.41		
4-Chloro-3-methylphenol	2.53	51	2.69	54	2.88	58	3.33	67	3.3	66	3.74	75	2.85	57	3.29	66	22-147	0.364457	1.09	3	
Hexachlorocyclopentadiene	1.21	24	1.62	32	1.64	33	2.1	42	2.19	44	2.2	44	2.05	41	2.4	48		0.293509	0.88		
2,4,6-trichlorophenol	2.54	51	2.99	60	2.83	57	3.5	70	3.69	74	3.75	75	2.85	57	3.25	65	37-144	0.389095	1.17	2.7	
2,4,5-Trichlorophenol	2.26	45	2.5	50	2.61	52	3.42	68	3.4	68	3.27	65	2.41	48	3.24	65		0.449794	1.35		
2-Chloronaphthalene	2.77	55	3.24	65	3.14	63	3.92	78	3.85	77	3.79	76	3.01	60	3.41	68	60-118	0.370945	1.11	1.9	
2-Nitroaniline	2.07	41	1.88	38	2.28	46	2.65	53	2.77	55	2.44	49	2.52	50	2.85	57		0.329892	0.99		
1,4-Dinitrobenzene	1.8	36	1.93	39	2.05	41	2.65	53	2.77	55	2.64	53	2.89	58	2.76	55		0.37792	1.13		
Dimethylphthalate-d6	2.82	56	3.41	68	3.38	68	3.85	77	4.08	82	3.88	78	2.94	59	2.96	59		0.452953	1.36		
Dimethylphthalate	3.25	65	3.76	75	3.78	76	4.32	86	4.67	93	4.39	88	3.62	72	3.8	76	D-112	0.403585	1.21	1.6	
1,3-Dinitrobenzene	2.01	40	2.24	45	2.33	47	2.64	53	2.94	59	3.09	62	2.94	59	2.91	58		0.331246	0.99		
Acenaphthylene-d8	2.59	52	3.15	63	3.09	62	3.65	73	3.78	76	3.58	72	2.56	51	2.77	55		0.46328	1.39		
Acenaphthylene	3.09	62	3.55	71	3.44	69	4.28	86	4.09	82	4.29	86	3.15	63	3.6	72	33-145	0.447863	1.34	3.5	
2,6-Dinitrotoluene	2.87	57	3.34	67	3.4	68	4.01	80	4.1	82	4.12	82	3.65	73	3.57	71	50-158	0.331533	0.99	1.9	
1,2-Dinitrobenzene	2.6	52	2.81	56	3.3	66	3.69	74	4.05	81	3.73	75	3.38	68	3.33	67		0.396779	1.19		
Acenaphthene	3.11	62	3.55	71	3.41	68	4.06	81	4.01	80	4.17	83	3.18	64	3.67	73	47-145	0.369974	1.11	1.9	
3-Nitroaniline	2.21	44	2.24	45	2.09	42	2.33	47	2.29	46	2.31	46	2.47	49	1.84	37		0.203786	0.61		
2,4-Dinitrophenol	1.65	33	1.79	36	1.75	35	1.53	31	1.98	40	1.97	39	2.04	41	1.07	21	D-191	0.340525	1.02	24	
Dibenzofuran	2.91	58	3.46	69	3.28	66	3.89	78	3.92	78	4.06	81	3.11	62	3.55	71		0.356371	1.07		
2,4-Dinitrotoluene	2.29	46	2.95	59	3.03	61	3.51	70	3.69	74	3.45	69	3.53	71	3.24	65	39-139	0.276207	0.83	5.7	
4-Nitrophenol-d4	1.35	27	1.61	32	1.77	35	1.94	39	1.83	37	2.03	41	1.77	35	2.09	42		0.167204	0.50		
4-Nitrophenol	1.68	34	1.75	35	1.69	34	1.82	36	1.97	39	1.81	36	1.53	31	1.49	30	D-132	0.168989	0.51	2.4	
2,3,4,6-Tetrachlorophenol	2.86	57	3.24	65	3.24	65	3.62	72	4.04	81	3.88	78	3.44	69	3.47	69		0.306509	0.92		
2,3,5,6-Tetrachlorophenol	2.48	50	2.97	59	3.17	63	3.36	67	3.8	76	3.65	73	3.23	65	3.17	63		0.292908	0.88		
Fluorene-d10	0.11	2	0.14	3	0.13	3	0.15	3	0.15	3	0.15	3	2.79	56	3.06	61		1.359245	4.08		
Fluorene	3.05	61	3.65	73	3.45	69	4.12	82	4.3	86	4.32	86	3.29	66	3.6	72	59-121	0.421404	1.26	1.9	
4-Chlorophenyl phenyl ether	3.12	62	3.7	74	3.45	69	4.14	83	4.24	85	4.29	86	3.19	64	3.51	70	25-158	0.435332	1.31	4.2	
Diethylphthalate	3.38	68	3.79	76	3.7	74	4.3	86	4.72	94	4.34	87	3.79	76	3.6	72	D-114	0.41912	1.26	1.9	
4-Nitroaniline	1.58	32	1.43	29	n.d. < 1.42	28	1.5	30	1.63	33	1.97	39	1.93	39	1.98	40		0.24996	0.75		
Phenanthrene-d10	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00		

Entry Identifier	MDL #1		MDL #2		MDL #3		MDL #4		MDL #5		MDL #6		MDL #7		MDL #8		EPA Limit	STD DEV	MDL	EPA
	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec	ug/L	% Rec				
4,6-Dinitro-2-methylphenol-d2	2.25	45	2.41	48	n.d. < 2.52	50	2.6	52	2.76	55	2.88	58	1.48	30	1.52	30		0.62083	1.86	
4,6-Dinitro-2-methylphenol	2.52	50	2.59	52	2.73	55	2.73	55	3.23	65	3.22	64	2.98	60	2.77	55	D-181	0.254343	0.76	42
Diphenylamine	2.87	57	3.36	67	3.37	67	3.74	75	4.24	85	3.77	75	3.46	69	3.4	68		0.322852	0.97	1.9
Azobenzene	3.34	67	3.79	76	3.84	77	4.16	83	4.5	90	4.28	86	3.74	75	3.91	78		0.286498	0.86	
4-Bromophenyl phenyl ether	2.85	57	3.33	67	3.3	66	3.7	74	3.86	77	3.95	79	3.31	66	3.44	69	53-127	0.276578	0.83	
Hexachlorobenzene	3.08	62	3.3	66	3.24	65	3.6	72	3.87	77	3.66	73	3.38	68	3.21	64	D-152	0.248452	0.74	1.9
Pentachlorophenol	1.94	39	1.91	38	2.11	42	2.26	45	2.3	46	2.03	41	2.28	46	2.29	46	14-176	0.153994	0.46	3.6
Phenanthrene	3.83	77	4.02	80	3.98	80	4.31	86	4.65	93	4.49	90	4.02	80	3.99	80	54-120	0.275525	0.83	5.4
Anthracene-d10	2.99	60	3.45	69	3.44	69	3.61	72	3.93	79	3.76	75	3.41	68	3.15	63		0.255725	0.77	1.9
Anthracene	3.73	75	3.89	78	3.86	77	4.24	85	4.64	93	4.49	90	4	80	4.2	84	27-133	0.297233	0.89	
Carbazole	3.71	74	3.75	75	3.76	75	4.05	81	4.69	94	4.11	82	4.45	89	4.13	83		0.342533	1.03	
Di-n-butylphthalate	3.93	79	4.02	80	4.15	83	3.98	80	5.15	103	4.54	91	4.46	89	4.05	81	1-118	0.420827	1.26	2.5
Fluoranthene	4.01	80	3.95	79	4.17	83	4.04	81	4.95	99	4.49	90	4.78	96	4.27	85	26-137	0.377201	1.13	2.2
Pyrene	4.14	83	3.95	79	4.13	83	4.16	83	4.88	98	4.56	91	4.65	93	4.4	88	52-115	0.329343	0.99	1.9
Benzyl butyl phthalate	4.5	90	4.11	82	4.28	86	4.37	87	5.2	104	5	100	5.01	100	4.28	86	D-152	0.444512	1.33	2.5
Benzidine	0.29	6	0.08	2	0.07	1	0.07	1	0.09	2	0.06	1	1.42	28	0.86	17		0.544776	1.63	44
Chrysene-d12	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00	
Benzo[a]anthracene	3.83	77	3.82	76	3.79	76	3.81	76	4.43	89	4.23	85	4.26	85	3.73	75	33-143	0.285832	0.86	7.8
Chrysene	4	80	4.04	81	3.96	79	4.04	81	4.77	95	4.43	89	4.86	97	4.05	81	17-168	0.379549	1.14	2.5
Bis(2-ethylhexyl)adipate	3.37	67	3.47	69	3.65	73	3.52	70	4.18	84	3.81	76	4.08	82	3.41	68		0.303284	0.91	
3,3-Dichlorobenzidine	1.09	22	1	20	0.89	18	0.84	17	0.89	18	0.87	17	1.34	27	1.23	25	D-262	0.197773	0.59	16.5
Bis(2-ethylhexyl)phthalate	3.64	73	4.53	91	4.03	81	3.67	73	4.43	89	4.35	87	3.86	77	3.79	76	8-158	0.34146	1.02	2.5
Di-n-octyl phthalate	3.28	66	3.54	71	3.61	72	3.54	71	4.17	83	3.83	77	3.86	77	3.5	70	4-146	0.244774	0.73	2.5
Perylene-d12	20	100	20	100	20	100	20	100	20	100	20	100	20	100	20	100		0	0.00	
Benzo[a]pyrene-d12	3.26	65	4.09	82	3.54	71	3.13	63	3.63	73	3.44	69	0.79	16	0.62	12		1.425546	4.27	
Benzo[a]pyrene	3.67	73	3.45	69	3.68	74	3.53	71	4.11	82	3.92	78	4.54	91	3.7	74	17-163	0.379109	1.14	2.5
Benzo[b]Fluoranthene	3.8	76	3.64	73	3.73	75	3.73	75	4.2	84	4.19	84	4.69	94	3.59	72	24-159	0.405821	1.22	4.8
Benzo[k]fluoranthene	4.08	82	3.91	78	4.15	83	3.99	80	4.58	92	4.5	90	4.77	95	4.08	82	11-162	0.330742	0.99	2.5
Indeno[1,2,3-cd]pyrene	3.16	63	4.07	81	4.29	86	3.3	66	3.84	77	3.64	73	3.76	75	2.93	59	D-171	0.458984	1.38	3.7
Dibenz[a,h]anthracene	3.43	69	3.26	65	3.37	67	2.88	58	3.74	75	3.79	76	3.62	72	3.16	63	D-227	0.332802	1.00	2.5
Benzo[g,h,i]perylene	3.7	74	3.45	69	3.43	69	3.22	64	3.69	74	3.65	73	4.27	85	3.05	61	D-219	0.393731	1.18	4.1

Expanded Methodology




EPA Methods (examples)



- EPA 8082
- EPA 8081
- EPA 1694
- EPA 1664
- EPA 1613
- EPA 8151
- EPA 1614
- EPA 8310
- EPA 8290
- EPA 8330
- EPA 8141
- Any EPA 500 series SPE method

Examples

Extraction Cartridge



Clean-up Cartridge



Examples

Validated Protocols with inline drying/cleaning

PAH's with drying and/or Silica Gel clean-up

**Organochlorine Pesticides with drying
and/or Florisil**

Organophosphorus Pesticides with Drying



Typical mid level OCP Recoveries

	<u>LCS</u>	<u>LCS</u>	<u>LCS</u>	<u>LCS</u>	<u>Mean</u>	
<u>Analyte</u>	<u>#1</u>	<u>#2</u>	<u>#3</u>	<u>#4</u>	<u>Rec</u>	<u>RSD</u>
α -BHC	102.5%	106.8%	106.3%	107.9%	106%	2.2%
γ -BHC (Lindane)	102.1%	105.9%	105.5%	107.0%	105%	2.0%
β -BHC	96.3%	98.0%	100.9%	100.0%	99%	2.1%
δ -BHC	107.7%	110.6%	113.0%	111.8%	111%	2.1%
Heptachlor	86.1%	106.4%	107.2%	101.6%	100%	9.8%
Aldrin	68.7%	80.2%	79.6%	77.0%	76%	6.9%
Heptachlor Epoxide	99.5%	100.0%	102.1%	103.6%	101%	1.9%
γ -Chlordane	95.5%	88.7%	92.8%	91.3%	92%	3.1%
Endosulfan I	93.2%	93.4%	96.8%	96.2%	95%	2.0%
α -Chlordane	98.9%	97.2%	100.7%	101.8%	100%	2.0%
4,4'-DDE	91.8%	92.3%	97.1%	94.5%	94%	2.6%
Dieldrin	101.2%	99.3%	103.0%	104.3%	102%	2.1%
Endrin	108.2%	107.6%	112.9%	113.3%	111%	2.7%
4,4'-DDD	104.6%	97.2%	105.5%	104.3%	103%	3.7%
Endosulfan II	102.5%	97.7%	102.7%	104.1%	102%	2.7%
Endrin Aldehyde	85.4%	77.6%	89.4%	84.1%	84%	5.8%
4,4'-DDT	105.2%	97.5%	106.9%	105.9%	104%	4.1%
Endosulfan Sulfate	104.3%	99.1%	104.5%	106.2%	104%	3.0%
Methoxychlor	104.8%	90.9%	101.2%	103.4%	100%	6.3%
Endrin Ketone	102.7%	93.1%	101.6%	101.8%	100%	4.5%
<u>Surrogate</u>						
TCMX	66%	73%	75%	75%	72.3%	5.7%
PCB-209	88%	69%	74%	84%	78.7%	11.2%

Typical mid level OPP Recoveries

<u>Compound</u>	<u>Position #1</u>	<u>Position #2</u>	<u>Mean Rec</u>	<u>RSD</u>	<u>Compound</u>	<u>Position #1</u>	<u>Position #2</u>	<u>Mean Rec</u>	<u>RSD</u>
Alachlor	4.72	4.61	93.3%	2.4%	Mevinphos	4.57	4.38	89.5%	4.2%
Ametryn	4.46	4.5	89.6%	-0.9%	MGK 264	4.13	3.59	77.2%	14.0%
Atraton	4.59	4.26	88.5%	7.5%	Molinate	4.24	3.97	82.1%	6.6%
Atrazine	4.49	4.42	89.1%	1.6%	Napropamide	4.49	4.28	87.7%	4.8%
Bromacil	4.18	5.07	92.5%	-19.2%	Norflurazon	4.74	4.67	94.1%	1.5%
Butachlor	4.08	4.09	81.7%	-0.2%	Pebulate	4.07	3.52	75.9%	14.5%
Butylate	3.62	3.54	71.6%	2.2%	Prometon	4.16	4.3	84.6%	-3.3%
Carboxin	4.3	3.64	79.4%	16.6%	Prometryn	4.4	4.57	89.7%	-3.8%
Chlorophram	4.61	4.17	87.8%	10.0%	Pronamide**	4.46	4.36	88.2%	2.3%
Cycloate	4.11	3.9	80.1%	5.2%	Propazine	4.72	4.55	92.7%	3.7%
Dichlorvos	4.23	3.91	81.4%	7.9%	Simazine	4.2	4.24	84.4%	-0.9%
Diphenamid	4.13	4.55	86.8%	-9.7%	Simetryn	4.79	4.73	95.2%	1.3%
Disulfoton*	4.03	3.32	73.5%	19.3%	Stirofos	4.33	4.4	87.3%	-1.6%
EPTC	3.74	3.41	71.5%	9.2%	Tebuthiuron	4.36	4.18	85.4%	4.2%
Ethoprop	4.64	4.44	90.8%	4.4%	Terbacil	4.81	5.33	101.4%	-10.3%
Fenamiphos	4.21	5.07	92.8%	-18.5%	Terbufos**	3.58	3.51	70.9%	2.0%
Fenarimol	5.27	4.82	100.9%	8.9%	Terbutryn	4.43	4.82	92.5%	-8.4%
Fluridone	5.51	5.06	105.7%	8.5%	Triademefon	5.69	4.82	105.1%	16.6%
Hexazinone	4.86	4.69	95.5%	3.6%	Tricyclazole	5.25	4.95	102.0%	5.9%
Methyl Paraoxon	5.35	5.56	109.1%	-3.8%	Vernolate	3.69	3.529	72.2%	4.5%
Metolachlor	4.5	4.41	89.1%	2.0%					

Typical mid level PAH Recoveries

<u>Compound</u>	<u>LCS</u>	<u>LCSD</u>	<u>RPD</u>
Naphthalene	92.0%	86.0%	-6.7%
2-Methylnaphthalene	93.6%	86.8%	-7.5%
Acenaphthylene	97.8%	85.2%	-13.8%
Acenaphthene	96.3%	107.2%	10.7%
Fluorene	101.6%	109.6%	7.6%
Phenanthrene	97.6%	106.8%	9.0%
Anthracene	107.2%	109.6%	2.2%
Fluoranthene	108.0%	106.8%	-1.1%
Pyrene	110.0%	108.4%	-1.5%
Benzo[a]anthracene	110.0%	108.0%	-1.8%
Chrysene	108.8%	108.4%	-0.4%
Benzo[b]fluoranthene	109.2%	106.8%	-2.2%
Benzo[k]fluoranthene	108.8%	106.4%	-2.2%
Benzo[a]pyrene	110.0%	108.0%	-1.8%
Indeno[1,2,3-c,d]pyrene	98.0%	94.4%	-3.7%
Dibenzo[a,h]anthracene	101.6%	105.2%	3.5%
Benzo[g,h,i]perylene	109.6%	106.4%	-3.0%

Summation

- Fully Automated 625/8270 extractions
- Expandable to meet any SPE method
- Capable of performing in line extract drying and/or Cartridge extract clean-ups
- In-line evaporation with direct to GC vial tubes.



Questions?

