



# EPA 8270D Optimized for Modern Instrumentation

---

NEMC 2018

Dale Walker, Melissa Churley, Bruce Quimby,  
Mike Szelewski, Anastasia Andrianova

Agilent

## The goal of EPA 8270 re-optimization

8270 criteria were initially developed on specific instrument platforms

- rigorous due to diverse mix of target compound classes and wide calibration range

The goal of re-optimization was to achieve the **widest calibration range in just one injection**

- **using standard column dimensions and specified upper limit**



## EPA 8270 initial calibration considerations

Initial calibration results are indicative of how long the continuing calibration will last

- a lower %RSD during initial calibration is indicative that the calibration will last longer
- decreased cost of operation results from less frequent intervention by the operator

Splitless sample injection often results in overloaded column capacity and manual integration is often required

- at calibration levels approaching 160 ppm and with a 0.25  $\mu\text{m}$  column film thickness to avoid long run times

Many practitioners have adopted a “ramped flow” approach

- reduces the effective amount injected

## Batch review: initial calibration 0.05 – 160 ppm (10 levels)

Calibration points were deleted from each end of range to meet one of the **required EPA 8270D criteria**:

1. Avg RF %RSD  $\leq 20$  (usually most desirable)
2. If not, linear curve fit  $R^2 \geq 0.990$  (usually desired to be minimal)
3. If not, quadratic fit (meets method criteria)

6 points needed for a curve fit

Accuracy for lowest point needs to be within  $\pm 30\%$

## Systems used for 8270 study

All systems: 30m x 250  $\mu$ m x 0.25  $\mu$ m DB-UI 8270D column

7890 - Inert Plus 5977 MSD    3, 6 mm lens (3 Units)

7890 - Inert Plus 5977 MSD    9 mm lens (4 Units)    Also with 0.50  $\mu$ m film (2 Units)

Intuvo GC - Inert 5977 MSD    6 mm lens

7890 – Inert Plus 5977 MSD    9 mm lens, with online H<sub>2</sub>



## Conditions tested

Tunes: DFTPP, Atune

Extractor lens diameters : 3 (standard), 6, and 9 mm

Inlet parameters: Pulsed splitless, Pulsed split, “ramped flow” \*

Coated (deactivated) or uncoated SSL

\*technical splitless injection but with brief initial flow ramp and Purge ON, which effectively reduces the amount of loaded sample. Mimics about a 6:1 split.



## Conditions tested

GC Liners: Agilent dual taper no wool (5190-3315)

split liner, low pressure drop, LPD (5190-2295)

single taper Glass Wool UI (5190-2293)

Other manufacturer splitless ST Glass Wool



## Conditions tested

Constant Flow at 1.2 mL/min was used throughout (Retention Time Locked)

Temperatures held constant: Inlet, transfer line, source (300°C) and quadrupole (150°C)

Oven ramps were consistent





## A note about GC liners with Ultra Inert SSL inlets

UI Universal Low Pressure Drop (LPD) Liner 5190-2295



We found this LPD 2295 liner to work better, especially in UI coated inlets. This liner has the smaller OD of a split liner and is not so tight a fit.

# Compounds in calibration standards

## Phenols

Phenol, 2-fluoro-  
Phenol-d6-  
Phenol  
2-chlorophenol  
2 methylphenol  
p-Cresol  
2,4- Dimethylphenol  
2,4-dichlorophenol  
4-chloro-3-methylphenol  
2,4,5-trichlorophenol  
2,4,6-trichlorophenol  
2-methyl-4,6-dinitrophenol  
2,4,6-tribromophenol  
Pentachlorophenol

## Bases

benzidine  
Pyridine  
Aniline  
4-Chloroaniline  
Carbazole  
3,3'-Dichlorobenzidine

## Nitros

Nitrobenzene-D5  
nitrobenzene  
2-nitrophenol  
2- nitroaniline  
2,6-dinitrotoluene  
3-Nitroaniline  
2,4-dinitrophenol  
4-nitrophenol  
2,4-dinitrotoluene  
4-Nitroaniline

## Acids

Benzoic acid

## Ethers

Bis(2-chloroethyl) ether  
Bis(2-chloro-1-methylethyl) ether  
bis(2-chloroethoxy)-methane  
4 chlorophenylphenyl ether  
4 bromophenylphenyl ether  
Dibenzofuran

## Nitrosos

N-Nitrosodimethylamine  
N-Nitrosodipropylamine  
N-Nitrosodiphenylamine

## Phthalates

Dimethyl phthalate  
Diethyl Phthalate  
Dibutyl phthalate  
Benzyl butyl phthalate  
Bis(2-ethylhexyl) phthalate  
Di-n-octyl phthalate

## ISTDs

1,4-Dichlorobenzene-D4  
Naphthalene-D8  
Acenaphthene-D10  
Phenanthrene-D10  
Chrysene-D12  
Perylene-D12


## PAHs

Naphthalene  
2-methylnaphthalene  
2-fluorobiphenyl  
2-chloronaphthalene  
Acenaphthylene  
Acenaphthene  
Fluorene  
Phenanthrene  
Anthracene  
Fluoranthene  
Pyrene  
p-Terphenyl-d14  
Benz[a]anthracene  
Chrysene  
Benzo[b]fluoranthene  
Benzo[k]fluoranthene  
Benzo[a]pyrene  
Indeno[1,2,3-cd]pyrene  
Dibenz[a,h]anthracene  
Benzo(g,h,i)perylene

## Misc

1,3-dichlorobenzene  
1,4-dichlorobenzene  
Benzyl alcohol  
1,2-dichlorobenzene  
Hexachloroethane  
Isophorone  
1,2,4-trichlorobenzene  
Hexachlorobutadiene  
Hexachlorocyclopentadien  
Azobenzene  
Hexachlorobenzene

## 77 Compounds and 6 ISTDs



TIC: 50ppm\_NewMix\_0.D\data.ms




## Instrument optimization – fit the range to a specific system

1. Inject the highest calibration standard, at a low gain setting (e.g., 1), adjust the split ratio (or ramped flow)
  - check performance of low responders
  - **look for overloaded peak shapes**
  - ensure that the separation of benzo[b] and benzo[k] fluoranthene isomers meets requirements (>50%)
2. **Properly set the gain factor** so that the BPC tallest peak is  $3-5 \times 10^6$  counts (SQ) or  $10^7$  (TQ)
  - Avoid detector saturation




## Instrument optimization – fit the range to a specific system



TIC: 50ppm\_NewMix\_0.D\data.ms

# Why is gain factor important? The detector ages over time



## Analysis of individual batches

Data points for each compound were deleted in some cases to meet method criteria (lower, upper end of range)

Results were entered into a spreadsheet columns (for each batch)

- average %RSD for the entire 77-compound list
- number of compounds that met average RF %RSD criteria
- those that required a linear curve fit



# Initial calibration batch review showing tracking of exceptions

0.05 – 160 ppm (10 points)

Avg RF %RSD: green if passes ( $\leq 20\%$ )

Blue if linear fit

Pulsed split 3:1 Inst 1, LPD liner, 9mm Analyst 1 Atune					Exceptions			
All curve fits are LINEAR unless there is an Exception								
Calibration is 0.05 - 160 ppm, 10 Levels = 0.05, 0.2, 0.5, 2, 5, 20, 50, 80, 120, 160 (4 ppm IS constant)								
Name	RT	Avg. RF RSD	CF R2	Cal Range	# points	Curve fit	pts deleted	
N-Nitrosodimethylamine	3.010	25.4	1.000					
Pyridine	3.052	9.8	1.000					
Phenol, 2-fluoro-	4.764	13.6	1.000					
Phenol-d6-	6.347	14.2	1.000					
Phenol	6.374	11.7	0.999					
Aniline	6.411	7.7	0.993					
Bis(2-chloroethyl) ether	6.561	2.8	1.000					
2-chlorophenol	6.615	12.4	1.000					

Exceptions – range, # points deleted, etc.



## Summary of all batches spreadsheet - sorted for best results

Excerpt, ranked by lowest average %RSD:

System	Inlet Coat.	Lens diam.	Injection Type	Liner	Tune Type	Pass DFTPP	Chemist	77 Cmpds	# Cmpds	# Cmpds	# Cmpds	CAL Range	# Cmpds Deleted	# Points Deleted	Presumptive		Compounds that are
								Avg %RSD	< 20% RSD	Linear	Quadr				Cal Range		
System 2	C	9	RF	LPD	A	Y	a	10.2	73	3	1	0.05-160	14	24	0.2	160	7
System 2	C	9	3:1 PSP	LPD	A	Y	a	11.2	69	7	1	0.05-160	9	13	0.2	160	4
System 3	C	9	RF	LPD	D	Y	b	11.3	70	6	1	0.05 - 160	34	51	0.5	160	6
System 2	C	6	RF	Split	A	Y	a	11.5	66	10	1	0.05 - 160	18	31	0.5	160	5

- number of **compounds** that required point deletion to meet acceptability criteria
- total number of **data points deleted** for all compounds (to meet criteria)
- presumptive working calibration range for each compound
  - useful range when the number of exceptions was < 8 of 77 compounds, or 10% of the target list
- number of compound exceptions to that range (acceptable, but outside the working range)

## Improved sensitivity for low responding compounds with 9mm lens

MDL ratios were calculated for easy comparison (3 vs. 6mm and 3 vs. 9 mm)

On the table on the next slide:

Higher ratios (>1) show that the MDL is lower (better sensitivity) using a larger diameter lens

- MDL Criterion: “MDL” is the amount (pg) required to achieve  $S/N \geq 5$ 
  - MH Quant – noise regions manually selected
  - peak-to-peak



# Improved sensitivity for low responding compounds with 9mm lens

	r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm		r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm
Acenaphthene	10.904	0.2	0.3	Phenol-d6-	6.337	0.8	0.8
Diethyl Phthalate	11.327	0.3	0.7	Fluorene	11.439	0.8	0.5
Phenol, 2-fluoro-	4.764	0.4	0.2	2,4,6-tribromophenol	11.685	0.9	1.0
4 chlorophenylphenyl ether	11.445	0.4	0.2	2-methylnaphthalene	9.733	0.9	0.7
bis(2-chloroethoxy)-methane	8.578	0.4	0.3	Benzo[a]anthracene	15.777	1.0	0.7
Di-n-octyl phthalate	17.344	0.4	0.6	Benzo[b]fluoranthene	17.975	1.0	0.9
Hexachlorobutadiene	9.054	0.5	0.5	1,2-dichlorobenzene	7.267	1.0	0.4
Chrysene	15.841	0.5	0.4	4 bromophenylphenyl ether	11.942	1.0	0.8
Carbazole	12.621	0.6	0.7	N-Nitrosodipropylamine	7.668	1.0	0.7
2-fluorobiphenyl	10.156	0.6	0.5	1,3-dichlorobenzene	6.893	1.0	0.4
Dimethyl phthalate	10.594	0.6	0.9	1,2,4-trichlorobenzene	8.786	1.0	0.3
Phenanthrene	12.413	0.6	0.3	Benzo[k]fluoranthene	18.029	1.0	1.0
Naphthalene	8.883	0.6	0.8	Bis(2-ethylhexyl) phthalate	15.9	1.1	1.3
Dibenzofuran	11.086	0.6	0.4	1,4-dichlorobenzene	7.027	1.1	0.4
2-chloronaphthalene	10.279	0.7	0.9	Benzo[g,h,i]perylene	20.923	1.1	0.5
Fluoranthene	13.648	0.7	1.0	2,4-dinitrophenol	10.937	1.1	1.5
Hexachlorobenzene	11.99	0.7	0.3	Pyrene	13.926	1.2	1.6
Hexachloroethane	7.792	0.7	0.5	Benzyl butyl phthalate	14.841	1.2	1.2
4-Chloroaniline	8.968	0.7	0.5	3-Nitroaniline	10.83	1.2	1.0
Bis(2-chloro-1-methylethyl) ether	7.471	0.8	0.3	2,4-dichlorophenol	8.674	1.2	0.6

- MDL ratios, small vs. large diameter lens: MDL is defined as the number of pg required to achieve S/N  $\geq 5$  (p2p)

- Ratios > 1 show that the MDL is lower (better sensitivity) using a larger diameter lens

- Many of the highest responders are seen in the top left of the table:

	r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm		r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm
p-Terphenyl -d14	14.124	1.2	0.7	Pentachlorophenol	12.188	1.9	2.6
Acenaphthylene	10.723	1.3	0.8	2,4,6-trichlorophenol	10.086	1.9	1.3
Anthracene	12.461	1.3	0.9	2,4,5-trichlorophenol	10.054	2.0	1.5
Pyridine	3.069	1.3	0.6	Benzoic acid	8.519	2.1	7.2
3,3'-Dichlorobenzidine	15.745	1.3	1.3	2-nitrophenol,	8.353	2.3	2.5
Hexachlorocyclopentadiene	9.915	1.3	0.7	2-methyl-4,6-dinitrophenol	11.487	2.3	5.3
4-chloro-3-methylphenol	9.567	1.4	1.7	4-Nitroaniline	11.455	2.3	2.1
Benzyl alcohol	7.235	1.4	0.3	2-nitroaniline	10.391	2.4	0.8
2-chlorophenol	6.615	1.4	1.3	Nitrobenzene -D5	7.872	2.4	0.8
Isophorone	8.251	1.5	1.1	Bis(2-chloroethyl) ether	6.556	2.5	1.2
p-Cresol	7.663	1.5	0.7	Aniline	6.412	2.9	1.0
N-Nitrosodimethylamine	3.015	1.6	1.0	Benzo[a]pyrene	18.58	3.0	1.3
N-Nitrosodiphenylamine	11.562	1.7	1.2	benzidine	13.809	3.1	4.9
2,4-dinitrotoluene	11.07	1.7	2.2	nitrobenzene	7.898	3.3	2.4
Phenol	6.363	1.7	1.0	Dibenz[a,h]anthracene	20.511	3.3	3.6
Azobenzene	11.605	1.8	1.9	2,6-dinitrotoluene	10.658	3.8	2.1
Dibutyl phthalate	12.964	1.8	1.2	2 methylphenol	7.417	3.9	0.7
Indeno[1,2,3-cd]pyrene	20.457	1.8	1.5	4-nitrophenol	10.995	4.4	2.4
2,4-Dimethylphenol	8.439	1.9	0.8				

- These are compounds for which sensitivity is much less of an issue and many have a ratio < 1

- Many challenging compounds (yellow) improved the most with a 9mm diameter lens (circled in red)

Low responding compounds

# Improved sensitivity for low responding compounds with 9mm lens

	r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm		r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm
Acenaphthene	10.904	0.2	0.3	Phenol-d6-	6.337	0.8	0.8
Diethyl Phthalate	11.327	0.3	0.7	Fluorene	11.439	0.8	0.5
Phenol, 2-fluoro-	4.764	0.4	0.2	2,4,6-tribromophenol	11.685	0.9	1.0
4 chlorophenylphenyl ether	11.445	0.4	0.2	2-methylnaphthalene	9.733	0.9	0.7
bis(2-chloroethoxy)-methane	8.578	0.4	0.3	Benzo[a]anthracene	15.777	1.0	0.7
Di-n-octyl phthalate	17.344	0.4	0.6	Benzo[b]fluoranthene	17.975	1.0	0.9
Hexachlorobutadiene	9.054	0.5	0.5	1,2-dichlorobenzene	7.267	1.0	0.4
Chrysene	15.841	0.5	0.4	4 bromophenylphenyl ether	11.942	1.0	0.8
Carbazole	12.621	0.6	0.7	N-Nitrosodipropylamine	7.668	1.0	0.7
2-fluorobiphenyl	10.156	0.6	0.5	1,3-dichlorobenzene	6.893	1.0	0.4
Dimethyl phthalate	10.594	0.6	0.9	1,2,4-trichlorobenzene	8.786	1.0	0.3
Phenanthrene	12.413	0.6	0.3	Benzo[k]fluoranthene	18.029	1.0	1.0
Naphthalene	8.883	0.6	0.8	Bis(2-ethylhexyl) phthalate	15.9	1.1	1.3
Dibenzofuran	11.086	0.6	0.4	1,4-dichlorobenzene	7.027	1.1	0.4
2-chloronaphthalene	10.279	0.7	0.9	Benzo(g,h,i)perylene	20.923	1.1	0.5
Fluoranthene	13.648	0.7	1.0	2,4-dinitrophenol	10.937	1.1	1.5
Hexachlorobenzene	11.99	0.7	0.3	Pyrene	13.926	1.2	1.6
Hexachloroethane	7.792	0.7	0.5	Benzyl butyl phthalate	14.841	1.2	1.2
4-Chloroaniline	8.968	0.7	0.5	3-Nitroaniline	10.83	1.2	1.0
Bis(2-chloro-1-methylethyl) ether	7.471	0.8	0.3	2,4-dichlorophenol	8.674	1.2	0.6

	r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm		r.t.	MDL ratio 3mm/6mm	MDL ratio 3mm/9mm
p-Terphenyl -d14	14.124	1.2	0.7	Pentachlorophenol	12.188	1.9	2.6
Acenaphthylene	10.723	1.3	0.8	2,4,6-trichlorophenol	10.086	1.9	1.3
Anthracene	12.461	1.3	0.9	2,4,5-trichlorophenol	10.054	2.0	1.5
Pyridine	3.069	1.3	0.6	Benzoic acid	8.519	2.1	7.2
3,3'-Dichlorobenzidine	15.745	1.3	1.3	2-nitrophenol,	8.353	2.3	2.5
Hexachlorocyclopentadiene	9.915	1.3	0.7	2-methyl-4,6-dinitrophenol	11.487	2.3	5.3
4-chloro-3-methylphenol	9.567	1.4	1.7	4-Nitroaniline	11.455	2.3	2.1
Benzyl alcohol	7.235	1.4	0.3	2-nitroaniline	10.391	2.4	0.8
2-chlorophenol	6.615	1.4	1.3	Nitrobenzene -D5	7.872	2.4	0.8
Isophorone	8.251	1.5	1.1	Bis(2-chloroethyl) ether	6.556	2.5	1.2
p-Cresol	7.663	1.5	0.7	Aniline	6.412	2.9	1.0
N-Nitrosodimethylamine	3.015	1.6	1.0	Benzo[a]pyrene	18.58	3.0	1.3
N-Nitrosodiphenylamine	11.562	1.7	1.2	benzidine	13.809	3.1	4.9
2,4-dinitrotoluene	11.07	1.7	2.2	nitrobenzene	7.898	3.3	2.4
Phenol	6.363	1.7	1.0	Dibenz[a,h]anthracene	20.511	3.3	3.6
Azobenzene	11.605	1.8	1.9	2,6-dinitrotoluene	10.658	3.8	2.1
Dibutyl phthalate	12.964	1.8	1.2	2 methylphenol	7.417	3.9	0.7
Indeno[1,2,3-cd]pyrene	20.457	1.8	1.5	4-nitrophenol	10.995	4.4	2.4
2,4-Dimethylphenol	8.439	1.9	0.8				

- MDL ratios, small vs. large diameter lens: MDL is defined as the number of pg required to achieve S/N  $\geq 5$  (p2p)

Sensitivity is better or equal using a 6 mm vs. a 3 mm lens (53 of 77 compounds)

- Ratios > 1 show that the MDL is lower (better sensitivity) using a larger diameter lens

- Many of the highest responders are seen in the top left of the table:


9 mm lens vs. 3 mm: 34 of 77 compounds show better or equal sensitivity

- These are compounds for which sensitivity is much less of an issue and many have a ratio < 1

- Many challenging compounds (yellow) improved the most with a 9mm diameter lens (circled in red)

Low responding compounds


## Sufficient peak resolution for benzo [b and k] fluoranthene using a 9mm lens



50 ppm, 1:3 split, LPD liner, 9 mm lens (17 ng injected)


From 8270D: Sufficient resolution is achieved if the height of the valley between two isomer peaks is less than 50% of the average of the two peak heights (at the mid-point concentration level)

## Calibration curve for benzoic acid using a 9mm lens




Calibration results for the challenging compound benzoic acid (0.05 – 120 ppm)

## Wide calibration range in one injection




## Wide calibration range in one injection





## Maximum number of compounds passing criteria and lowest %RSD

(conditions, in order of significance)

- ✓ ion source temperature  $\geq 300^{\circ}\text{C}$
- ✓ 9mm ion source drawout lens diameter
- ✓ PSP, pulsed split injection - best peak shape, widest calibration range
- ✓ Ultra Inert LPD split liner with wool  Also, the single taper Ultra Inert with glass wool
- ✓ either DFTPP tune or Autotune



## Summary of recommended conditions

<b>Injection type</b>	<b>Pulsed split (PSP)</b>	<b>Adjust to avoid peak overload Meet peak resolution criteria</b>
<b>GC liner</b>	<b>Low pressure drop split (LPD)</b>	<b>Best peak shape</b>
<b>column dimensions</b>	<b>30m x 0.25 mm x 0.25 <math>\mu</math>m DB - UI 8270D</b>	<b>Optimize peak separation</b>
<b>ion source draw out lens diameter</b>	<b>9 mm</b>	<b>Wide calibration range Improves sensitivity for active compounds</b>



---

Thank you very much