

Direct Injection Method for LC-MS/MS Detection of Polyfluoroalkyl Substances in Water

Douglas Stevens


Kari Organtini
Waters Corporation



Overview

- **PFAS Overview**
- **Drinking water analysis**
 - SPE sample prep
- **Other water type analyses**
 - Large volume direct injection (LVI/DI)

- 
- A polar bear is shown standing in a snowy, open landscape. The bear is facing the camera, with its head slightly turned to the left. It has thick, white fur and dark eyes. The ground is covered in snow with some dry grass and small plants visible. The background is a vast, flat, snowy expanse under a pale sky.



PFOS

METHOD 537. DETERMINATION OF SELECTED PERFLUORINATED ALKYL ACIDS IN DRINKING WATER BY SOLID PHASE EXTRACTION AND LIQUID CHROMATOGRAPHY/TANDEM MASS SPECTROMETRY (LC/MS/MS)

Version 1.1
September 2009

UCMR3 - PFASs

Third Unregulated Contaminant Monitoring Rule

Six Perfluorinated Compounds

Contaminant	CAS Registry Number ¹	Minimum Reporting Level	Sampling Points ²	Analytical Methods
perfluorooctanesulfonic acid (PFOS)	1763-23-1	0.04 µg/L 40 ppt	EPTDS	EPA 537 Rev 1.1
perfluorooctanoic acid (PFOA)	335-67-1	0.02 µg/L 20 ppt	EPTDS	EPA 537 Rev 1.1
perfluorononanoic acid (PFNA)	375-95-1	0.02 µg/L 20 ppt	EPTDS	EPA 537 Rev 1.1
perfluorohexanesulfonic acid (PFHxS)	355-46-4	0.03 µg/L 30 ppt	EPTDS	EPA 537 Rev 1.1
perfluoroheptanoic acid (PFHpA)	375-85-9	0.01 µg/L 10 ppt	EPTDS	EPA 537 Rev 1.1
perfluorobutanesulfonic acid (PFBS)	375-73-5	0.09 µg/L 90 ppt	EPTDS	EPA 537 Rev 1.1

Unregulated Contaminant Monitoring Rule

https://www.epa.gov/sites/production/files/2016-02/documents/ucmr4_stakeholderwebinarslides_160113_508.pdf



UCMR History

- UCMR 1 (2001-2005, 26 contaminants)
- UCMR 2 (2007-2011, 25 contaminants)
- UCMR 3 (2012-2016, 30 contaminants)
 - Monitoring concluded in 2015
 - Data review will occur in 2016
- UCMR 4 (2017-2021, 30 contaminants)
 - Proposed in the FR on December 11, 2015
 - Final publication anticipated in late 2016/early 2017

National occurrence data publically available:
<http://www.epa.gov/dwucmr/occurrence-data-unregulated-contaminant-monitoring-rule>


January 2016 U.S. Environmental Protection Agency 12

PFAS's in the Literature

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Hu, Xindi C., et al. "Detection of Poly-and Perfluoroalkyl Substances (PFASs) in US Drinking Water Linked to Industrial Sites, Military Fire Training Areas, and Wastewater Treatment Plants." *Environmental Science & Technology Letters* (2016)



PFAS Analysis Sensitivity Requirements

SPE enrichment

- Sample prep allows for use of low- to mid-level sensitivity for mass spectrometer



Large Volume Injection

- High sensitivity TQ required



Increasing Sensitivity

Comparison of analysis methods

	EPA Method 537	ISO 21675	ASTM
Sample Prep	SPE	SPE	Dilute and Filter aka No Prep
Injection Volume	10 µL	10 µL	30 µL
Number of compounds	14	27	38

UPLC/TQD Conditions


- System: ACQUITY UPLC/TQD with PFC Kit
- Analytical column: ACQUITY UPLC BEH C18 2.1x 50 mm
- Mobile Phase A: 2mM ammonium acetate in water/MeOH [95:5]
- Mobile Phase B: MeOH
- Column Temp: 50 °C
- Injection volume: 10 µL (full loop)

Time (min)	Flow (mL/min)	%B	Curve
0.0	0.40	25	
0.5	0.40	25	6
5.0	0.40	85	6
5.1	0.40	100	6
5.6	0.40	100	6
7.0	0.55	100	1
9.0	0.40	25	1



537

PFAS's Interferences/Contamination



PFC Kit

PEEK tubing
for solvent lines

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®



Stainless steel tubing
Isolator column to injector

PFC Isolator Column


Solvent mixer

Coiled stainless steel tube
Mixer to
PFC isolator column

PEEK tubing
for autosampler
lines




MRM Chromatograms of PFOA With PFC Kit Installed



■ Generates all the necessary methods for routine analysis

- Methods -> Sample List -> Data Processing



Method files and a sample batch list are required to run this analysis.
Enter the name to use for all the files created, the project in which they will be created and specify which files to create.

Name:

Project:

Files to create:


- ☒ MS method
- ☒ Chromatography method
- ☒ TargetLynx method
- ☒ Sample list

#	Name	Sample Text	RT	Area	Response	Reverse Fit	Forward Fit
1	0732013_46	Matrix match spiked 10 ppb (eq to sample) 39.55 ppb (in vial)	8.11	39517.477	147.774	873	726
2	0732013_50	Matrix match spiked 20 ppb (eq to sample) 79.09 ppb (in vial)	8.11	71666.023	251.299	806	619
3	0732013_54	Matrix match spiked 25 ppb (eq to sample) 98.87 ppb (in vial)	8.10	85963.734	359.276	937	743
4	0732013_58	Matrix match spiked 50 ppb (eq to sample) 197.74 ppb (in vial)	8.08	112297.180	507.213	942	749

Typical Calibration Curve

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

PFOA calibration curve (0.09 ppb to 9.4 ppb) $r^2 > 0.997$




537

Recovery Study

Target 10 PFCs	1 ppt spiked water (n=6)	
	Recovery (%)	RSD (%)
PFBuS	105	6
PFHxA	112	11
PFHpA	124	15
PFHxS	103	5
PFOA	122	17
PFNA	107	9
PFOS	104	6
PFDA	109	8
PFUnA	100	9
PFDoA	105	12

537

PFOS in Bottled Water



Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

PFOS in the bottled water

Injection	sample 1	Sample 2	sample 3	sample 4
1	0.42	0.39	0.44	0.29
2	0.37	0.44	0.43	0.35
3	0.48	0.52	0.39	0.35
4	0.48	0.46	0.39	0.38
5	0.44	0.52	0.42	0.35
6	0.36	0.43	0.36	0.36
Average (ppt)	0.42	0.46	0.40	0.35
RSD (%)	12.7	11.1	7.7	8.1
Mean = 0.41 ppt				537

Instrument LOQ <25ppt; **SPE enrichment factor =250; Method LOQ <0.1ppt**

Direct (large volume) Injection Approach ASTM 7979-17

ASTM 7979-17 Direct Injection Method Sample Prep

1. 5 mL MeOH + 5 mL sample in 15mL tube, vortex
2. Filter: Glass fiber + GHP (polypropylene)
3. Add 10 ul acetic acid
4. 500 ul to LC vial

- Reagent Water
- Surface (River) Water
- Ground Water
- Influent Waste Water
- Effluent Waste Water

June 2018 Application Note: **720006329EN**

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®



Instrument Methods

Source Parameters

- Instrument: **Xevo TQ-XS**
- Ion Mode: ESI-
- Capillary Voltage: 1.0 kV
- Desolvation Temperature: 500°C
- Desolvation Flow: 1100 L/hr
- Cone Flow: 150 L/hr

MS Method

- Developed using **QuanOptimize**
 - MRMs, CV, CE
- Divert flow to waste from 15 – 21 mins

LC Method

- Acquity I-Class with PFC Kit
- CSH Phenyl Hexyl 2.1mm x 100 mm, 1.7 μ m
- A: 95:5 H₂O:MeOH + 2 mM am ac
- B: MeOH + 2 mM am ac
- **Inj Vol: 30 μ l**
- Gradient:

Time (min)	Flow (mL/min)	%A	%B
0	0.3	100	0
1	0.3	80	20
6	0.3	55	45
13	0.3	20	80
14	0.4	5	95
17	0.4	5	95
18	0.3	100	0
22	0.3	100	0

7979-17

Compounds Included in Method



- 39 Native Compounds
- 27 Isotope Labels – Recovery only!

Abbreviation	Molecular Formula	CAS Number	PFAS Class
PFBA	C ₄ HF ₇ O ₂	375-22-4	carboxylate
PFPeA	C ₅ HF ₉ O ₂	2706-90-3	carboxylate
PFHxA	C ₆ HF ₁₁ O ₂	307-24-4	carboxylate
PFHpA	C ₇ HF ₁₃ O ₂	375-85-9	carboxylate
PFOA	C ₈ HF ₁₅ O ₂	335-67-1	carboxylate
PFNA	C ₉ HF ₁₇ O ₂	375-95-1	carboxylate
PFDA	C ₁₀ HF ₁₉ O ₂	335-76-2	carboxylate
PFUnDA	C ₁₁ HF ₂₁ O ₂	2058-94-8	carboxylate
PFDoDA	C ₁₂ HF ₂₃ O ₂	307-55-1	carboxylate
PFTriDA	C ₁₃ HF ₂₅ O ₂	72629-94-8	carboxylate
PFTreDA	C ₁₄ HF ₂₇ O ₂	376-06-7	carboxylate
PFHxDA	C ₁₆ HF ₃₁ O ₂	67905-19-5	carboxylate
PFOcDA	C ₁₈ HF ₃₅ O ₂	16517-11-6	carboxylate
PFBS	C ₄ HF ₉ SO ₃	29420-49-3	sulfonate
PFPeS	C ₅ HF ₁₁ SO ₃	2706-91-4	sulfonate
PFHxS	C ₆ HF ₁₃ SO ₃	3871-99-6	sulfonate
PFHpS	C ₇ HF ₁₅ SO ₃	375-92-8	sulfonate
PFOS	C ₈ HF ₁₇ SO ₃	1763-23-1	sulfonate

PFNS	C ₉ HF ₁₉ SO ₃		sulfonate
PFDS	C ₁₀ HF ₂₁ SO ₃	335-77-3	sulfonate
FOSA	C ₈ H ₂ F ₁₇ NO ₂ S	754-91-6	sulfonamide
N-EtFOSAA	C ₁₂ H ₈ F ₁₇ NO ₄ S	2991-50-6	sulfonamidoacetic acid
N-MeFOSAA	C ₁₁ H ₆ F ₁₇ NO ₄ S	2355-31-9	sulfonamidoacetic acid
N-MeFOSA	C ₉ H ₄ F ₁₇ NO ₂ S	31506-32-8	sulfonamide
N-EtFOSA	C ₁₀ H ₆ F ₁₇ NO ₂ S	4151-50-2	sulfonamide
FHUEA	C ₈ H ₂ F ₁₂ O ₂	70887-88-6	unsaturated telomer acid
FOUEA	C ₁₀ H ₂ F ₁₆ O ₂	70887-84-2	unsaturated telomer acid
8:2 diPAP	C ₂₀ H ₈ F ₃₄ PO ₄	678-41-1	phosphate ester
4:2 FTS	C ₈ H ₂ F ₉ SO ₃	757124-72-4	telomer sulfonate
6:2 FTS	C ₈ H ₂ F ₁₃ SO ₃	29420-49-3	telomer sulfonate
8:2 FTS	C ₁₀ H ₂ F ₁₇ SO ₃	39108-34-4	telomer sulfonate
PFecHS	C ₈ HF ₁₅ SO ₃	67584-42-3	cyclic
FHEA	C ₈ H ₂ F ₁₃ O ₂	53826-12-3	telomer acid
FOEA	C ₁₀ H ₂ F ₁₇ O ₂	27854-31-5	telomer acid
FDEA	C ₁₂ H ₂ F ₂₁ O ₂	N/A	telomer acid
FHpPA	C ₁₀ H ₂ F ₁₉ O ₂	812-70-4	other
ADONA	C ₇ H ₂ F ₁₂ O ₄	958445-44-8	other
9CI-PF3ONS	C ₈ HF ₁₆ ClSO ₄	73606-19-6	other
11CI-PF3OUDS	C ₁₀ HF ₂₀ ClSO ₄	73606-19-6	other

7979-17

Overall Method Summary




7979-17

Spiked Surface Water

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Sulfonates




7979-17

PFBA and PFPeA at 300 ng/L; 4:2, 6:2, and 8:2 FTS at 1200 ng/L; all other compounds at 60 ng/L.
(*) this compound shown off scale.

Data Independent Acquisition (DIA) Concurrent MRM + Full Scan (aka RADAR)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®



Method Detection Limits (MDL)

Below reporting limits described by ASTM 7979





Compound	Sample spike (ng/L)	MDL (ng/L)	Reporting range (ng/L)*	R ²
PFBA	100	25.20	50–2000	0.993
PFPeA	10	1.04	50–2000	0.999
PFHxA	10	1.33	10–400	0.999
PFHpA	10	0.91	10–400	0.999
PFOA	10	1.42	10–400	0.999
PFNA	10	1.22	10–400	0.999
PFDA	10	0.84	10–400	0.998
PFUnDA	10	2.52	10–400	0.996
PFDoDA	10	1.76	10–400	0.993
PFTriDA	10	2.34	10–400	0.991
PFTreDA	10	1.99	10–400	0.993
PFHxDA	200	25.41	–	0.984
PFOcDA	400	41.99	–	0.983
PFBS	10	1.21	10–400	0.999
PFPeS	10	1.07	10–400	0.999
PFHxS	10	1.41	10–400	0.999
PFHpS	10	1.57	10–400	0.999
PFOS	10	1.61	10–400	0.999
PFNS	10	1.67	10–400	0.999
PFDS	10	1.44	10–400	0.997

Compound	Sample spike (ng/L)	MDL (ng/L)	Reporting range (ng/L)*	R ²
FOSA	10	1.29	10–400	0.999
N-Et-FOSAA	10	1.90	10–400	0.997
N-Me-FOSAA	10	1.59	10–400	0.999
N-Et-FOSA	10	1.45	–	0.997
N-Me-FOSA	10	1.19	–	0.999
FHUEA	10	1.53	10–400	0.999
FOUEA	10	1.36	–	0.999
8:2 diPAP	300	50.16	–	0.988
4:2 FTS	10	1.50	10–400	0.999
6:2 FTS	10	N/A	10–400	0.999
8:2 FTS	10	2.62	10–400	0.997
PFecHS	10	1.17	10–400	0.998
FHEA	200	42.19	300–8000	0.994
FOEA	200	50.38	200–8000	0.997
FDEA	200	79.48	200–8000	0.993
FHpPA	10	1.47	10–400	0.999
ADONA	10	0.82	–	0.999
9CI-PF3ONS	10	1.06	–	0.999
11CI-PF3OUdS	10	1.45	–	0.998



7979-17

Linearity and Sensitivity

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®



**2.5 ng/L injection
(0.5x Reporting
Limit)**



7979-17

Recovery and Repeatability


Matrix	% Recovery PFOS
Reagent Water	92.8
Ground Water	90.6
Surface Water	94.6
Influent Water	94.1
Effluent Water	93.5

Overall Recoveries: 80 – 125%

*PFTreDA, PFTriDA ~ 200% - matrix effects

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

RSD: 5.5625



Overall RSD: < 10 %

7979-17

Water Sample Types and Recoveries

- Reagent Water
- Surface (River) Water
- Ground Water
- Influent Waste Water
- Effluent Waste Water

7979-17

	Average Percent Recovery in Matrix					Recovery Compound
	Reagent Water	Ground Water	Surface Water	Influent Water	Effluent Water	
PFBA	82.7	82.1	80.9	80.8	85.5	¹³ C-PFBA
PFPeA	89.1	87.7	90.2	88.1	91.4	¹³ C ₅ -PFPeA
PFHxA	89.7	90.1	91.7	91.3	93.3	¹³ C ₅ -PFHxA
PFHpA	90.6	89.8	92.6	91.3	91.9	¹³ C ₄ -PFHpA
PFOA	92.5	92	94.2	94.7	94.3	¹³C₈-PFOA
PFNA	93	92.2	94.3	94.8	95.2	¹³ C ₉ -PFNA
PFBS	94.6	92.1	96.8	93.8	96.1	¹³ C ₃ -PFBS
PFPeS	94.6	92.1	96.8	93.8	96.1	¹³ C ₃ -PFBS
PFHxS	89.8	88.1	91.3	91.9	93.5	¹³ C ₃ -PFHxS
PFHpS	92.8	90.6	94.6	94.1	93.5	¹³ C ₈ -PFOS
PFOS	92.8	90.6	94.6	94.1	93.5	¹³C₈-PFOS
PFNS	92.8	90.6	94.6	94.1	93.5	¹³ C ₈ -PFOS

Water Sample Reproducibility

- Reagent Water
- Surface (River) Water
- Ground Water
- Influent Waste Water
- Effluent Waste Water

	% RSD in Matrix			
	Ground Water		Surface Water	
	Low spike (n=3)	High spike (n=3)	Low spike (n=3)	High spike (n=3)
PFBA	15.6	2.7	2.1	2.8
PFPeA	7.4	1.6	2.2	0.4
PFHxA	6.6	1.7	1.6	1.9
PFHpA	4.4	2.4	8.9	2.1
PFOA	5.9	1.5	9.5	2.5
PFNA	6.3	2.5	8.2	1.0
PFDA	8.7	1.3	5.2	1.9
PFBS	6.2	1.2	1.5	1.5
PFPeS	2.5	0.3	1.4	0.1
PFHxS	2.6	1.7	2.5	4.1
PFHpS	6.2	3.8	2.6	1.6
PFOS	7.9	3.1	7.2	1.5
PFNS	4.0	2.8	5.6	1.3

7979-17

PFAS's Summary

- Occurrence studies have led to increased PFAS's testing needs
- Drinking water analysis is increasing steadily due to increased public awareness
- Other water types (non-drinking water) of increasing interest as well
- SPE enrichment allows entry level TQ to achieve ppt detection. HP TQ required for LVI/DI water analysis
- MRM database, analyte optimization tool and isolator column facilitate implementation



Thank you for your attention!

Questions???