

Emerging Per/Polyfluoroalkyl Substance (PFAS) Analysis in Water, Soil & Biota


NEMC Conference 2018

Tarun Anumol, Ph.D.
Global Environmental Industry Manager
Agilent, Wilmington, DE, USA




Terminology

Perfluoroalkyl substance



ALL H atoms linked to C in alkyl chain are substituted with F

Polyfluoroalkyl substance



SOME (but not all) H atoms linked to C in alkyl chain are substituted with F

Per/Polyfluoroalkyl substances

Unique Properties

- Thermal & Chemical stability: grease-proof food packaging, stain repellents
- Zwitterionic properties: surfactants
- Surface-tension lowering: fire-fighting foams



Contamination in Food and Water

Concern Grows Over Tainted Drinking Water

Vermont, New Hampshire and New York expand efforts to find out how much of a potentially toxic chemical is in drinking water

England since August 2014, when a resident of Hoosick Falls, N.Y., near the Vermont border, tested his drinking water and found high levels of the acid. The man was concerned because his father, a former employee of the town's plastics plant that used PFOA, died of cancer. The state in March sampled PFOA levels up to 620 parts per trillion in private Litchfield wells, well above the 100-parts-per-trillion level at which New Hampshire officials start to consider the amount unsafe. Tests in Merrimack measured as high as 1,600 parts per trillion.

environment

Williamtown water contamination highlights dangers of PFOS and PFOA

Qantas faces \$180,000 fine over toxic foam spill at Brisbane Airport

21 April 2017 - 04:56pm - First published 21 April 2017 - 11:02am
By Ruth McCosker

Qantas has been hit with an investigation notice while residents continue to be warned not to eat seafood following a toxic spill into Brisbane's waterways.

On Monday, April 10, 22,000 litres of a firefighting foam containing perfluorinated compounds was spilled from a Qantas hangar at Brisbane Airport.

f t e A A A



Volume 94 Issue 20 | pp. 20-22
Issue Date: May 16, 2016 | Web Date: May 11, 2016

Perfluorinated chemicals taint drinking water

Are there toxins in your fast food packaging?



QLD News

Brisbane River seafood warning still in effect

Chris Honnery, The Courier-Mail
April 20, 2017 8:49pm

Regulatory Guidance for PFAS

Compound	Agency/State	Concentration (ng/L)	Year
PFOA + PFOS	EPA DWHA	70	2016
PFOA	MN	35	2017
PFOS	MN	27	2017
PFOA	NJ	14 (40)	2017 (2007)
GenX	NC	140	2017
PFOA/PFOS	UK	300 (Reg. 10)	2009
PFOS/PFHxS	AUS	70	2016
PFOS	VIC EPA	0.23	2016

Standard Methods

EPA Method 537

DRINKING WATER

EPA Document #: EPA/600/R-08/092

METHOD 537. DETERMINATION OF SELECTED PERFLUORINATED ALKYL
ACIDS IN DRINKING WATER BY SOLID PHASE EXTRACTION
AND LIQUID CHROMATOGRAPHY/TANDEM MASS
SPECTROMETRY (LC/MS/MS)

Solid Phase Extraction

10 µL injection (96/4: MeOH/Water)

14 analytes

J.L. Shonkoff US EPA, Office of Research and Development, National Exposure
Research Laboratory
P.E. Grimmett US EPA, Office of Research and Development, National Exposure
Research Laboratory
B.K. Boutin The National Council on Aging, Senior Environmental Employment
Program

MDLs: 0.5 – 6.5 ng/L

ASTM Method 7979



SURFACE & WASTEWATER, SLUDGE

Standard Test Method for
Determination of Perfluorinated Compounds in Water,
Sludge, Effluent and Wastewater by Liquid
Chromatography Tandem Mass Spectrometry (LC/MS/MS)¹

LLE with 5 mL MeOH + pH adjustment

30 µL injection (50/50: MeOH/Water)



21 analytes

21 min run time

MDLs: 1.3 – 107 ng/L (1.3 – 11.6 ng/L for 537)

PFAS Analysis – LC Instrument Setup

Eliminate Background Contamination



Potential Contamination Sources

- Solvents
- Filtration apparatus
- Teflon lined tubing



[Anumol et al.](#), Recommended Plumbing Configurations for Reduction in Per/Polyfluoroalkyl Substance Background with Agilent 1260/1290 Infinity (II) LC's, Agilent Application Note (5991-7863EN)



Extracting PFAS from matrix

Water

Weak Anion Exchange (WAX) Cartridge





WAX SPE with sorbent mass of 30, 60 and 150 mg available.

For PFAS in water, the 6 cc, 150 mg with 30 μ m particle size is recommended.

Foods & Biologicals


Enhanced Matrix Removal-Lipid (EMR-L)



EMR-L specifically retains lipids while letting PFAS pass through

Available in powder and flow-through cartridge (1, 3 and 6 cc) format

Solid Phase Extraction USEPA 537 vs Modified



Experimental: Procedure: Phospholipid Removal Evaluation:

Protein Precipitation:

Add 400 μ L of Acetonitrile with 1% Formic Acid into a test tube



Add 100 μ L of blank plasma, pre-spun



Vortex on a Heidolph Multi Reax® at 800-1000 rpm, 5 minutes



Centrifuge at 5000 rpm, 5 min



Pipette the supernant to an autosampler vial for analysis



Inject directly into LC/MSMS

Captiva EMR-Lipid:

Add 400 μ L of Acetonitrile with 1% Formic Acid to **Captiva EMR-Lipid 1 mL cartridge**



Add 100 μ L of blank plasma, pre-spun



In-well mixing. Allow PPT.



Pull low vacuum, 2-4 psi, 1 drop/3-5 sec



Collect in polypropylene test tubes



Inject directly into LC/MSMS, using polypropylene* autosampler vials

* Polypropylene collection tubes and autosampler vials are recommended since PFASs will stick to glass

Results and Discussion: Recovery and RSD





Figure 3. Recovery and RSD for the PFASs evaluated in the study at 5 and 20 ng/mL

The recovery and relative standard deviation for the 22 PFASs were determined at 5 and 20 ng/mL shown in Figure 3. The overall recovery was between 75-125% for both 5 and 20 ng/mL. Relative standard deviation was 0.8-14% for 5 and 20 ng/mL.

Results and Discussion:

Phospholipid Removal by Captiva EMR-Lipid versus Protein Precipitation:




- Captiva EMR-Lipid, a novel phospholipid removing sorbent available in a SPE cartridge format allows for in-situ protein precipitation and phospholipid removal as the extract passes through the sorbent during elution.
- Figure 4 shows the overall phospholipid removal of Captiva EMR-Lipid when compared to protein precipitation, overlay.
- Captiva EMR-Lipid removes ~99% of the phospholipids based on peak area determined from the LC/MSMS MRM method for 11 phospholipids.

PFOA – 1 μ L injection of 24 fg/ μ L




RSD = 9%

IDL = 6.4 fg




No inject



Methanol
blank

52 PFAS in one LC-MS/MS method (5 ng/mL chromatogram)



Compound
classes:

PFCAs

PFSAs

ADONA

PFPAs

n:2 FTS

FASAs

FASAAs

FASEs


diPAPs

diSAmPAP

PFPis



Calibration curves PFOA and PFOS




PFOA:

2 μ L injection in MeOH

linear ($r^2 = 0.999$)

Range 0.1 to 100 ng/mL




PFOS:

2 μ L injection in MeOH

linear ($r^2 = 0.999$)

Range 0.1 to 100 ng/mL

Calibration curves 6:2 FTS and 6:2 diPAP




6:2 FTS:

2 µL injection in MeOH

linear ($r^2 = 0.999$)

Range 0.1 to 25 ng/mL



6:2 diPAP:

2 µL injection in MeOH

linear ($r^2 = 0.999$)

Range 0.1 to 100 ng/mL

Fish Extract Recovery & Data

Fish Extract Recovery & Data

Acknowledgements

Agilent Technologies

Peter Abrahamsson
Sheher Mohsin
Michael Woodman
James Pyke
Courtney Milner
Craig Marvin

RMIT University

Bradley Clarke
Timothy Coggan