

# Low-cost, rapid and *in situ* accurate quantification of chloramines and ammonia





National Environmental  
Monitoring Conference 2018


**Merwan Benhabib, PhD**  
VP Engineering




# Chlorine + Ammonia



Monochloramine




Dichloramine



Trichloramine or nitrogen trichloride

**Rate of formation depends on pH, temperature, time, and initial Cl<sub>2</sub> : NH<sub>3</sub> ratio**

# Chloramination for sanitation



Appearance	Name	Molecular Weight	Preferred pH Value	Biocidating Effect
$NH_2Cl$	Monochloramine	52	> 7	Good
$NHCl_2$	Dichloramine	85	4-7	Tolerable
$NCI_3$	Trichloramine	119	1-3	Average
$RNHCl$	Organic chloramines	Varies	Unknown	Bad

# Chloramination for sanitation


Inorganic chloramines are **alternate disinfectants to chlorine** that

- Are recognized as a **safer disinfectant**
- forms significantly **lower amount of toxic disinfectant byproducts**
- does not **contribute to taste or odor**
- **more stable in** solution than free residual chlorine

Appearance	Name	Molecular Weight	Preferred pH Value	Biocidating Effect
<i>NH<sub>2</sub>Cl</i>	Monochloramine	52	> 7	Good
<i>NHCl<sub>2</sub></i>	Dichloramine	85	4-7	Tolerable
<i>NCI<sub>3</sub></i>	Trichloramine	119	1-3	Average
<i>RNHCl</i>	Organic chloramines	Varies	Unknown	Bad

# Why optimize

- Control the ratio chlorine : nitrogen
- Avoid nitrification
- Guarantee treatment performance
- Respect EPA guidelines of 4 mg/L for chloramine



=> A method capable of continually, accurately and easily measure chloramine and ammonia is required

# Current monitoring procedure

## GC/MS :

- Expensive
- Skilled labor required
- Dedicated system and team
- Time consuming
- Not deployable



## Colorimetric (EPA method 350.1):

using ferrous diethyl-p-phenylenediamine (DPD) or leucocrystal violet (LCV) dyes:

- Interferences (Ca, Mg, organic amines, color absorbing compound)
- Lower performance than GC
- Speciation is technique dependant





# Raman spectroscopy

1920s



1990s



2000s



Optical telecommunications drove technology needed for portable, in-line, and compact Raman spectroscopy


But Raman is weak

- $\sigma_{NR} \sim 10^{-30} \text{ cm}^2/\text{molecule}$
- **1 in 10 million photons**



# Surface-Enhanced Raman Scattering (SERS)


- At a rough metal surface
  - Increased field intensity
  - Which means increased Raman signal
- SERS activity quantified by Enhancement Factor
  - EF range: 1 -  $10^{10}$





# Multi-particle effects


Nanoparticles (gold/silver) enable **ppb-level** detection



“Hot Spots”

# Analyte / substrate interaction

- The SERS effect requires an interaction at the surface – within a couple nanometers
- Analyte must
  - Interact with the substrate
  - Interact with a linker molecule
  - Change the properties of another SERS-active molecule
- Gold/Silver nanoparticles: many options to control surface properties



# Why not widely used?

- Achilles' heel: reproducibility
  - Variations in substrate properties
  - Stochastic nanoparticle alignments
- One reviewer: "SERS doesn't work"

**OndaVia has made SERS a quantitative, repeatable method using:**


- Internal standards
- Nanoparticle structure
- Surface modifications
- Intelligent software

# Pioneering colorimetric SERS

- Most colorimetric dyes have an amine group -> affinity with SERS substrate
- SERS coupled with internal standards can be applied to many colorimetric tests:
  - Improve detection limits, accuracy
  - Eliminate sample blanks, interferences
  - Add speciation capability
- Examples:
  - Selenium, ammonia, chlorinated solvents, lead, alcohols



# Colorimetric SERS, Berthelot reaction




# Colorimetric SERS, Berthelot reaction

- Addition of Bleach => Ammonia measurement
- No Bleach => Chloramine measurement

Berthelot dye:

Salicylate instead of phenol


-> Not hazardous





# Colorimetric-SERS enables trace-level lab-grade measurement

**Internal standard:**  $^{15}\text{NH}_3 \Rightarrow$  self-calibration



## A 15-min standard operating procedure

1. Mix 0.5 ml sample with 0.5 ml I.S.
2. Add 50 ul of reagents (mixture of sodium salicylate, nitroprusside , citrate and base)
3. Add Bleach for Ammonia measurement
4. Heat for 10 min in 70C water bath
5. Mix with Nanoparticles
6. Measure

# Raman spectrometer equipment


## Instrumentation

- 785-nm, 60-mW at substrate
- Cooled (-20°C) CCD detector
- 200-2000-cm<sup>-1</sup>, 4-cm<sup>-1</sup> resolution




# Lab Standard Calibration Curve

- Range of interest 0 - 15ppm (as N)
- Accurate quantification RSD<5%. (Internal standard => self-calibration)



# Quantitative analysis in tap water



# Ultra compact solid state spectrometer


## Instrumentation

- 785-nm, 50-mW, 10  $\mu\text{m}$  slit, No temp control
- APS-CMOS Detector
- 200-2000- $\text{cm}^{-1}$ , 5- $\text{cm}^{-1}$  resolution






Relative intensity at  
 $1530\text{cm}^{-1}$



Known concentration (ppm as N)



# Summary & Future Directions

- The successful results obtained in lab settings for standards and field samples demonstrates:
  - No pretreatment necessary
  - Low-cost, rapid and deployable
  - Repeatable accuracies similar to IC
  - Robust to interferents (Ca, Mg, amines)

Next :

- Pilot: looking for users
- Continue to work on the portability of the analysis
- Develop an In-line system
- Measure of Di-Trichloramines (CDC SBIR, chlorine in pool => asthma)

# Special thanks to...



Kristle Cruz Garcia  
Mayra Zaragoza  
Melissa Yao  
Mark Peterman, PhD  
George Janssen

## And to...



Questions?



Thanks to...

**NEMC**

Merwan Benhabib

merwan@ondavia.com

<http://www.ondavia.com>