

Implementing EPA Method 6020 with the Agilent ICP-MS Portfolio

NEMC August 2018

Craig Jones
Applications Specialist
ICP-MS, ICP-MS/MS




The Agilent Atomic Spectroscopy Lineup

An Instrument for Any Application!



55B AA



240FS/280FS Flame AA



240Z/280Z Furnace AA



4200 MP-AES



5110 ICP-OES



7800 ICP-MS



7900 ICP-MS



NEW! 8900 ICP-QQQ

Leading the way in atomic spectroscopy innovation

www.agilent.com/chem/atomic

Agilent's History in ICP-MS

30 Years of ICP-MS Innovation

Enabling high
sensitivity
metal analysis
PMS series



First computer-
controlled ICP-MS

1987

Enabling
routine robust
ICP-MS analysis
4500



First benchtop ICP-MS
Cool plasma

1994

Enabling control
of common
interferences
7500



9 orders detector
ORS cell

2000

Enabling ease
of use and
productivity
7700



HMI
ISIS-DS
MassHunter SW

2009

Enabling controlled
reaction chemistry
8800 ICP-QQQ



World's first ICP-QQQ

A new era in ICP-
MS performance
7900



UHMI
ODS detector
ISIS 3

2012

Flexible, high
performance MS/MS
New 8900 ICP-QQQ



Second generation
ICP-QQQ

Enabling simplified
ICP-MS workflows
7800



Solution ready
Method automation

2015

2016

#1 selling ICP-MS !

Most Compact Instruments on the Market



Agilent 7800



Agilent 7900



Agilent 8800



Agilent 8900

Some of the Benefits of ICP-MS

- Multi-elemental technique
- High sensitivity, ppq/ppt for most elements
- Short analysis time (~ 3 min to <60sec)
- Extremely linear with wide dynamic range
- Minimum number of interferences
- High Productivity



Most Important Performance Consideration for Environmental Analysis

- Sensitivity
- Interference Removal
- Matrix Tolerance
- Linear Dynamic Range



4 Key Benefits of Agilent ICP-MS


1) Matrix Tolerance



Ultra High Matrix Introduction System (UHMI)

Analyze Higher TDS samples than any other instrument on the market.


2) Interference Removal



Octopole Reaction System (ORS⁴) He Mode

Simple, effective removal of all common interferences.


3) Dynamic Range



Widest Dynamic Range

Up to 11 orders dynamic range:
0.1ppt(DL) to 10,000ppm

4) Mass Hunter Software



4 Key Benefits of Agilent ICP-MS


1) Matrix Tolerance



Ultra High Matrix Introduction System (UHMI)

Analyze Higher TDS samples than any other instrument on the market.


2) Interference Removal



Octopole Reaction System (ORS⁴) He Mode

Simple, effective removal of all common interferences.


3) Dynamic Range








Widest Dynamic Range


Up to 11 orders dynamic range:
0.1ppt(DL) to 10,000ppm

4) Mass Hunter Software






Ionization Suppression as a Function of Oxides



HMI – Effect on Matrix Suppression in Undiluted Seawater (~3% NaCl)



HMI dilutes aerosol density & water vapor, as well as sample matrix.

Gives much higher plasma temp; much better matrix decomposition. Matrix suppression is almost eliminated.

Plot shows % recovery in undiluted seawater vs aqueous calibrations.

With HMI, ALL results are within +/- 15% recovery (shaded area)

Without HMI, ALL results are below 60% recovery

Sample Introduction – HMI



US007671329B2

(12) **United States Patent**
Sakata et al.

(10) **Patent No.:** **US 7,671,329 B2**
(45) **Date of Patent:** **Mar. 2, 2010**

(54) **INDUCTIVELY COUPLED PLASMA MASS
SPECTROMETER**

(75) Inventors: **Kenichi Sakata**, Tokyo (JP); **Noriyuki
Yamada**, Tokyo (JP)

(73) Assignee: **Agilent Technologies, Inc.**, Santa Clara,
CA (US)

(*) Notice: Subject to any disclaimer, the term of this
patent is extended or adjusted under 35
U.S.C. 154(b) by 231 days.

(21) Appl. No.: **11/825,702**

(22) Filed: **Jul. 9, 2007**

(65) **Prior Publication Data**

US 2008/0035844 A1 Feb. 14, 2008

(30) **Foreign Application Priority Data**

Aug. 11, 2006 (JP) 2006-219520

5,185,523 A *	2/1993	Kitagawa et al.	250/281
5,308,977 A *	5/1994	Oishi et al.	250/288
5,495,107 A *	2/1996	Hu et al.	250/281
5,519,215 A *	5/1996	Anderson et al.	250/288
6,265,717 B1 *	7/2001	Sakata et al.	250/289
6,936,787 B2 *	8/2005	Tao et al.	219/121.51
2004/0195218 A1 *	10/2004	Tao et al.	219/121.51
2006/0024199 A1 *	2/2006	Tao et al.	422/54
2006/0087651 A1 *	4/2006	Montaser et al.	356/338
2008/0099671 A1 *	5/2008	Sakata et al.	250/288

FOREIGN PATENT DOCUMENTS

JP	10-188877	7/1998
JP	10-208691	8/1998
JP	11-006788	1/1999

* cited by examiner

Primary Examiner—David A Vanore

Assistant Examiner—Michael Maskell

(74) *Attorney, Agent, or Firm*—Marc Bobys

Introducing the Ground-Breaking Agilent UHMI



- UHMI (Ultra High Matrix Interface) – much more than just a simple T-piece
- UHMI uses optimized gas mixing geometry and sophisticated plasma/gas-flow tuning algorithm to set reproducible conditions for predictable aerosol dilution rate

“Big Four” Spiked into Different Salt Matrices

NaCl Amount	75 As [25 ppb]	114 Cd [50 ppb]	208 Pb [50 ppb]	201 Hg [1 ppb]
0%	26.9	49.2	49.7	0.85
0.5%	24.2	49.0	50.1	0.99
1%	24.8	51.5	50.2	0.93
1.5%	25.5	50.0	50.5	0.88
2%	24.6	50.0	49.7	1.03
5%	25.4	48.7	50.7	0.89
10%	22.8	46.1	49.8	0.91
25%	26.2	45.4	49.0	0.96
Average	25.1	48.7	50.0	0.93
% Recovery	100%	97%	100%	93%
% RSD	5%	4%	1%	6%



Data supplied by Wim Proper, Eurofins
Analytico, NL
Implementing EPA Method 6020 with the Agilent ICP-MS Portfolio



Interfered Elements Spiked into Different Salt Matrices

NaCl Amount	51 V [50 ppb]	52 Cr [50 ppb]	60 Ni [50 ppb]	63 Cu [50 ppb]
0%	49.2	49.1	49.9	49.6
0.5%	47.3	50.3	48.3	48.6
1%	49.5	49.3	48.8	48.8
1.5%	50.5	50.3	49.9	49.6
2%	49.7	49.1	49.4	48.7
5%	48.9	50.3	47.2	47.7
10%	47.8	50.3	46.3	47.7
25%	48.0	48.7	50.9	50.5
Average	48.9	49.7	48.8	48.9
% Recovery	98%	99%	98%	98%
% RSD	2%	1%	3%	2%



Data supplied by Wim Proper, Eurofins
Analytico, NL

Implementing EPA Method 6020 with the Agilent ICP-MS Portfolio



4 Key Benefits of Agilent ICP-MS


1) Matrix Tolerance



Ultra High Matrix Introduction System (UHMI)

Analyze Higher TDS samples than any other instrument on the market.


2) Interference Removal



Octopole Reaction System (ORS⁴) He Mode

Simple, effective removal of all common interferences.


3) Dynamic Range



Widest Dynamic Range

Up to 11 orders dynamic range:
0.1ppt(DL) to 10,000ppm

4) Mass Hunter Software



4 Key Benefits of Agilent ICP-MS


1) Matrix Tolerance



Ultra High Matrix Introduction System (UHMI)

Analyze Higher TDS samples than any other instrument on the market.


2) Interference Removal



Octopole Reaction System (ORS⁴) He Mode

Simple, effective removal of all common interferences.


3) Dynamic Range



Widest Dynamic Range


Up to 11 orders dynamic range:
0.1ppt(DL) to 10,000ppm

4) Mass Hunter Software



What about Interferences?


- Three main types of interferences
 - Spectroscopic
 - Physical
 - Memory
- How to control them?
 - Cell technologies
 - HMI / UHMI
 - Discrete sampling techniques




Troublesome Region of the Periodic Table

H																	He
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	L	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	A															
		L	La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
		A	Ac	Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr


Polyatomic Interference Formation - Matrix




Polyatomic Interference Formation - Argon




Polyatomic Interference Formation - Both



ICP-MS with ORS Technology



Processes of Interference Removal in Collision/Reaction Cell



Collisional Dissociation

- Limited in ICP-MS, as collision energy must be higher than bond dissociation energy


Reaction

- Can be very efficient – up to 9 orders reduction – but can also be non-selective. Highly reactive gases may react with analytes, matrix components and residual cell contamination, giving analyte loss and the formation of complex cluster ions

Energy Discrimination


- Useful in complex, variable and unknown matrices, as interference removal occurs, regardless of the level, source and chemistry of the interfering species. Can use inert cell gas, so no reaction with analytes and no formation of new cluster ions

Principle of Cell Gas Mode and KED for removing polyatomic interferences*



Polyatomic Interferences in No Gas Mode


Color of spectrum indicates which matrix gave each interfering peak



No Gas Mode

Polyatomic Interferences in He Mode


Color of spectrum indicates which matrix gave each interfering peak




He Mode

Matrix Mix with Spike (10ppb) in He Mode


Consistent sensitivity and perfect template match for all elements



Cal 0, 1, 10 ppb Arsenic in 1% HNO₃/0.5% HCl (ArCl interference on m/z 75)



Helium Mode Eliminates ArCl Interference on As



	Rjct	Conc.	Calc Conc.	CPS	Ratio	Det.	RSD
1	☐	0.000	0.000	48.33	0.0007	P	12.
2	☐	0.100	0.100	762.35	0.0106	P	7.
3	☐	1.000	1.017	7487.12	0.1021	P	1.
4	☐	5.000	5.041	36672.57	0.5033	P	0.
5	☐	10.000	9.978	73250.29	0.9954	P	1.
6	☐						

$$y = 0.0997 * x + 6.6377E-004$$

$$R = 1.0000$$

$$DL = 0.002449$$

$$BEC = 0.006658$$

Weight: <None>

Min Conc: <None>

Calibration standards were prepared in 2% HNO₃ / 1% HCl (in order to stabilize Hg). Notice complete elimination of ArCl.

Largest Analytical Range of any ICP-MS

Calibration ranges

Hg (10 – 200ppt) – NoGas Mode

As (10 – 200 ppt) – He Mode

Se (10 – 200 ppt) – He Mode

Na (0.05 – 1000 ppm) – He Mode

Overall calibration range 10ppt (Hg, As, Se)
to 1000 ppm (Na) in a single method
- *without* attenuating ion transmission to
increase working range

Na

Typically, ICP-MS cannot measure
above 200ppm Na without changing quad
resolution or ion lens settings





Hg

Hg LOD on 7800 is about 2ppt – 7800 can
QUANTITATE at 10ppt!

7800: 10 orders (1000's ppm)


7900: 11 orders (% Level)

These 4 plots were obtained under the same analytical conditions on
the 7800 – only the gas mode (NoGas for Hg) changed



Far Wider Dynamic Range Than Any Other ICP-MS 11 orders - low and high level calibrations in one run


Cd (1ppt - 1ppb) and Na (100ppb - 10,000ppm (1%)) in the same run



Analyte Information (111 Cd [No Gas])						
Current Sample	Calc Conc.	CPS	Ratio	Det.	Conc. RSD	
1ppbBeAsCdHgPb	1000.006	50630.23		P	0.4	
Calibration						
Level	Rjct	Conc.	Calc Conc.	CPS	Ratio	RSD
1		0.000	0.000	3.00	P	47.1
2		1.000	0.960	51.60	P	14.9
3		10.000	10.078	513.21	P	3.8
4		50.000	49.851	2526.83	P	2.9
5		100.000	100.008	5066.11	P	2.2
6		1000.00	1000.006	50630.23	P	0.4
7						
8						
9						
10						
11						
12						
13						
14						

Both calibrations are linear.

Total concentration range covered from Cd blank (BEC of <0.1ppt) to Na top standard (1%) is 11 orders



Analyte Information (23 Na [HMI He] ISTD: 115 In [HMI He])						
Current Sample	Calc Conc.	CPS	Ratio	Det.	Conc. RSD	
rinse	55.072	3281161.01	52.9048	A	31.7	
Calibration						
Level	Rjct	Conc.	Calc Conc.	CPS	Ratio	RSD
1		0.000	0.000	1944.19	0.0477	5.8
2		0.100	0.109	6393.55	0.1528	2.9
3		1.000	1.049	43220.26	1.0541	1.3
4		10.000	10.192	421823.28	9.8296	3.2
5		100.000	99.124	4004216.88	95.1848	1.7
6		1000.00	1089.335	40455458.40	1,045.570	1.5
7		10000.0	9991.075	328778051.1	9,589.288	3.7
8						
9						
10						
11						
12						
13						
14						

Enviro: 7800 HMI / Estuarine Samples: 40-75% Seawater

		~1745 ppm Na in Solution	~860 ppm Na in Solution	~430 ppm Na in Solution	
Acq. Date-Time	Sample Name	5	10	20	10/20/2009 5:59 PM
Sample Name Dilution	Dilution				Sample C 20x 20
9 Be [NoGas]	9 Be [NoGas]	0.0000	0.0000	0.0000	0.0418
23 Na [He]	23 Na [He]	8726108	8611994	8628625	4432966
24 Mg [He]	24 Mg [He]	1040282	1031898	1058318	544663
27 Al [He]	27 Al [He]	1.6166	0.6615	3.7364	1.6384
39 K [He]	39 K [He]	380963	371819	365013	183981
44 Ca [He]	44 Ca [He]	414310	409728	404271	214749
51 V [HEHe]	51 V [HEHe]	1.4565	1.2708	1.2975	2.9053
52 Cr [He]	52 Cr [He]	0.2165	0.3066	0.3074	0.1845
55 Mn [He]	55 Mn [He]	1.0943	1.1272	0.9813	1.7931
56 Fe [He]	56 Fe [He]	<2.5000	<5.0000	<10.000	<10.000
59 Co [He]	59 Co [He]	0.0221	0.0234	0.0266	0.0515
60 Ni [He]	60 Ni [He]	0.2081	0.2032	<0.6000	1.1988
63 Cu [He]	63 Cu [He]	0.2909	0.2419	0.2371	1.4446
66 Zn [He]	66 Zn [He]	0.0254	0.1200	0.4770	1.3007
75 As [He]	75 As [He]	1.3844	1.4452	1.4770	1.7367
78 Se [HEHe]	78 Se [HEHe]	0.3003	0.3113	<0.4000	<0.4000
95 Mo [NoGas]	95 Mo [NoGas]				5.4514
107 Ag [NoGas]	107 Ag [NoGas]				<0.1000
111 Cd [NoGas]	111 Cd [NoGas]				0.0497
121 Sb [NoGas]	121 Sb [NoGas]				0.4872
137 Ba [NoGas]	137 Ba [NoGas]				25.616
201 Hg [NoGas]	201 Hg [NoGas]				0.1662
205 Tl [NoGas]	205 Tl [NoGas]				<0.1000
208 Pb [NoGas]	208 Pb [NoGas]				<0.5000
232 Th [NoGas]	232 Th [NoGas]				<0.0600
238 U [NoGas]	238 U [NoGas]				1.4673
		~0.275 ppb As in Solution	~0.144 ppb As in Solution	~0.074 ppb As in Solution	

7900 Heavy Matrix Analysis (HMI-25)

-Major constituents reported in ppm

-Trace elements reported in ppb

General Seawater Composition


Element	At.weight	ppm
Hydrogen H ₂ O	1.00797	110,000
Oxygen H ₂ O	15.9994	883,000
Sodium NaCl	22.9898	10,800
Chlorine NaCl	35.453	19,400
Magnesium Mg	24.312	1,290
Sulfur S	32.064	904
Potassium K	39.102	392
Calcium Ca	10.08	411
Bromine Br	79.909	67.3

Most analytes acquired using Helium mode.


H2 Mode used for select elements to provide the best DL's (Si, Ca, Se)

6020 Unknown Results ppb (Na,Mg,S,K,Ca,Fe ppm)			
Acq. Date-Time		9/8/2017 13:24	9/8/2017 13:27
Type		Sample	Sample
Level			
Sample Name		SAMPLE A	SAMPLE B
Total Dil.		1.11	1.11
9 Be [No Gas]	Conc. [ppb]	86.3314	83.1840
11 B [He]	Conc. [ppb]	15661.9829	16122.3051
23 Na [He]	Conc. [ppm]	1167.3130	1153.2386
24 Mg [He]	Conc. [ppm]	2812.7092	2786.6283
27 Al [He]	Conc. [ppb]	356519.5441	349253.6402
28 Si [H2]	Conc. [ppb]	48360.4526	47745.2191
34 S [He]	Conc. [ppm]	553.3986	537.7674
39 K [He]	Conc. [ppm]	61.2055	60.5476
40 Ca [H2]	Conc. [ppm]	483.6432	491.6683
51 V [He]	Conc. [ppb]	271.6009	247.0839
52 Cr [He]	Conc. [ppb]	23758.5895	22959.0748
55 Mn [He]	Conc. [ppb]	74152.6376	72835.0402
56 Fe [He]	Conc. [ppm]	544.9229	530.3951
59 Co [He]	Conc. [ppb]	1129.1326	1132.3477
60 Ni [He]	Conc. [ppb]	21512.9884	21352.7933
63 Cu [He]	Conc. [ppb]	5.3540	2.0105
66 Zn [He]	Conc. [ppb]	3902.0441	4234.4300
75 As [He]	Conc. [ppb]	13.1571	12.4169
78 Se [H2]	Conc. [ppb]	96.1427	114.5517
88 Sr [He]	Conc. [ppb]	800.8821	771.5266
95 Mo [He]	Conc. [ppb]	3.5913	2.3013
107 Ag [He]	Conc. [ppb]	0.1056	0.0166
111 Cd [He]	Conc. [ppb]	30.5723	39.1334
118 Sn [He]	Conc. [ppb]	0.3993	0.3039
121 Sb [He]	Conc. [ppb]	0.6053	0.3508
137 Ba [He]	Conc. [ppb]	11.3219	10.4092
205 Tl [He]	Conc. [ppb]	0.5324	0.6342
208 Pb [He]	Conc. [ppb]	0.1570	0.0861

Enviro: 7900 UHMI / Very High Matrix Samples (HMI)



7900 UHMI: Samples Containing ~12% TDS



*Calibration Standards in 2% HNO3/1% HCl

4 Key Benefits of Agilent ICP-MS


1) Matrix Tolerance



Ultra High Matrix Introduction System (UHMI)

Analyze Higher TDS samples than any other instrument on the market.


2) Interference Removal



Octopole Reaction System (ORS⁴) He Mode

Simple, effective removal of all common interferences.


3) Dynamic Range



Widest Dynamic Range

Up to 11 orders dynamic range:
0.1ppt(DL) to 10,000ppm

4) Mass Hunter Software



Mass Hunter 4.4: Truly Automated Software

Agilent INTERNAL USE ONLY!!! Online ICP-MS MassHunter - Pro Caps (UW 100mg HNO3 and Hot Plate) 020414.b

File Edit View Instrument Hardware Startup Batch Queue User Tune Tools Help Script

Hardware Plasma Batch Tune Queue Data Analysis Report Auto Sampler Settings

Instrument Status

Error: 2/2/2014 10:44:32 AM 2152. Warning: Analyzer pressure is too high.

Standby EMP

IF/BK Press 2.62E+0 Pa Analyzer Press 1.14E-5 Pa Water Temp 16.8 °C Internal Temp 37.2 °C Inlet Temp 26.1 °C

Hardware -> Mainframe -> Maintenance Log

Mainframe Autosampler(ALS) Sample Introduction Plasma Ion Lenses Collision/ Reaction Cell Quadrupole Detector

Add Maintenance Log Delete Maintenance Log

	Date and Time	Operator	Maintenance Log	
1	1/14/2014 5:06:43 PM	Matt Reuer	Miscellaneous	Installation of 7900 at Santa Clara COE

Startup

Add to Queue Select Custom Settings

Hardware Settings	On	Vial #	Monitored Masses	PeriPump
Torch Axis	☒	1307	7 89 205	Set
EM	☒	1307	80	Set
Plasma Correction	☐	1307		Set
Standard Lenses Tune	☒	1307	7 89 205	Set
Resolution/Axis	☒	1307	7 89 205	Set
Performance Report	☒	1307	Oxide: 140; Doubly Charged: 140	Set
Full Spectrum	☐	1307		Set
P/A Factor	☒	1306	7 89 205 6 9 72 23 24 27 45 51 52 55 59 60 6	Set

Standard Setting

Vial# 1307 Monitored Masses: 7 89 205 Set PeriPump/ISIS

[Standby] Load Option: Standard Tune

Daily Performance Log

747 565 118 - TeamViewer - Free license (non-commercial use only)

[Agilent INTERNAL USE ONLY!!!] Online ICP-MS MassHunter - Pro Caps (UW 100mg HNO₃ and Hot Plate) 020414.b

File Edit View Instrument Hardware Startup Batch Queue User Tune Tools Help Script

Hardware Plasma Batch Tune Queue Data Analysis Report Auto Sampler Settings

Instrument Status

Error: 2/2/2014 10:44:32 AM 2152, Warning: Analyzer pressure is too high.

Standby EMF

IF/BK Press 2.65E+0 Pa Analyzer Press 1.14E-5 Pa Water Temp 16.9 °C Internal Temp 37.4 °C Inlet Temp 26.3 °C

Hardware -> Mainframe -> Performance Report

Mainframe Autosampler(ALS) Sample Introduction Plasma Ion Lenses Collision/ Reaction Cell Quadrupole Detector


Run Performance Report History View Select Data Displayed in History View Export Clear Performance Report

Created Date	Sensitivity			Background			Standard Tune Parameters Ion Lenses			Standard T. Cell
	Channel 1 Count	Channel 2 Count	Channel 3 Count	Channel 1 Count	Channel 2 Count	Channel 3 Count	Extract 1	Extract 2	Omega Bias	
1/31/2014 9:01:55 AM	7848.41	57429.06	46811.01	1.00	0.55	3.35	0.0 V	-230.0 V	-100 V	0.0 mL/min
1/28/2014 12:38:36 PM	6033.13	57663.31	56772.38	1.20	0.25	2.45	0.0 V	-200.0 V	-100 V	0.0 mL/min
1/18/2014 7:43:58 AM	5234.87	46754.61	44780.73	1.45	0.30	1.25	0.0 V	-200.0 V	-100 V	0.0 mL/min
1/18/2014 7:34:49 AM	897.95	564.17	11943.48	0.00	0.00	0.00	0.0 V	-200.0 V	-80 V	5.0 mL/min
1/17/2014 11:19:59 AM	6196.21	45472.99	42122.55	1.55	0.25	1.90	0.0 V	-200.0 V	-100 V	0.0 mL/min
1/16/2014 9:38:52 AM	6015.78	41287.56	38027.58	0.85	0.45	1.25	0.0 V	-200.0 V	-100 V	0.0 mL/min
1/14/2014 6:09:42 PM	6169.79	41334.88	39337.92	1.00	0.55	1.70	0.0 V	-200.0 V	-100 V	0.0 mL/min

[Standby] Load Option: Standard Tune

10:45 PM 2/4/2014

Detail View of Autotune



Pre-Configured Methods for Common Applications

Select Prest Method

Application Method	Generic Method
Drinking Water (with He)	7900 Application Method for Drinking Water (with ORS)
ChP	7900 Application Method for Elemental Impurities in Pharma, using China Ph
USP<232>	7900 Application Method for Elemental Impurities in Pharma, using USP<232>
EPA200.8	7900 Application Method for EPA200.8 (No minerals, without ORS)
EPA6020	7900 Application Method for EPA6020

Compatible Sample Types:
Simple aqueous matrices and acid digested samples up to 0.4% dissolved solids.

Pre-Defined Analytes:
V, Cr, Ni, Cu, As, Mo, Ru, Rh, Pd, Cd, Os, Ir, Pt, Hg, Pb

Comment:
Method uses **General Purpose** plasma conditions and He mode only, but requires the **7900**. User Access Control option software is required for regulated plasma analysis. The selection of an appropriate internal standard should consider the analyte and the sample matrix.

Required Hardware:
7900, x-Lens, Micro Mist Nebulizer

OK Cancel



Select Elements on Periodic Table

Tune Mode: He

Selected Elements: 33

Mass	Element
45	Sc
51	V
52	Cr
53	Cr
60	Ni
62	Ni
63	Cu
65	Cu
72	Ge
75	As
95	Mo
97	Mo
101	Ru
103	Rh
105	Pd
111	Cd

Go to Mass Scale Element Information...

Correction Equations: 1

OK Cancel

Acq Parameters PenPump/ISIS Tune

Sample Introduction: General Tune Vial: <Manual>

Sample Acquisition	Time [sec]	Speed [ps] Nebulizer Pump	Vial#
Pre Run			
Sample Uptake	30	0.30	Sample
Stabilize	40	Tune Parameter	Sample
Acquisition			
Speed		Tune Parameter	Sample
Post Run			
Probe Rinse (Sample)	10	0.30	Rinse Port
Probe Rinse (Std)	10		Rinse Port
Rinse 1	60	0.10	<Manual>
Probe Rinse 1			<Manual>
Rinse 2			<Manual>
Probe Rinse 2			<Manual>
Rinse 3			<Manual>

Intelligent Rinse

☒ Preemptive Rinse

Start rinse before acquisition complete: 40 sec

☐ Move to the next rinse step at acquisition complete



Mass	Element Name	Monitor	Integ Time /Mass [sec]
45	Sc	☒	0.10
51	V	☒	0.50
52	Cr	☒	0.50
53	Cr	☒	0.50
60	Ni	☒	0.50
62	Ni	☒	3.00
63	Cu	☒	0.50
65	Cu	☒	0.50
72	Ge	☒	0.10
75	As	☒	3.00
95	Mo	☒	0.50
97	Mo	☒	0.50
101	Ru	☒	0.50
103	Rh	☒	0.50
105	Pd	☒	0.50
111	Cd	☒	1.00
114	Cd	☒	1.00
115	In	☒	0.10

Automated P/A Calibration

[Agilent INTERNAL USE ONLY!!!] Offline ICP-MS MassHunter - 022714_a.b

File Edit View Instrument Hardware Startup Batch Queue Tools Help

Hardware Plasma Batch Tune Queue Data Analysis Report

Batch - 022714_a.b

Save Batch Add to Queue Validate Method Select Elements Tune Mode: <All>

Acq Method Data Analysis Method Sample List

Acq Parameters PeriPump/ISIS Tune

Acq Mode

Spectrum

Spectrum Mode Option

Peak Pattern: 1 Point

Replicates: 3

Sweeps/Replicate: 50

Acq Option

☐ Auto/Semi Auto Tune before Batch


☐ Print Tune Report

☒ P/A Factor Adjustment

Total Acq Time

Tune Mode		#1: No Gas	#2: H2	#3: He				
Quick Scan		☐	☐	☒				
Independent P/A Factors		☒	☒	☒				
Stabilization Time [sec]		0	5	5				
Resolution		Standard	Standard	Standard				
Mass	Element Name	Monitor	+0.5 u	IntegTime /Mass [sec]	+0.5 u	IntegTime /Mass [sec]	+0.5 u	IntegTime /Mass [sec]
6	Li	☒	☐	0.30	☐	N/A	☐	0.30
7	[Li]	☐	☐	0.30	☐	N/A	☐	0.30
9	Be	☐	☐	1.00	☐	N/A	☐	1.00
11	B	☐	☐	1.00	☐	N/A	☐	1.00
27	Al	☐	☐	0.50	☐	N/A	☐	0.30
35	Cl	☐	☐	N/A	☐	N/A	☐	0.30
39	K	☐	☐	0.30	☐	N/A	☐	0.30
40	Ca	☐	☐	N/A	☐	0.30	☐	N/A

Seamless Pulse to Analog Transition




EPA Compliance Software

200.8 Tune Check Sample Report

Batch Folder C:\Agilent\ICPMH\1\DATA\epa tune check.b
Report Comment American West Tune 110923
Instrument Name G3281A JP10420651

[No Gas]



Mass	Peak Height	Axis (Actual)	Axis (Required)	Axis (Flag)	Width-X% (Actual)	Width-X% (Required)	Width-X% (Flag)
9	49907	9.00	8.9 - 9.1		0.742	0.900	
24	14962090	23.95	23.9 - 24.1		0.742	0.900	
25	2081180	25.00	24.9 - 25.1		0.769	0.900	
26	2515438	25.95	25.9 - 26.1		0.740	0.900	
59	502004	59.00	58.9 - 59.1		0.733	0.900	
115	5395881	115.00	114.9 - 115.1		0.734	0.900	
206	409162	205.95	205.9 - 206.1		0.737	0.900	
207	341574	206.95	206.9 - 207.1		0.728	0.900	
208	845273	207.95	207.9 - 208.1		0.730	0.900	

X% = 5 Integration Time [sec] = 0.1 Acquisition Time [sec] = 235 Y Axis = Linear

Individual Sample QC Limits (ie: ICSAB)

22	LLCCV	Analyte <table><thead><tr><th colspan="6"></th><th colspan="4">Main Criteria1</th></tr><tr><th>Tune Mode</th><th>Mass</th><th>Name</th><th>QC Check</th><th>Report</th><th>Lower Exp.</th><th>Lower Error Lim.</th><th>Upper Exp.</th><th>Upper Error Lim.</th></tr></thead><tbody><tr><td>11</td><td>3: He</td><td>52</td><td>Cr</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>12</td><td>1: No Gas</td><td>55</td><td>Mn</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>13</td><td>3: He</td><td>55</td><td>Mn</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>14</td><td>2: H2</td><td>56</td><td>Fe</td><td>☐</td><td>☐</td><td>250000</td><td>0.8</td><td>250000</td><td>1.2</td></tr><tr><td>15</td><td>3: He</td><td>56</td><td>Fe</td><td>☐</td><td>☐</td><td>250000</td><td>0.8</td><td>250000</td><td>1.2</td></tr><tr><td>16</td><td>1: No Gas</td><td>59</td><td>Co</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>17</td><td>3: He</td><td>59</td><td>Co</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>18</td><td>1: No Gas</td><td>60</td><td>Ni</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>19</td><td>3: He</td><td>60</td><td>Ni</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>20</td><td>1: No Gas</td><td>63</td><td>Cu</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>21</td><td>3: He</td><td>63</td><td>Cu</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>22</td><td>1: No Gas</td><td>65</td><td>Cu</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr><tr><td>23</td><td>1: No Gas</td><td>66</td><td>Zn</td><td>☒</td><td>☒</td><td>200</td><td>0.8</td><td>200</td><td>1.2</td></tr></tbody></table>							Main Criteria1				Tune Mode	Mass	Name	QC Check	Report	Lower Exp.	Lower Error Lim.	Upper Exp.	Upper Error Lim.	11	3: He	52	Cr	☒	☒	200	0.8	200	1.2	12	1: No Gas	55	Mn	☒	☒	200	0.8	200	1.2	13	3: He	55	Mn	☒	☒	200	0.8	200	1.2	14	2: H2	56	Fe	☐	☐	250000	0.8	250000	1.2	15	3: He	56	Fe	☐	☐	250000	0.8	250000	1.2	16	1: No Gas	59	Co	☒	☒	200	0.8	200	1.2	17	3: He	59	Co	☒	☒	200	0.8	200	1.2	18	1: No Gas	60	Ni	☒	☒	200	0.8	200	1.2	19	3: He	60	Ni	☒	☒	200	0.8	200	1.2	20	1: No Gas	63	Cu	☒	☒	200	0.8	200	1.2	21	3: He	63	Cu	☒	☒	200	0.8	200	1.2	22	1: No Gas	65	Cu	☒	☒	200	0.8	200	1.2	23	1: No Gas	66	Zn	☒	☒	200	0.8	200	1.2
						Main Criteria1																																																																																																																																																	
Tune Mode	Mass		Name	QC Check	Report	Lower Exp.	Lower Error Lim.	Upper Exp.	Upper Error Lim.																																																																																																																																														
11	3: He		52	Cr	☒	☒	200	0.8	200	1.2																																																																																																																																													
12	1: No Gas		55	Mn	☒	☒	200	0.8	200	1.2																																																																																																																																													
13	3: He		55	Mn	☒	☒	200	0.8	200	1.2																																																																																																																																													
14	2: H2		56	Fe	☐	☐	250000	0.8	250000	1.2																																																																																																																																													
15	3: He		56	Fe	☐	☐	250000	0.8	250000	1.2																																																																																																																																													
16	1: No Gas		59	Co	☒	☒	200	0.8	200	1.2																																																																																																																																													
17	3: He		59	Co	☒	☒	200	0.8	200	1.2																																																																																																																																													
18	1: No Gas		60	Ni	☒	☒	200	0.8	200	1.2																																																																																																																																													
19	3: He		60	Ni	☒	☒	200	0.8	200	1.2																																																																																																																																													
20	1: No Gas		63	Cu	☒	☒	200	0.8	200	1.2																																																																																																																																													
21	3: He		63	Cu	☒	☒	200	0.8	200	1.2																																																																																																																																													
22	1: No Gas		65	Cu	☒	☒	200	0.8	200	1.2																																																																																																																																													
23	1: No Gas		66	Zn	☒	☒	200	0.8	200	1.2																																																																																																																																													
23	LLICV																																																																																																																																																						
24	ICB																																																																																																																																																						
25	CCB																																																																																																																																																						
26	AllRef																																																																																																																																																						
27	SpkDup																																																																																																																																																						
28	Dup																																																																																																																																																						
29	PDS																																																																																																																																																						
30	ICSA																																																																																																																																																						
31	ICSB																																																																																																																																																						
32	LFM																																																																																																																																																						
33	LFMDup																																																																																																																																																						
34	LRB																																																																																																																																																						
35	LFB																																																																																																																																																						
36	QCS1																																																																																																																																																						
37	LRS																																																																																																																																																						
38	Dil																																																																																																																																																						
39	CCV2																																																																																																																																																						
40	Blank																																																																																																																																																						
41	PB																																																																																																																																																						
42	LCS																																																																																																																																																						
43	FRB																																																																																																																																																						
44	SpikeRef																																																																																																																																																						

ISTD (Method Name : IS)									
						Main Criteria1			
	Tune Mode	Mass	Name	QC Check	Report	Lower Exp.	Lower Error Lim.	Upper Exp.	Upper Error Lim.
1	1: No Gas	45	Sc	☒	☒	100	0.6	100	1.25
2	2: H2	45	Sc	☒	☒	100	0.6	100	1.25

Each QC Sample is Customized for its' individual QC requirements

ICSAB, CRDL, ICV, CCV, CCB, etc.

Associated Sample Report may also be Customized to Display all required Information

Unique Performance of the 7800 & 7900

Better matrix tolerance than any other ICP-MS

Higher plasma temperature (lower CeO/Ce ratio) under standard conditions than any other system

- 7800 with HMI: **~2-3% TDS** & 7900 with UHMI: **~25% TDS**

Best performance with Helium cell gas – eliminates need for reaction gases in all common applications

- 7800/7900 ORS⁴ – removes polyatomic interferences in complex matrices while maintaining excellent sensitivity.

Wider dynamic range than any other quadrupole ICP-MS

- 7800: 10 orders of Magnitude (~1000s ppm) (3 ms Integration)
- 7900: 11 orders of Magnitude (~Percent) (0.1 ms Integration spNP)



Maximizing Productivity for High Matrix Sample Analysis Using the Agilent 7900 ICP-MS with ISIS 3 Discrete Sampling System

Publication number: 5991-5208EN

Publication date: October 2014

Environmental

Kazuo Yamanaka¹ and Steve Wilbur²

¹Agilent Technologies, Japan

²Agilent Technologies, USA

Introduction

Faster analysis, improved ease-of-use, superior matrix tolerance

Combining the unique benefits from each of the following provides much faster analysis, improved ease-of-use and superior matrix tolerance:

- Agilent 7900 ORS⁴ He mode collision cell
 - For interference removal
- Ultra High Maxtrix Introduction (UHMI)
 - Aerosol dilution
- ISIS 3
 - Discrete sampling to achieve the best possible productivity

Experimental

Operating conditions

Table 1. Agilent 7900 ICP-MS operating conditions

ICP-MS parameters	No gas mode	He mode
RF power (W)	1600	
Carrier gas flow (L/min)	0.77	
Dilution gas flow (L/min)	0.28	
Lens tune	Autotune	Autotune
Cell gas flow (mL/min)	0.0	4.3
Energy discrimination (V)	5.0	
Number of elements	1 analyte, 1 ISTD	25 analytes, 6 ISTDs
Total acquisition time (3 reps) (seconds)	41	

Experimental

ISIS-DS operation




Figure 1. Overview of ISIS-DS operation. Valve in load position (upper left) and inject position (upper right). Actual analyte and internal standard signals during ISIS 3 operation are shown in lower plot with annotation (shown on next page).

Experimental ISIS-DS operation




Figure 1 continued. Actual analyte and internal standard signals during ISIS 3 operation are shown with annotation.

Results

This study used He cell mode for most elements and no gas mode for low mass elements.

383 samples were measured in 9 hours and 35 minutes with a run-to-run time of 90 seconds per sample.

Results

Meeting EPA QA/QC requirements




Figure 4. Internal standard recoveries normalized to the calibration blank for all samples. Due to limited space, not all sample names are shown. No internal standard failures occurred.

Results and Discussion

Recovery of certified reference values

Table 4. Mean measured values, recoveries, percent Relative Standard Deviations (%RSDs) for all certified elements in the six CRMs analyzed. Blank cells indicate no certified value.

* Concentration of 1/10 diluted solution. ND: less than detection limit

Isotope/ Element	NIST 1640a (n=7)			River Sediment A (1/10, n=14)*			River Sediment B (1/10, n=14)*		
	Mean conc. (ppb)	%RSD	Mean recovery (%)	Mean conc. (ppb)	%RSD	Mean recovery (%)	Mean conc. (ppb)	%RSD	Mean recovery (%)
9 Be	2.80	2.4	93	0.026	9.8		ND		
23 Na	3,112	2.0	107	5,326	1.8	107	5,610	1.6	112
24 Mg	1,062	1.7	101	7,375	1.8	105	13,130	1.3	109
27 Al	57.8	4.4	109	25,180	1.3	101	61,460	0.6	102
39 K	561	5.4	97	14,730	2.4	98	20,220	2.3	101
44 Ca	5,373	3.2	97	29,050	1.8	97	28,960	1.6	97
51 V	15.4	1.5	103	25.5	1.1	102	104	0.7	104
52 Cr	40.5	1.6	101	30,950	1.1	103	1,592	0.9	106
55 Mn	40.0	1.7	100	786	2.2	98	587	1.1	98
56 Fe	44.2	5.7	121	122,000	0.9	102	41,610	0.6	104
59 Co	20.7	2.1	103	12.1	1.5	121	16.3	1.3	109
60 Ni	25.7	1.8	103	52.2	1.4	104	52.8	1.4	106
63 Cu	87.9	2.0	103	101.9	1.7	102	104	1.4	104
66 Zn	56.1	1.3	102	1,454	0.9	97	493	0.8	99
75 As	7.99	1.4	100	60.2	1.3	100	20.4	1.4	102
78 Se	18.5	5.1	93	2.28	5.7	114	0.95	10.5	95
95 Mo	45.8	2.0	101	0.074	7.0		0.17	4.8	
107 Ag	8.25	2.8	103	ND			0.17	2.3	
111 Cd	3.97	1.8	100	10.2	1.1	102	3.18	2.4	106
121 Sb	5.07	1.8	100	50.7	0.7	102	4.22	0.9	106
137 Ba	147	1.1	98	49.3	0.7	99	392	0.5	98
201 Hg	ND			ND			0.023	10.9	
205 Tl	1.58	3.0	98	0.99	1.4	99	1.15	1.6	115
208 Pb	12.5	4.4	104	742	1.8	106	212	1.7	106
232 Th	0.002	40		2.04	2.6	102	9.93	2.3	99
238 U	25.2	3.8	100	1.02	3.1	102	3.02	2.9	101

Conclusions

Long sequences of high Total Dissolved Solids (TDS) samples can be analyzed with high accuracy, precision, and long term stability using the Agilent 7900 ICP-MS.

This was shown by combining the benefits of He mode for removing polyatomic interferences with the UHMI for highly robust plasma conditions, and the ISIS 3 discrete sampling system.

Questions?

Craig Jones

Craig_Jones@agilent.com

