

The Analysis of Soils and Waters in Accordance with U.S. EPA Method 6020B using the NexION[®] 2000 ICP-MS



Anthony Palermo



Introduction

- Liquid and soil wastes enter the environment
 - Quantity and diversity of wastes increasing
 - Contaminate soils and waters
- U.S. EPA Method 6020 developed to address contaminated soils and waters
 - Updated several times throughout the years
 - Applies to both waters and soils
 - Latest iteration is 6020B, which includes
 - New elements – 23 in total
 - New detection limit criteria
 - New QC criteria

Elements

Analyte	Symbol	Recommended Mass (amu)
Aluminum	Al	27
Antimony	Sb	121
Arsenic	As	75
Barium	Ba	137
Beryllium	Be	9
Boron*	B	11
Cadmium	Cd	111
Calcium	Ca	44
Chromium	Cr	52
Cobalt	Co	59
Copper	Cu	63
Iron	Fe	56
Lead [#]	Pb	208

Analyte	Symbol	Recommended Mass (amu)
Magnesium	Mg	24
Manganese	Mn	55
Mercury	Hg	202
Molybdenum*	Mo	95
Nickel	Ni	60
Potassium	K	39
Selenium	Se	78
Silver	Ag	107
Sodium	Na	23
Thallium	Tl	205
Vanadium	V	51
Zinc	Zn	66

- B, Mo: Not part of 6020B; analyzed for “information only”
- Pb: Sum Pb206 + Pb207 + Pb 208 isotopes

Why Sum the Pb Isotopes?


- Pb occurs naturally with 4 isotopes:

Pb Isotope	Natural Abundance
204	1.4%
206	22.4%
207	22.1%
208	52.4%

- Other elements undergo radioactive decay and turn into Pb
 - Th 232 \rightarrow Pb 208 \rightarrow Decay time = 1.4 billion years
 - U 235 \rightarrow Pb 207
 - U 238 \rightarrow Pb 206

} Decay Time = 4.46 billion years
- For highest accuracy, the Pb isotopes must be summed

Daily Work Flow



Samples & Sample Preparation

- Samples
 - Water & soil reference materials
 - NIST 1640a Trace Elements in Natural Water
 - HPS Soil Solution B
 - HPS River Sediment B
- Sample Preparation
 - Waters: Acidified to 2% HNO₃ + 1% HCl
 - Soils: 10x dilution with 2% HNO₃+ 1% HCl
 - 200 ppb Au added to aid with Hg washout

A Word on Soil Reference Materials

- Solid soil reference materials exist
 - For example NIST 2711a Montana II Soil
- The analytical results depend on sample preparation
 - Since soils contain Si, complete digestion requires HF
 - EPA methods for leaching soils, without the use of HF
 - Silicates not dissolved → not total metals
 - Various sample preparation methods available

HCl (mL)	HNO ₃ (mL)
0	10
3	9
9	3

Will give different results

A Word on Soil Reference Materials

- NIST 2711 Montana II Soil as an example
 - Drying Samples
 - Recommended at 110°C for 2 hours
 - But don't dry for volatile elements (i.e. As, Hg, Se)
 - Need to correct for moisture content
 - Total Results
 - Many from solid-sampling techniques (i.e. XRF, INAA, etc.)
 - Analyses done by NIST and USGS
 - Leach Results
 - Determined with ICP-OES and Method 3050b (hot plate digestion)
 - 6 commercial labs
 - % Recoveries relative to total composition: 1.5% (Na) to 100% (Hg)

**What is the correct, accepted concentration
in solid soil reference materials?**

Calibration Standards

Analytes	Standard 1 (µg/L)	Standard 2 (µg/L)	Standard 3 (µg/L)	Standard 4 (µg/L)
Be, B, V, Cr, Mn, Co, Ni, Cu, Zn, As, Se, Mo, Ag, Cd, Sb, Ba, Tl, Pb	1	10	50	100
Hg	0.1	1	5	10
Na, Mg, Al, K, Ca, Fe	100	1000	5000	10000
Ga, Ge, In, Ir	Internal Standards			

- Prepared in 2% HNO₃ + 1% HCl
- 200 ppb Au added to aid with Hg washout
- External calibration curves – linear through zero

Instrumental Conditions: NexION[®] 2000 ICP-MS


Parameter	Description/Value
Sample Delivery Rate	350 µL/min
Nebulizer	MEINHARD [®] plus Glass Type C
Spray Chamber	Glass cyclonic (baffled), 2 °C
Injector	2.0 mm i.d.
Nebulizer Flow	Optimized for < 2% oxides
RF Power	1600 W
Cones	Ni
Replicates	3
Aerosol Dilution	Set to 2.5x



- All elements measured in Collision mode
 - Increases confidence in results when matrices & interferences can vary from sample to sample

Instrumental Conditions: NexION[®] 2000 ICP-MS

- On-line aerosol dilution
 - Decreases matrix suppression
 - Set to 2.5x




Instrumental Conditions: NexION® 2000 ICP-MS

- 7-port switching valve with 1 mL sample loop
 - Increases sample throughput
 - Minimizes sample uptake and washout
- On-line addition of internal standards through switching valve
 - Increased productivity
 - Less chance of error
 - Enhanced mixing




Results: Reference Materials



All Recoveries within 10% of Certified Values

Results: Stability – CCV

- 8 hour analysis of soil samples
- Calibration check standard analyzed every 10 samples
- No recalibration, no extra rinsing




Variations within $\pm 10\%$ over 8 hours

Stability Established

Results: Stability – Internal Standards

- 8 hour analysis of soil samples
- Monitor internal standard intensities for all samples
 - Recoveries relative to calibration blank




Well within method-defined criteria over 8 hours

Results: MDLs & LLOQs

- Method Detection Limits
 - 7 replicate measurements of a blank
 - Multiply standard deviation by 3.14
- Lower Limit of Quantitation
 - Lowest concentration standard that can be read within 35% of its true value against the calibration curve

Results: MDLs & LLOQs



Summary

- Criteria of U.S. EPA Method 6020B met
 - Accuracy
 - Stability
 - Matrix effects
- Same methodology applied to both waters and soils
- Everything run in Collision mode
 - More confidence in results as samples can vary significantly