




Overcoming the Challenges of Trace Elemental Analysis in Environmental Samples Using ICP-MS

Daniel Kutscher, Fergus Keenan
Thermo Fisher Scientific


The world leader in serving science

Method Comparison for EPA Method 200.8 and SW846 Method 6020B

EPA 200.8




21 elements



18 common elements

EPA 6020B



23 elements

	Method 200.8 v. 5.4 (1994)	Method 6020B (2014)
Applications	Drinking and wastewater	Solid waste and soils and sediments
Method type	Method for regulatory compliance; QC acceptance criteria for lab accreditation	For guidance purpose, performance-based method
Approved Elements	21	23
Collision cell usage	Approved for wastewater, not for drinking water	Approved


Challenges

- Sample Matrix

- Instrumental drift
- Internal standard suppression higher than defined method limits

- Interferences

- Mathematical corrections may amplify errors
- False positive results if not resolved




Physical Interferences – ICP-MS Only

Dissolved Solids in Sample Matrix

- **Dissolved solids deposit onto interface cones**
 - Continued deposition can upset stability of instrument
 - Some sample matrix components deposit more than others
- **Too many ions in ion beam can cause spread of ion beam focusing**
 - Degrades recovery of internal standards
- **Solution – dilute the sample**
 - Dilute the sample to reduce the total matrix present
 - Use an internal standard to correct for changes in ionization, transport, etc.



Customized Interfaces for Maximum Matrix Tolerance




- Development research to model various skimmer cone base and tip temperatures
- Studies indicated:
 - Hotter skimmer cone tip leads to reduced deposition
 - Improved signal stability
 - Reduced cleaning/maintenance frequency

		
Robust 4.5 mm	High Matrix 3.5 mm	High Sensitivity 2.8 mm

Improved Matrix Tolerance

- Analysis of seawater samples diluted 1:10
 - 10 hours of continuous analysis
 - Data collected using standard and iCAP RQ skimmer cones for comparison




• Standard Skimmer Cone



• iCAP RQ Skimmer Cone

Handling Matrix – Customized Interfaces for Maximum Flexibility



Unique interface designed for maximum coverage of sensitivity and dynamic range:

			
High Sensitivity 2.8 mm	High Matrix 3.5 mm	Robust 4.5 mm	Argon Gas Dilution

STDS
KEDS
CCTS

STD
KED
CCT

STDR
KEDR
CCTR

STD AGD
KED AGD
CCT AGD

Easily implemented through dedicated grouping and autotunes in the software

He KED Performance – Precision and Sensitivity in High Matrix Samples

- Continuous high throughput analysis of brine samples
- Sample matrices contained 25% NaCl
- Sub-ppb DLs for 22 analytes and 5 internal standards
- Precision of all measurements $\leq 5\%$

Isotope	Detecton Limit (µg/L)	Measurement Precision (%RSD)
⁹ Be	0.56	5.3
⁵² Cr	0.61	2.2
⁵⁵ Mn	0.38	3.1
⁵⁶ Fe	0.42	2.5
⁵⁹ Co	0.12	2.0
⁶⁵ Cu	0.82	2.1
⁷⁵ As	0.25	1.9
⁷⁷ Se	0.55	5.3
⁸⁸ Sr	0.09	2.2
¹¹¹ Cd	0.05	2.4
¹³⁷ Ba	0.14	2.6
²⁰² Pb	0.14	4.2
²⁰⁸ Pb	0.04	1.8

Matrix Handling – Auto-dilution

- Supported by Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) Software

- Automated *prescriptive* dilution for preparation of:

- Samples
- Standards

- Automated *intelligent* dilution:

- Internal standard range auto-dilution
- Over calibration range auto-dilution



ESI
prepFAST



TELEDYNE CETAC SDX_{HPLD}

- Close coupling of the discrete sampling valves to ICP for minimized uptake and washout – increases throughput and reduces contamination

Matrix Handling - Intelligent Auto-dilution

• Calibration Over-Range

- User defines:
 - The maximum calibration limit as a percentage of the top calibration standard, e.g. 110%
 - The target concentration of the analyte after auto-dilution

Calibration Range

☒ Enable

Limit [%] 110

Target [%] 60

Action on Failure Abort LabBook

	Label	Status	Sample Type	Rack	Vial	prepFAST DF	Standard	Total Dilution Factor
1	Blank	●	BLK	1	1	1		
2	Level 1	●	STD	1	2	100	Tune B	
3	Level 2	●	STD	1	2	10	Tune B	
4	Level 3	●	STD	1	2	2	Tune B	
5	Over-range	●	UNKNOWN	1	2	1		
6	Over-range	●	UNKNOWN	1	2	3.367		
7	Washout	●	QC	1	1			1

Calibration line generated by auto-dilution from a stock solution in the autosampler rack

Automatically added Analysis

Automatically defined dilution factor

• Internal Standard Recovery

- User defines:
 - The range of acceptable Internal Standard Recovery, e.g. 60 - 125%
 - The auto-dilution factor step and # of possible dilutions

Internal Standard

☒ Enable

Upper Limit 125 [%] of Internal Standard Recovery

Lower Limit 60 [%] of Internal Standard Recovery

Autodilution Factor 5

Max. # of Autodilutions 2


Action on Failure Wash and Continue

No	Time	Sample Type	Label	89Y	115In	175Lu
2	6/17/2014 7:37:54 PM	BLK		100.0%	100.0%	100.0%
3	6/17/2014 7:39:10 PM	STD				
22	6/17/2014 8:00:57 PM	UNKNOWN	3% NaCl + 10 ppb Spike	13.5%	10.2%	7.1%
23	6/17/2014 8:01:46 PM	UNKNOWN	3% NaCl + 10 ppb Spike	76.0%	68.2%	62.8%
26	6/17/2014 8:04:32 PM	UNKNOWN	3% NaCl + 10 ppb Spike	90.9%	84.7%	85.2%

A sample of 3% NaCl is analyzed directly

The 3% NaCl sample is automatically diluted until the Internal Standard Recovery is within the defined limits


A Typical Day in the Busy Lab without Auto-dilution



Improvement Opportunities:

- Operator time wasted for simple tasks
- Risk of contamination
- Manual interaction may be error prone

A Typical Day in the Busy Lab with Auto-dilution



Action	Time	#	Total
Dilution	20s / sample	300	1.7 h
Preparation of Calibration/QC solutions	30 Minutes		0.5 h
Performance Verification	30 Minutes		0.5 h
Re-run failed samples	2 minutes	10%	1.0 h
Remaining Time			4.3 h



Improvement Opportunities:

- Operator time wasted for simple tasks
- Risk of contamination
- Manual interaction may be error prone

Gain: 3 h Operator time per day!

EPA Method 6020 B using SQ-ICP-MS




Parameter	Value
Forward Power	1550 W
Nebulizer	Burgener MiraMist®
Interface	Ni cones, Robust (4.5mm) insert
QCell Settings	KED using 100% He at 4.5 mL·min ⁻¹ , 3V KED voltage
Scan settings	0.01-0.3s

- Calibration

Standard #	Concentration [µg·L ⁻¹]	Elements
1	2	Ag, Al, As, Ba, Be, Ca, Cd, CO, Cr, Cu, Fe, K, Mg, Mn, Na, Ni, Pb, Sb, Se, Tl, V, Zn
2	100	
3	1,000	
4	5,000	
5	100,000	Ca, Fe, K, Mg, Na
6	40,000, 50	Al, Hg

- Sample – similar to NIST SRM 2781 Domestic Sludge
- QC Protocol as per EPA method 6020
 - CCV low and high, DUP, SER, LCS

Results – EPA Method 6020 B



- ✓ Internal Standard Recovery well within regulatory limits
- ✓ Consistent QC recovery
- ✓ Excellent long term stability over 12 h

Results – EPA Method 6020 B

	Sim. Sample conc.	Measured Sim. Sample	Recovery	Spiked conc.	Measured Spike	Recovery	IDL	MDL
	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$	%	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$	%	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$
⁹ Be	N/A	N/A		50	45	90	0.13	0.18
²⁴ Mg	20000	20619	103	50			0.77	1.3
²⁷ Al	50000	51604	103	50			0.44	1.5
³⁹ K	15000	16018	107	50			3.5	5.4
⁴⁴ Ca	120000	125986	105	50			8.3	7.4
⁵¹ V	250	295	118	50	59	119	0.02	0.02
⁵² Cr	500	463	93	50	41	81	0.03	0.04
⁵⁵ Mn	2000	2003	100	50	46	92	0.05	0.07
⁵⁷ Fe	85000	86863	102	50			0.36	1.6
⁵⁹ Co	N/A	N/A		50	59	118	0.003	0.01
⁶⁰ Ni	250	280	112	50	50	99	0.001	0.04
⁶³ Cu	2000	2003	100	50	47	93	0.07	0.12
⁶⁶ Zn	4000	3727	93	50	55	109	0.12	0.23
⁷⁵ As	25	29	116	50	58	117	0.01	0.04
⁸² Se	50	44	88	50	45	90	0.01	0.45
¹⁰⁷ Ag	300	257	86	50	43	87	0.02	0.3
¹¹¹ Cd	50	50	100	50	52	104	0.005	0.02
¹²¹ Sb	N/A	N/A		50	52	105	0.001	0.07
¹³⁷ Ba	2000	1824	91	50	60	119	0.03	0.06
²⁰² Hg	10	10	99				0.002	0.01
²⁰⁵ Tl	N/A	N/A		50	53	105	0.002	0.04
²⁰⁸ Pb	600	629	105	50	55	109	0.002	0.02

Handling Interferences – What Options are Available?

• But those aren't sufficient!

Some typical interferences in ICP-MS:



- Alternative Sample Introduction Techniques
 - E.g. aerosol desolvation to minimize the population of the precursor (parent) ions in the ICP ion source

NOT A UNIVERSAL SOLUTION

- Cold Plasma:
 - Reduce plasma power to reduce amount of Ar ionization


NOT A UNIVERSAL SOLUTION

- Mathematical Correction Equations:
 - Measure the isotope of interest, an un-interfered isotope, a polyatomic isotope and mathematically de-convolute to return an interference free value

AMPLIFIES ERRORS, DOESN'T ACCOUNT FOR UNKNOWN INTERFERENCES


Handling Interferences with Collision Reaction Cell Technology

- Thermo Scientific™ iCAP™ QCell Technology
- He Kinetic Energy Discrimination (KED)




Quadrupole set to filter out exact mass of target analyte

QCell in collision mode with pure He uses energy discrimination




Handling Interferences with Collision Reaction Cell Technology

- ...and Low Mass Cut Off




Quadrupole set to filter out exact mass of target analyte

QCell flatapole dynamically applies LMCO relative to target analyte




Interference Removal Collision Cell: Example ^{56}Fe

He KED
w/o CUT-OFF



Interference Removal Collision Cell: Example ^{56}Fe

He KED
with CUT-OFF



Handling Interferences – Low Mass Cut Off

- Low Mass Cut Off (LMCO)
 - Prevents the lower mass ions that contribute to interferences from traveling through the QCell...

Mass	Interferences	Precursors
⁵¹V	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, Cl
⁵⁶Fe	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁶³Cu	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar
⁷⁵As	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar

Mass	LMCO	Interferences	Precursors
⁵¹V	35	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, Cl
⁵⁶Fe	39	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁶³Cu	45	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar
⁷⁵As	47	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar

Handling Interferences – Low Mass Cut Off

- Low Mass Cut Off (LMCO)
 - Prevents the lower mass ions that contribute to interferences from traveling through the QCell...


Mass	Interferences	Precursors
⁵¹V	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, Cl
⁵⁶Fe	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁶³Cu	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar
⁷⁵As	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar

Mass	LMCO	Interferences	Precursors
⁵¹V	35	³⁵ Cl ¹⁶ O, ³⁷ Cl ¹⁴ N, ³⁴ S ¹⁶ OH	H, N, O, S, Cl
⁵⁶Fe	39	⁴⁰ Ar ¹⁶ O, ⁴⁰ Ca ¹⁶ O	O, Ar, Ca
⁶³Cu	45	⁴⁰ Ar ²³ Na, ¹² C ¹⁶ O ³⁵ Cl, ³¹ P ³² S	C, N, O, Na, P, S, Cl, Ar
⁷⁵As	47	⁴⁰ Ar ³⁵ Cl, ⁴⁰ Ca ³⁵ Cl, ⁴⁰ Ar ³⁴ SH, ³⁷ Cl ² H	H, S, Cl, Ca, Ar

- ...and reduces BECs even further than He KED alone


He KED Performance – Detection Limits in the Presence of Interferences

^{51}V KED



LOD = 7 ppt BEC = 12 ppt

^{75}As KED



LOD = 5 ppt BEC = 17 ppt

Sample Matrix: 2.5% HNO_3 + 1.5% HCl

EPA Method 6020 B using SQ-ICP-MS




Parameter	Value
Forward Power	1550 W
Nebulizer	Burgener MiraMist®
Interface	Ni cones, Robust (4.5mm) insert
QCell Settings	KED using 100% He at 4.5 mL·min ⁻¹ , 3V KED voltage
Scan settings	0.01-0.3s

- Calibration

Standard #	Concentration [µg·L ⁻¹]	Elements
1	2	Ag, Al, As, Ba, Be, Ca, Cd, CO, Cr, Cu, Fe, K, Mg, Mn, Na, Ni, Pb, Sb, Se, Tl, V, Zn
2	100	
3	1,000	
4	5,000	
5	100,000	Ca, Fe, K, Mg, Na
6	40,000, 50	Al, Hg

- Sample – similar to NIST SRM 2781 Domestic Sludge
- QC Protocol as per EPA method 6020
 - CCV low and high, DUP, SER, LCS

Results – EPA Method 6020 B





- ✓ Internal Standard Recovery well within regulatory limits
- ✓ Consistent QC recovery
- ✓ Excellent long term stability over 12 h


Results – EPA Method 6020 B

	Sim. Sample conc.	Measured Sim. Sample	Recovery	Spiked conc.	Measured Spike	Recovery	IDL	MDL
	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$	%	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$	%	$\mu\text{g}\cdot\text{L}^{-1}$	$\mu\text{g}\cdot\text{L}^{-1}$
⁹ Be	N/A	N/A		50	45	90	0.13	0.18
²⁴ Mg	20000	20619	103	50			0.77	1.3
²⁷ Al	50000	51604	103	50			0.44	1.5
³⁹ K	15000	16018	107	50			3.5	5.4
⁴⁴ Ca	120000	125986	105	50			8.3	7.4
⁵¹ V	250	295	118	50	59	119	0.02	0.02
⁵² Cr	500	463	93	50	41	81	0.03	0.04
⁵⁵ Mn	2000	2003	100	50	46	92	0.05	0.07
⁵⁷ Fe	85000	86863	102	50			0.36	1.6
⁵⁹ Co	N/A	N/A		50	59	118	0.003	0.01
⁶⁰ Ni	250	280	112	50	50	99	0.001	0.04
⁶³ Cu	2000	2003	100	50	47	93	0.07	0.12
⁶⁶ Zn	4000	3727	93	50	55	109	0.12	0.23
⁷⁵ As	25	29	116	50	58	117	0.01	0.04
⁸² Se	50	44	88	50	45	90	0.01	0.45
¹⁰⁷ Ag	300	257	86	50	43	87	0.02	0.3
¹¹¹ Cd	50	50	100	50	52	104	0.005	0.02
¹²¹ Sb	N/A	N/A		50	52	105	0.001	0.07
¹³⁷ Ba	2000	1824	91	50	60	119	0.03	0.06
²⁰² Hg	10	10	99				0.002	0.01
²⁰⁵ Tl	N/A	N/A		50	53	105	0.002	0.04
²⁰⁸ Pb	600	629	105	50	55	109	0.002	0.02


Determination of Arsenic and Selenium in Food Samples



Example: REE Mass Range for Apple Leaves



Accurate Analysis of Arsenic and Selenium in REE Matrix




Accurate Analysis of Arsenic and Selenium in REE Matrix

Single Quad ICP-MS: KED

Typically enhances M^{2+} Interferences

Solution:

Mass shift As and Se using O_2



Accurate Analysis of Arsenic and Selenium in REE Matrix


Single Quad ICP-MS: CCT O₂

Typically enhances M²⁺ Interferences

Solution:

Mass shift As and Se using O₂

Other interferences: ⁹¹Zr⁺, ^{94,96}Mo⁺, if present in the sample



Same Hardware Platform – Familiar Look and Feel

iCAP RQ ICP-MS

Innovative collision cell

Reaction Finder Software

Simplified power connections

Built-in safety for handling reactive gases



Intuitive user-friendly software

4 mass flow controllers:
He, O₂, H₂, NH₃

Quick connect and push-fit sample intro components

Bench-level easy-access interface

iCAP TQ ICP-MS

Compact footprint

Additional quadrupole for superior interference removal

Robust RF generator

Analysis with SQ and TQ in a single sample run




If you know how to use our single quadrupole, you already know how to use our triple quadrupole!

Accurate Analysis of Arsenic and Selenium in REE Matrix

- Solution is Triple Quadrupole ICP-MS:
 - Control ions entering Q2 using **Q1**

Type of Interference	^{75}As	SQ Mode	^{78}Se	SQ Mode
Polyatomic	$^{40}\text{Ar}^{35}\text{Cl}$	KED	$^{40}\text{Ar}^{38}\text{Ar}$	KED, H_2
	$^{40}\text{Ca}^{35}\text{Cl}$			
Isobaric	$^{150}\text{Nd}^{2+}$	O_2	$^{156}\text{Gd}^{2+}$	O_2
	$^{150}\text{Sm}^{2+}$		$^{156}\text{Dy}^{2+}$	

Q1 set to
analyte mass




Accurate Analysis of Arsenic and Selenium in REE Matrix

- Solution is Triple Quadrupole ICP-MS:
 - Control ions entering Q2 using **Q1**
 - Use O₂ to efficiently convert As and Se to AsO and SeO in **Q2**
 - REE⁺⁺ species don't react

Type of Interference	⁷⁵ As	SQ Mode	⁷⁸ Se	SQ Mode
Polyatomic	⁴⁰ Ar ³⁵ Cl	KED	⁴⁰ Ar ³⁸ Ar	KED, H ₂
	⁴⁰ Ca ³⁵ Cl			
Isobaric	¹⁵⁰ Nd ²⁺	O ₂	¹⁵⁶ Gd ²⁺	O ₂
	¹⁵⁰ Sm ²⁺		¹⁵⁶ Dy ²⁺	

Q2 filled with reactive gas (O₂)


Q1 set to analyte mass



Accurate Analysis of Arsenic and Selenium in REE Matrix

- Solution is Triple Quadrupole ICP-MS:
 - Control ions entering Q2 using **Q1**
 - Use O₂ to efficiently convert As and Se to AsO and SeO in **Q2**
 - REE⁺⁺ species don't react
 - Selectively detect AsO and SeO free from REE⁺⁺ interference, using **Q3**


Type of Interference	⁷⁵ As	SQ Mode	⁷⁸ Se	SQ Mode
Polyatomic	⁴⁰ Ar ³⁵ Cl	KED	⁴⁰ Ar ³⁸ Ar	KED, H ₂
	⁴⁰ Ca ³⁵ Cl			
Isobaric	¹⁵⁰ Nd ²⁺	O ₂	¹⁵⁶ Gd ²⁺	O ₂
	¹⁵⁰ Sm ²⁺		¹⁵⁶ Dy ²⁺	




Accurate Analysis of Arsenic and Selenium in REE Matrix

- Solution is Triple Quadrupole ICP-MS:
 - Control ions entering Q2 using **Q1**
 - Use O₂ to efficiently convert As and Se to AsO and SeO in **Q2**
 - REE⁺⁺ species don't react
 - Selectively detect AsO and SeO free from REE⁺⁺ interference, using **Q3**

Type of Interference	⁷⁵ As	TQ Mode	⁷⁸ Se	TQ Mode
Polyatomic	⁴⁰ Ar ³⁵ Cl	O ₂	⁴⁰ Ar ³⁸ Ar	O ₂
	⁴⁰ Ca ³⁵ Cl			
Isobaric	¹⁵⁰ Nd ²⁺		¹⁵⁶ Gd ²⁺	
	¹⁵⁰ Sm ²⁺		¹⁵⁶ Dy ²⁺	



Kick Interference to the Curb



Reaction Finder – eliminate the complexity of triple quadrupole ICP-MS analysis

Eliminate the Complexity of Triple Quadrupole ICP-MS

- Reaction Finder for Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ Software



Step 1: Select your element/s or isotope/s

Step 2: You're done!

- **Reaction Finder proposes the most appropriate gas/scan settings**
- **Settings for both single quad mode and triple quad mode are suggested, for reference**

Identifier	Q3 Analyte	SQ / TQ	CR Gas	Dwell time (s)	Channels	Spacing (u)	
▶ 78Se 78Se.160	78Se.160 (93.912 ▾)	TQ	O ₂	0.1	1	0.1	Normal
80Se 80Se.160	80Se.160	TQ	O ₂	0.1	1	0.1	Normal



Redefining triple quadrupole technology with unique ease of use

Simplest Method Development Using Reaction Finder

Without Reaction Finder

Select

- Select the Analytes to be measured

Select

- Select the internal standard element

Decide


- Are the suggested settings ok? If not, update them

Analyze

- Enter sample names and positions or import from LIMS and start the LabBook

★ Operator skills required

Analyte



Result

With Reaction Finder

Select

- Select the Analytes to be measured

Select

- Select the internal standard element

Decide


- Are the suggested settings ok? If not, update them

Analyze

- Enter sample names and positions or import from LIMS and start the LabBook

Less than 20 Minutes until a method is set up and the samples are ready to run!

Multi-element results: River water reference material (NRC SLRS-5)




Analyte	Measurement Mode	Measured $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
^{23}Na	KED	5085	5380	95
^{24}Mg	KED	2665	2540	105
^{27}Al	KED	55.3	49.5	112
^{39}K	KED	863	839	103
^{56}Fe	KED	93.2	91.2	102
^{59}Co	KED	0.05	0.05	107
^{60}Ni	KED	0.52	0.48	110
^{63}Cu	KED	18.2	17.4	105
^{75}As	TQ-O2	0.43	0.41	104
^{78}Se	TQ-O2	0.10	N/A	N/A
^{208}Pb	KED	0.08	0.081	99
^{238}U	KED	0.10	0.09	109

Multi-element results: River water reference material (CRM 1643F)

Analyte	Measurement Mode	Measured $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
51V	TQ-O2	36.03	36.07 ± 0.28	100
52Cr	TQ-O2	18.06	18.50 ± 0.10	98
55Mn	SQ-KED	39.50	37.14 ± 0.60	106
57Fe	SQ-KED	106.11	93.44 ± 0.78	114
58Ni	SQ-KED	55.97	59.8 ± 1.4	94
59Co	TQ-O2	24.22	25.30 ± 0.17	96
63Cu	SQ-KED	21.32	21.66 ± 0.71	98
66Zn	TQ-O2	87.78	74.4 ± 1.7	118
75As	TQ-O2	57.10	57.42 ± 0.38	99
80Se	TQ-O2	11.97	11.70 ± 0.081	102
98Mo	TQ-O2	123.8	115.3 ± 1.7	107
107Ag	SQ-KED	0.94	0.97 ± 0.0055	97
111Cd	SQ-KED	5.91	5.89 ± 0.13	100
121Sb	SQ-KED	54.59	55.45 ± 0.40	98
125Te	SQ-KED	0.95	0.977 ± 0.0084	97
205Tl	SQ-KED	6.38	6.892 ± 0.035	93
208Pb	SQ-KED	18.19	18.488 ± 0.084	98

Multi-element results: Sea Water reference material (NRC CASS-5)



Analyte	Measurement Mode	Measured $\mu\text{g}\cdot\text{L}^{-1}$	Certified Value $\mu\text{g}\cdot\text{L}^{-1}$	Recovery %
51V	TQ-O2	1.37	1.32 ± 0.14	104
52Cr	TQ-O2	0.069	0.106 ± 0.013	65
55Mn	SQ-KED	2.56	2.62 ± 0.20	98
57Fe	TQ-O2	1.49	1.44 ± 0.11	103
63Cu	SQ-KED	0.378	0.380 ± 0.028	99
75As	TQ-O2	1.17	1.24 ± 0.09	94
98Mo	TQ-O2	10.80	9.82 ± 0.72	110
111Cd	SQ-KED	0.0227	0.0215 ± 0.0018	106

Summary and Conclusion

- Productivity is key for most laboratories
 - Matrix tolerance allows to analyze more samples without internal standards going out of limits
 - Productivity tools may even free up more time from the operators

- Sometimes there is more than meets the eye!
 - Advanced interference removal combined with modern sample introduction systems allow high throughput analysis
 - Modern software makes the technology accessible for all laboratories
 - Right first time results straight away accessible in almost any matrix