


# A Direct Injection Approach for Analysis of Legacy and Emerging Perfluoroalkyl Substances in Environmental Water and Soil Samples

Kari Organtini  
Waters Corporation  
August 8, 2019

## Perfluoroalkylated Substances (PFAS)

- PFAS = PFC = AFFF
- First **created** in the **1930s**
- **Widespread applications**
  - Non stick **coatings**, **surfactants**, food **packaging**, firefighting **foams**
  - **Polymerization aid** for polytetrafluoroethylene (PTFE) and other fluoropolymers – how PFOS and PFOA became famous
- **Stable** and **persistent** in the environment (POP)
  - Bio-accumulative
- **Identified** in environmental samples **worldwide**
  - Found in arctic **polar bears**
  - Most humans have a range of **PFAS** in their **blood**




## Perfluoroalkylated Substances (PFAS)

- The **most-studied** PFAS chemicals are PFOA and PFOS.
- Studies in laboratory animals indicate that PFOA and PFOS **can cause effects** in
  - Reproductive and developmental systems
  - Liver and kidney
  - Immune system in laboratory animals
  - Tumors
- But there are **thousands** of **PFAS** that have been created...



Perfluoro Sulfonic Acid  
(PFOS)



Perfluoro Carboxylic Acid  
(PFOA)



Thousands...

Waters


THE SCIENCE OF WHAT'S POSSIBLE.™



Perfluoro Carboxylic Acid  
(PFOA)



Perfluoro Sulfonic Acid  
(PFOS)




Sulfonamidoacetic  
acid  
(N-EtFOSAA)



Perfluoro Telomer Acid  
(FHEA or 6:2 FTA)




Perfluoro Telomer Sulfonate  
(6:2 FTS)




Emerging PFAS  
(GenX)

Cyclic  
(PFecHS)



## How Are People Exposed to PFAS?



## A Word on Regulations...

- Currently there are **no federal level regulations** requiring monitoring **for PFAS** (worldwide)
- There are **advisory limits** and established testing methods
  - EPA drinking water series – UCMR 3
  - European water framework directive/drinking water directive
  - ISO and ASTM methods
- EPA advisory limit currently set at **70 ng/L (ppt)** for total PFOA/PFOS
- Individual state proposals and regulations (18 states)
- EPA released **PFAS Action Plan** in early 2019



## A Word on Regulations...

- Currently there are **no federal level regulations** requiring monitoring for **PFAS** (e)
- There are **advisory limits** for testing methods
  - EPA drinking water
  - EPA surface water
  - Is
- EPA **action level** set at **70 ng/L (ppt)** for total P
- Individual **proposals and regulations** (18 states)
- EPA released **PFAS Action Plan** in early 2019

STAY  
TUNED!



## Global Interest in PFAS

NATIONAL QUEENSLAND

**Chemicals in Brisbane Airport spill have sparked fears 'worldwide'**

**Report says incidents of GenX, other PFAS, nearly doubled in last year**

**White House, EPA headed off chemical pollution study**

The intervention by Scott Pruitt's aides came after one White House official warned the findings would cause a 'public relations nightmare.'

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



**NEW TEFLON TOXIN FOUND IN NORTH CAROLINA DRINKING WATER**

**EU project assesses 'critical' PFAS use in textiles**

**'New' perfluoroalkyl substances found in fish in China**

**PFAS costs Europe more than €50 billion a year in health problems**



## Options for the Analysis of Water Samples

### SPE enrichment prior to injection

- Sample prep allows for use of mid-level sensitivity for mass spectrometer



Xevo TQ-S micro

### Large Volume Injection

- High sensitivity mass spectrometer required



Xevo TQ-XS

## Comparison of PFAS water analysis methods

	<b>EPA Method 537*</b>	<b>ISO 25101</b>	<b>ASTM 7979</b>
<b>Sample Prep</b>	SPE (SDVB)	SPE	Dilute and Filter
<b>Injection Volume</b>	1 - 10 µL	1 - 10 µL	30 µL
<b>Instrument sensitivity required</b>	Low to mid	Low to mid	High
<b>Number of compounds</b>	18	2	24


\*EPA 537 is for drinking water analysis only



# PFAS Contamination




## PFAS's Interferences/Contamination






## The Waters PFAS Kit

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.™



## MRM Chromatograms of PFOA With PFAS Kit Installed



# **Direct Large Volume Injection for**

## **Water (ASTM 7979)**

**and**

## **Soil (ASTM 7968)**

# Sample Pre- Treatment

**WATER**



5 mL sample +  
5 mL Methanol

**SOIL**



2 g sample (soil) + 10 mL 1:1  
Water:Methanol  
Adjust pH to ~ 9-10  
Shake 1 hour

Centrifuge

Water  
Surface  
Ground  
Influent  
Effluent



Syringe filter entire  
sample



Acidify and transfer to  
polypropylene vial

Soils  
Sand  
Silt  
Lean Clay  
Fat Clay

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.™



## Instrument Methods

### Source Parameters

- Instrument: Xevo TQ-XS
- Ion Mode: ESI-
- Capillary Voltage: 1.0 kV
- Desolvation Temperature: 500° C
- Desolvation Flow: 1100 L/hr
- Cone Flow: 150 L/hr

### MS Method

- Developed using QuanOptimize
  - MRMs, CV, CE
- Divert flow to waste from 15 – 21 mins


### LC Method

- Instrument: Acquity I Class with PFAS Kit
- Column: CSH Phenyl Hexyl 2.1mm x 100 mm, 1.7 µm
- Mobile Phase A: 95:5 H<sub>2</sub>O:MeOH + 2 mM ammonium acetate
- Mobile Phase B: MeOH + 2 mM ammonium acetate
- Injection Volume: 30 µL
- Gradient:

Time (min)	Flow (mL/min)	%A	%B
0	0.3	100	0
1	0.3	80	20
6	0.3	55	45
13	0.3	20	80
14	0.4	5	95
17	0.4	5	95
18	0.3	100	0
22	0.3	100	0


## Compounds Included in Methods

### Carboxylates



C4 – C14, C16, C18

### Sulfonates



C4 – C10

### Emerging



GenX	ADONA
11CI-PF3OUdS	9CI-PF3ONS
PFEESA	NFHDA
PFMBA	

### Others

4:2/6:2/8:2 FTS	FHEA/FOEA/FDEA
FOSA	PFecHS
FHUEA	FHpPA
FOUEA	diPAP
NMeFOSAA/NEtFOSAA	
NMeFOSA/NEtFOSA	

## Overall Method Summary

ERA **certified** QC sample injected as **control** throughout analysis



## Method Detection Limits (MDL) for Water

★ Well **below** the necessary **reporting limits** for the compounds defined in **ASTM 7979** ★

Compound	Sample spike (ng/L)	MDL (ng/L)	Reporting range (ng/L)*	R <sup>2</sup>	Compound	Sample spike (ng/L)	MDL (ng/L)	Reporting range (ng/L)*	R <sup>2</sup>
PFBA	100	25.20	50–2000	0.993	FOSA	10	1.29	10–400	0.999
PFPeA	10	1.04	50–2000	0.999	N-Et-FOSAA	10	1.90	10–400	0.997
PFHxA	10	1.33	10–400	0.999	N-Me-FOSAA	10	1.59	10–400	0.999
PFHpA	10	0.91	10–400	0.999	N-Et-FOSA	10	1.45	–	0.997
PFOA	10	1.42	10–400	0.999	N-Me-FOSA	10	1.19	–	0.999
PFNA	10	1.32	10–400	0.999	FHUEA	10	1.53	10–400	0.999
PFDA	10	0.84	10–400	0.998	FOUEA	10	1.36	–	0.999
PFUnDA	10	2.52	10–400	0.996	8:2 diPAP	300	50.16	–	0.988
PFDoDA	10	1.76	10–400	0.993	4:2 FTS	10	1.50	10–400	0.999
PFTriDA	10	2.34	10–400	0.991	6:2 FTS	10	N/A	10–400	0.999
PFTreDA	10	1.99	10–400	0.993	8:2 FTS	10	2.62	10–400	0.997
PFHxDA	200	25.41	–	0.984	PFecHS	10	1.17	10–400	0.998
PFOcDA	400	41.99	–	0.983	FHEA	200	42.19	300–8000	0.994
PFBS	10	1.21	10–400	0.999	FOEA	200	50.38	200–8000	0.997
PFPeS	10	1.07	10–400	0.999	FDEA	200	79.48	200–8000	0.993
PFHxS	10	1.41	10–400	0.999	FHpPA	10	1.47	10–400	0.999
PFHpS	10	1.57	10–400	0.999	ADONA	10	0.82	–	0.999
PFOS	10	1.61	10–400	0.999	9CI-PF3ONS	10	1.06	–	0.999
PFNS	10	1.67	10–400	0.999	11CI-PF3OUdS	10	1.45	–	0.998
PFDS	10	1.44	10–400	0.997					

Reporting Range is the required range the method must be able to cover as defined in ASTM 7979 method









## Lower Limits of Quantitation (LLOQ) for Soil

Compound	LLOQ (ng/L)	LLOQ (ng/kg)
	In vial	In sample
PFBA	25	125
PFPeA	< 1	< 5
PFHxA	< 1	< 5
PFHpA	< 1	< 5
PFOA	< 1	< 5
PFNA	1	5
PFDA	< 1	< 5
PFUnDA	5	25
PFDoDA	5	25
PFTriDA	5	25
PFTreDA	5	25
PFBS	< 1	< 5
PFPeS	< 1	< 5
PFHxS	1	5
PFHpS	1	5


Compound	LLOQ (ng/L)	LLOQ (ng/kg)
	In vial	In sample
PFOS	1	5
PFNS	5	25
PFDS	5	25
FOSA	< 1	< 5
N-Et-FOSAA	5	25
N-Me-FOSAA	5	25
4:2 FTS	< 1	< 5
6:2 FTS	10	50
8:2 FTS	5	25
ADONA	< 1	< 5
9CI-PF3ONS	5	25
11CI-PF3OUdS	5	25
GenX	10	50
PFMBA	< 1	< 5
PFEESA	< 1	< 5

# Linearity and Sensitivity

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.™




# Recovery in Water



## Recovery in Soils


Percent Recovery in Soil





## Instrument Robustness

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.™




20 injections of  
surface water  
sample

Overall RSD: <  
10 %

# Spiked Surface Water

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™




PFBA and PFPeA at 300 ng/L; 4:2, 6:2, and 8:2 FTS at 1200 ng/L; all other compounds at 60 ng/L. (\*) this compound shown off scale.

# Spiked Lean Clay Sample

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



(\*) compound shown off scale (\*\*) compound shown zoomed

## Conclusions

- PFAS are ubiquitous environmental contaminants that are detected around the globe
- Direct Injection approach can be applied to both water and soil samples
- Both ASTM 7979 (water) and ASTM 7968 (soil) are suitable for the basic range of PFAS normally monitored
- Limited sample preparation allows for higher throughput of samples
- New emerging PFAS can easily be incorporated into these methods



## Acknowledgements

- EPA
  - Larry Zintek
  - Mark Strynar
  - Raj Singhvi
- Waters
  - Ken Rosnack
  - Lauren Mullin
  - Doug Stevens
  - Euan Ross
  - Simon Hird
  - Gareth Cleland



# Waters

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.™