

# A New Method for Measuring Ammonia and TKN in Wastewater


*FIA*lab<sup>®</sup>

*Fluidics Intelligently Automated*



## Table of Contents

- Current automated ammonia & TKN methods
- Why a new method?
- Operation
- Performance
- Approval
- Conclusions and questions





# HELLO!

**I am Dr. Ilkka Lahdesmaki**


Chief Scientist

You can find me at [ilkka@flowinjection.com](mailto:ilkka@flowinjection.com)

# 1

## CURRENT METHODS

Automated methods for  
ammonia and TKN






## CURRENT METHODS

**FIAlab**

Fluidics Intelligently Automated


- Flow methods (FIA, SFA) achieve a high degree of automation





## CURRENT METHODS

### ■ Use of gas diffusion to replace distillation



### ■ Detection of diffusion-separated ammonia

- ▷ Phenolate / salicylate + photometric detection (EPA 351.2, SM 4500-N<sub>org</sub> D.)
- ▷ Conductimetric detection (EPA Timberline-Ammonia-001)

# 2

## NEW METHOD

What is the new method?


Why would a new method be needed?



## NEW METHOD – HOW DOES IT WORK?

Based on:

- ▶ o-phthalaldehyde indicator
- ▶ fluorescence detection



*Felix et al.*  
*J. Braz. Chem. Soc.*, 23 (1), 142-147, 2012






# NEW METHOD – HOW DOES IT WORK?

**FIAlab**

Fluidics Intelligently Automated

**PUMP**  
Speed: 30%

**INJECTION VALVE**  
Sample loop: 35  $\mu$ L





## NEW METHOD – WHY INTRODUCE IT?

### ■ Extensive prior use in oceanography

#### 8.7 DETERMINATION OF AMMONIUM

##### 8.7.1 OVERALL DESCRIPTION OF THE METHOD

The fluorometric method offers many advantages over indophenol blue colorimetry: simplicity, very high sensitivity, stable reagents with low toxicity, no refractive index blank, mostly insignificant fluorescence blank from natural substances, almost unaltered by sample turbidity, and a low salt effect. Both methods produce equivalent results (Kérouel and Aminot, 1997). While a detection

*A. Aminot et al. in O. Wurl (ed.) "Practical Guidelines for the Analysis of Seawater", p. 166-169. CRC Press, Boca Raton, 2009.*



## NEW METHOD – WHY INTRODUCE IT?

■ Used by customers for **ammonia** and **TKN** measurements

- ▷ Oceanography research groups
- ▷ International environmental laboratories
- ▷ Very positive feedback: no phenol, more sensitive, stable baseline, **consistent**



■ Bring method available to accredited labs in the U.S.

- ▷ EPA Alternative Test Procedure (ATP) application
- ▷ Validation for use in wastewater analysis



## NEW METHOD – APPROVAL PROCESS

### ■ Recommendation letter from EPA – 10 Apr 2018

I have reviewed FIALab Method 100 (ATP Case No. N15-0001), “Determination of Inorganic Ammonia by Continuous Flow Gas Diffusion and Fluorescence Detector Analysis”, and the supporting validation data in ATP Case No. N15-0001. I determined that this method meets all requirements for measurement of ammonia in wastewater and total Kjeldahl nitrogen (TKN) digestates. That is, the performance of this method is substantially similar to methods listed at 40 CFR Part 136 for measurement of ammonia in wastewater and TKN digestates.

### ■ Inclusion in 2019 Methods Update Rule proposal – 11 Jun 2019

contained as supporting documents within the docket for this proposed rule. These proposed new methods include: FIALab Method 100, “Determination of Inorganic Ammonia by Continuous Flow Gas Diffusion and Fluorescence Detector Analysis,” MACHEREY-NAGEL GmbH and

# 3


## PRACTICAL USE

How to implement on FIA instrument?

Performance metrics



## PRACTICAL USE – INSTRUMENT CONFIGURATION



Pictured:

**FIAlyzer-1000  
with PMT Detector**




## PRACTICAL USE – DATA OUTPUT


**FIAlab**

Fluidics Intelligently Automated

Calibration 50 – 10,000  $\mu\text{g N / L}$



Single injection 50  $\mu\text{g N / L}$






# PRACTICAL USE – DATA OUTPUT

**FIALab**

Fluidics Intelligently Automated



Fit model: 2nd order polynomial

Weighting: none

Coeff A: 2526.62741

Coeff B: 19.53941

Coeff C: -0.00018

R<sup>2</sup>: 0.99995


RSE (%): 4.53106

☐ Apply Drift Correction

	Name	Peak Response	Known Concentration	Calculated Concentration	% Error	Enabled
▶	Blank	2464.9000	0	-3.159	NA	<input checked="" type="checkbox"/>
	Std 1	3442.4000	50	46.8881	6.22	<input checked="" type="checkbox"/>
	Std 2	4376.8000	100	94.7716	5.23	<input checked="" type="checkbox"/>
	Std 3	11951.6000	500	484.5086	3.10	<input checked="" type="checkbox"/>
	Std 4	21934.5000	1000	1002.4788	0.25	<input checked="" type="checkbox"/>
	Std 5	41811.5000	2000	2049.0256	2.45	<input checked="" type="checkbox"/>
	Std 6					<input checked="" type="checkbox"/>
	Std 7					<input checked="" type="checkbox"/>

Relative Standard Error

2nd order fit 50 – 10,000 µg N / L



Fit model: 1st order polynomial

Weighting: 1/response<sup>1/2</sup>

Coeff A: 2464.87341

Coeff B: 19.52027

Coeff C: 0

R<sup>2</sup>: 0.99991

RSE (%): 2.06344

☐ Apply Drift Correction

	Name	Peak Response	Known Concentration	Calculated Concentration	% Error	Enabled
▶	Blank	2464.9000	0	0.0014	NA	<input checked="" type="checkbox"/>
	Std 1	3442.4000	50	50.0775	0.16	<input checked="" type="checkbox"/>
	Std 2	4376.8000	100	97.9457	2.05	<input checked="" type="checkbox"/>
	Std 3	11951.6000	500	485.9936	2.80	<input checked="" type="checkbox"/>
	Std 4	21934.5000	1000	997.4056	0.26	<input checked="" type="checkbox"/>
	Std 5	41811.5000	2000	2015.6804	0.78	<input checked="" type="checkbox"/>
	Std 6	95201.4000	5000	4750.7809	4.98	<input type="checkbox"/>
	Std 7	180112.1000	10000	9100.6541	8.99	<input type="checkbox"/>

1st order fit 50 – 2,000 µg N / L





## PRACTICAL USE – PERFORMANCE METRICS

**FIAlab**


Fluidics Intelligently Automated

Detection Limit	12 µg N/L	
Reporting Limit	50 µg N/L	
Range Upper Limit	10 000 µg N/L	
Spike recovery	99.3%	POTW (Anaerobic digester sludge)
	97.3%	Industrial discharge (Food process)
	109%	Industrial discharge (Metal finish)
	96.9%	River water
	102%	POTW (Final effluent, pre-UV)
	105%	POTW (Primary clarifier effluent)
Throughput	50 samples / h	



## PRACTICAL USE – PERFORMANCE METRICS

- Above performance with 35  $\mu$ L sample loop
- ▶ "Standard" size that works both for ammonia and TKN
- ▶ Ability to adjust range by varying sample loop size



# 4

## CONCLUSIONS

What you should take home from this talk






## CONCLUSIONS

### ■ New method for measuring **ammonia** and **TKN**

- ▷ Applicable to wastewater and natural waters
- ▷ Included in 2019 MUR, signed and being proposed
- ▷ Pending inclusion in 40 CFR Part 136

### ■ Method highlights

- ▷ High sensitivity
- ▷ Quick calibration and analysis
- ▷ Robust – gas diffusion
- ▷ Simple - one-reagent chemistry



**FIAlyzer-1000**



## CREDITS

- Askew Scientific Consulting
  - ▷ Dr. Ed Askew
- EPA Clean Water Act ATP Coordination Team
  - ▷ Lem Walker
  - ▷ Kevin Roberts
- Laboratories participating in the ATP validation experiments



**Askew**  
Scientific Consulting LLC





# THANKS!

**Any questions?**

You can find me at  
[sales@flowinjection.com](mailto:sales@flowinjection.com)