

Moving beyond monitoring legacy per and polyfluoroalkyl substances (PFAS)

Screening strategies for the growing list.

James Pyke, Jerry Zweigenbaum, Tarun Anumol
Agilent Technologies, Inc.



PFAS related products

Common household products and industrial uses

Non-stick
surfaces



Grease-proof
food packaging



Surfactants and
lubricants

Fire-fighting
foams



Stain guards

Water repellents



PFAS related products

Common laboratory materials

Solvent caps,
filters and
tubing

Caps on
samples vials

Pump seals,
frits, degasser
materials




Air conditioning
filters

Gloves/Coats

Sample
preparation
consumables


PFAS

Perfluoroalkyl substance



ALL H atoms linked to C in alkyl chain are substituted with F

Polyfluoroalkyl substance



SOME (but not all) H atoms linked to C alkyl chain are substituted with F

- Thermal & Chemical stability: grease-proof food packaging, stain repellents
- Zwitterionic properties: surfactants
- Surface-tension lowering: fire-fighting foams



PFAS

Concern Grows Over Tainted Drinking Water

Vermont, New Hampshire and New York expand efforts to find out how much of a potentially toxic chemical is in drinking water

England since August 2014, when a resident of Hoosick Falls, N.Y., near the Vermont border, tested his drinking water and found high levels of the acid. The man was concerned because his father, a former employee of the town's plastics plant that used PFOA, died of cancer. The state in March sampled PFOA levels up to 620 parts per trillion in private Litchfield wells, well above the 100-parts-per-trillion level at which New Hampshire officials start to consider the amount unsafe. Tests in Merrimack measured as high as 1,600 parts per trillion.

environment

Williamstown water contamination highlights dangers of PFOS and PFOA

Qantas faces \$180,000 fine over toxic foam spill at Brisbane Airport

21 April 2017 - 04:56pm - First published 21 April 2017 - 11:02am
By Ruth McCosker

f t e A A A

Qantas has been hit with an investigation notice while residents continue to be warned not to eat seafood following a toxic spill into Brisbane's waterways.

On Monday, April 10, 22,000 litres of a firefighting foam containing perfluorinated compounds was spilled from a Qantas hangar at Brisbane Airport.



Volume 94 Issue 20 | pp. 20-22
Issue Date: May 16, 2016 | Web Date: May 11, 2016

Perfluorinated chemicals taint drinking water

Are there toxins in your fast food packaging?



QLD News


Brisbane River seafood warning still in effect

Chris Honnery, The Courier-Mail
April 20, 2017 8:49pm

Regulated PFAS monitoring

Common target compounds

Compound	Formula	Abbreviation
Perfluoroalkylcarboxylic acids (PFCAs)		
Perfluoro-n-butyric acid	C ₄ HF ₇ O ₂	PFBA
Perfluoro-n-pentanoic acid	C ₅ HF ₉ O ₂	PFPeA
Perfluoro-n-hexanoic acid	C ₆ HF ₁₁ O ₂	PFHxA
Perfluoro-n-heptanoic acid	C ₇ HF ₁₃ O ₂	PFHpA
Perfluoro-n-octanoic acid	C ₈ HF ₁₅ O ₂	PFOA
Perfluoro-n-nonanoic acid	C ₉ HF ₁₇ O ₂	PFNA
Perfluoro-n-decanoic acid	C ₁₀ HF ₁₉ O ₂	PFDA
Perfluoro-n-undecanoic acid	C ₁₁ HF ₂₁ O ₂	PFUdA
Perfluoro-n-dodecanoic acid	C ₁₂ HF ₂₃ O ₂	PFDoA
Perfluoro-n-tridecanoic acid	C ₁₃ HF ₂₅ O ₂	PFTTrDA
Perfluoro-n-tetradecanoic acid	C ₁₄ HF ₂₇ O ₂	PFTeDA
Perfluoro-n-pentadecanoic acid	C ₁₅ HF ₂₉ O ₂	PFPeDA
Perfluorinated sulfonates (PFSA)		
Perfluoro-1-butanedisulfonate	C ₄ F ₉ SO ₃ ⁻	PFBS
Perfluoro-1-hexanedisulfonate	C ₆ F ₁₃ SO ₃ ⁻	PFHxS
Perfluoro-1-octanedisulfonate	C ₈ F ₁₇ SO ₃ ⁻	PFOS
Perfluoro-1-decanedisulfonate	C ₁₀ F ₂₁ SO ₃ ⁻	PFDS
Perfluorinated sulfonamides (FOSA)		
Perfluoro-1-octansulfonamide	C ₈ H ₂ F ₁₇ NO ₂ S	FOSA
N-Methylperfluoro-1-octansulfonamide	C ₉ H ₄ F ₁₇ NO ₂ S	N-MeFOSA
N-Ethylperfluoro-1-octansulfonamide	C ₁₀ H ₆ F ₁₇ NO ₂ S	N-EtFOSA
Perfluorinated sulfonamidoethanols (FOSE)		
2-(N-methylperfluoro-1-octansulfonamido)-ethanol	C ₁₁ H ₈ F ₁₇ NO ₃ S	N-MeFOSE
2-(N-ethylperfluoro-1-octansulfonamido)-ethanol	C ₁₂ H ₁₀ F ₁₇ NO ₃ S	N-EtFOSE
Fluorinated Telomer Sulfonates (FTS)		
1H,1H,2H,2H-Perfluorohexanesulfonic acid	C ₆ H ₅ F ₉ SO ₃	4:2 FTS
1H,1H,2H,2H-perfluorooctane sulfonate	C ₈ H ₅ F ₁₃ SO ₃	6:2 FTS
1H,1H,2H,2H-Perfluorodecanesulfonic acid	C ₁₀ H ₅ F ₁₇ SO ₃	8:2 FTS
Perfluorinated sulfonamidoacetic acids (FOSAA)		
Perfluorooctane sulfonamidoacetic acid	C ₁₀ H ₄ F ₁₇ NO ₄ S	FOSAA
N-methylperfluorooctane sulfonamidoacetic acid	C ₁₁ H ₆ F ₁₇ NO ₄ S	N-MeFOSAA
N-ethylperfluorooctane sulfonamidoacetic acid	C ₁₂ H ₈ F ₁₇ NO ₄ S	N-EtFOSAA




>4000 PFAS compounds in commerce

Common Acronyms

PFCa	Perfluoroalkylcarboxylic acid
PFOA	Perfluorooctanecarboxylic acid
PFAS	Perfluoroalkylsulfonate
PFOS	Perfluorooctanesulfonate
PFASi	Perfluoroalkylsulfinate
FOSA	Perfluorooctanesulfonamide
FOSAA	Perfluorooctanesulfonamidoacetic acid
FOSE	Perfluorooctanesulfonamidoethanol
FTOH	Fluorinated telomer alcohol (-OH functional group)
FTA	Fluorinated telomer acid
FTUA	Fluorinated telomer unsaturated acid
FTS	Fluorinated telomer sulfonate
PFAPA	Perfluoroalkylphosphonic acid
PFPI	Perfluoroalkylphosphinate
PAP	Mono-substituted polyfluoroalkylphosphate ester
diPAP	Di-substituted polyfluoroalkylphosphate ester
PFAl	Perfluoroalkyl iodide
SFA	Semifluorinated alkane
FTI	Fluorinated telomer iodide
FTO	Fluorinated telomer olefin
FTAC	Fluorinated telomer acrylate

https://www.well-labs.com/docs/pfc_reference_handling_guide.pdf



Wang, Z et al. (2017). Environ. Sci. Technol. 51, 2508-2518.

Aim

Simplified and simultaneous Target quantitation and Suspect screening

MassHunter Acquisition software with Study Manager


MassHunter Quantitation software

Data acquisition

Target Quantitation


Targeted screen

Report results




Sampling for screening

Traditional SPE is selective





Sampling for screening

Move to direct injection to broaden the screen



Agilent 1290 Infinity II UHPLC conditions



High Speed Pump (G7120A)

Solvent A: 5mM ammonium acetate in water
Solvent B: Neat acetonitril
Flow Rate: 0.4 mL/min
Max Pressure Limit: 1300 bar

Gradient:

Time (min)	%B
0.00	3.00
1.00	25.00
9.00	85.00
10.00	97.00
12.00	97.00
14.00	3.00
15.00	3.00

Stop time: 15.00 min

Multisampler (G7167B)

Injection Volume: 30 μ L
Multiwash: Seat back flush and needle wash with 5s each of 100% isopropanol, then 100% acetonitrile, then 100% water

Multicolumn Thermostat (G7116B)

Column Temperature: 30°C
Column: Agilent InfinityLab Poroshell HPH-C18
2.1 $\times$ 100 mm, 1.9 μ m (p/n 695675-702)

Agilent 6546 LC/Q-TOF LC/MS system (G6546A) conditions

Agilent Jet-Stream Ion Source

Drying Gas Temp: 320 °C
Drying Gas Flow: 8 L/min
Nebulizer: 35 psi
Sheath Gas Temp: 350 °C
Sheath Gas Flow: 11 L/min
Capillary Voltage: 3500 V
Nozzle Voltage: 0 V
Fragmentor Voltage: 115 V



Tune mode

Ion Polarity: Negative
Mass Range: Low (1700m/z)
Slicer mode: High resolution

Acquisition mode


50-1100 m/z
Rate: 6 spectra/sec
Collision Energy: 0, 10, 20 V

Reference Mass Correction: Enabled using bottle A
119.03632 (M-H)⁻ adduct of purine
980.016375 (M+Ac)⁻ adduct of HP-0921




Quantitative analysis results

Separation of the target compounds




Quantitative analysis results


Lower Limit of Detection



PFBA @ 200 ng/L




PFOS @ 5 ng/L




PFOA @ 20 ng/L

Quantitative analysis results


Dynamic range



PFBA 200 – 20000 ng/L



PFOS 5 – 2000 ng/L



PFOA 20 – 2000 ng/L

Monitoring suspect PFAS

PFAS/EPA: ToxCast Chemical Inventory

Chemistry Dashboard | EPA/PAFAS

United States Environmental Protection Agency

PFAS|EPA: ToxCast Chemical Inventory

Search EPA/PAFAS Chemicals

Identifier substring search

List Details

Description: Per- and Polyfluoroalkyl Substances (PFAS) included in EPA's expanded ToxCast chemical inventory and available for testing. These PFAS chemicals were successfully procured from commercial suppliers (with a small number provided by National Toxicology Program partners) and were deemed suitable for testing (i.e., solubilized in DMSO above 5mM, and not gaseous or highly reactive). All or portions of this inventory are being made available to EPA researchers and collaborators to be analyzed and tested in various high-throughput screening (HTS) and high-throughput toxicity (HTT) assays.

The https://comptox.epa.gov/dashboard/chemical_lists/EPAPFAS751 list is a prioritized subset of this larger chemical inventory.

The https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASINSOL list were chemicals procured, but found to be insoluble in DMSO above 5mM.

Number of Chemicals: 430

430 chemicals

Select all Download Send to Batch Search Default TOXICO CASRN TOXCAST Hide chemicals that are: Order by Name

2H-Perfluoro-2-propanol
DTXSID: DTXSID0021134
CASRN: 920-66-1
TOXCAST: -

Perfluorooctanesulfonyl fluoride
DTXSID: DTXSID0027140
CASRN: 307-35-7
TOXCAST: -

N-Ethyl-N-(2-hydroxyethyl)perfluorooct...
DTXSID: DTXSID0027426
CASRN: 1691-99-2
TOXCAST: -

How to rapidly screen for PFAS with more confidence, without standards?


- Focus search criteria:
 - on mass accuracy
 - on expected retention time
- Use all data possible to make a putative identification

https://comptox.epa.gov/dashboard/chemical_lists/EPAPFASINV

Monitoring suspect PFAS

Adding focus through Retention Time (RT) Projection

Compound Name	RT (min)
Perfluorobutanoic acid	2.13
Perfluoropentanoic acid	2.81
4:2 Fluorotelomer sulfonic acid	3.21
Perfluorohexanoic acid	3.40
Perfluorobutanesulfonic acid	3.50
Perfluoroheptanoic acid	4.00
6:2 Fluorotelomer sulfonic acid	4.34
Perfluorooctanoic acid	4.60
Perfluorohexanesulfonic acid	4.90
Perfluorononanoic acid	5.20
Perfluoroheptanesulfonic acid	5.54
8:2 Fluorotelomer sulfonic acid	5.40
Perfluorodecanoic acid	5.72
Perfluorooctanesulfonic acid	6.10
Perfluoroundecanoic acid	6.20
Perfluorododecanoic acid	6.70
Perfluorotridecanoic acid	7.20
Perfluorotetradecanoic acid	7.70
Perfluorooctanesulfonamide	7.37
2-(N-Methylperfluorooctanesulfonamido)acetic acid	7.38
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	7.81
4,8-Dioxo-3H-perfluorononanoic acid	4.56
(2E)-3,4,4,5,5,6,6,7,7,8,8,8-Dodecafluoro-2-octenoic acid	4.44
2-Perfluorooctyl ethanoic acid	5.66
2-Perfluorodecyl ethanoic acid	6.83
2H-Perfluoro-2-octenoic acid	4.40
2H-Perfluoro-2-decenoic acid	5.64
2H-Perfluoro-2-dodecenoic acid	6.79
Ethyl heptafluorobutyrate	3.20
2H,2H,3H,3H-Perfluorooctanoic acid	4.41
3-Perfluoroheptylpropanoic acid	5.68
Perfluoropentanesulfonic acid	3.83
Perfluorononanesulfonic acid	6.07
Perfluorodecanesulfonic acid	6.60
Perfluorododecanesulfonic acid	7.52
Perfluoro(2-((6-chlorohexyl)oxy)ethanesulfonic acid)	5.84
2-((8-Chloro-1,1,2,2,3,3,4,4,5,5,6,6,7,7,8,8-hexadecafluorooctyl)oxy)-1,1,2,2-tetrafluoroethanesulfonic acid	6.92
10:2 Fluorotelomer sulfonic acid	7.11
N-Methylperfluorooctanesulfonamide	8.15
8:2 Fluorotelomer sulfonamide	8.52
N-((1,1,2,2,3,3,4,4,5,5,6,6,7,7,8,8,8-Heptafluorooctyl)sulfonyl)glycine	6.14
N-Methyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide	8.17
N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide	8.51
Phosphonic acid, (tridecafluorohexyl)-	3.22
PFOPA	4.35
Perfluorodecylphosphonic acid	5.57
6:2 Fluorotelomer phosphate diester	7.92
6:2/8:2 Fluorotelomer phosphate diester	8.48
8:2 Fluorotelomer phosphate diester	8.90
Bis(tridecafluorohexyl)phosphonic acid	7.39
(Heptafluorooctyl)(tridecafluorohexyl)phosphonic acid	8.07
Bis(heptafluorooctyl)phosphonic acid	8.60
Ammonium bis(N-ethyl-2-perfluorooctylsulfonaminoethyl)phosphate	9.35



21 measured RT's
32 projected RT's

Analytical and Bioanalytical Chemistry (2019) 411:3507–3520
<https://doi.org/10.1007/s00216-019-01829-8>

RESEARCH PAPER

A single analytical method for the determination of 53 legacy and emerging per- and polyfluoroalkyl substances (PFAS) in aqueous matrices




Timothy L. Coggan¹ • Tarun Anumol² • James Pyke² • Jeff Shimeta¹ • Bradley O. Clarke¹

Received: 30 November 2018 / Revised: 27 February 2019 / Accepted: 3 April 2019 / Published online: 9 May 2019
© Springer-Verlag GmbH Germany, part of Springer Nature 2019




Monitoring suspect PFAS

Adding focus through Retention Time (RT) **Prediction**



Monitoring suspect PFAS

Adding focus through Retention Time (RT) **Prediction**



Mansouri et al. J Cheminform (2018) 10:10
<https://doi.org/10.1186/s13321-018-0263-1>

Journal of Cheminformatics

RESEARCH ARTICLE

Open Access




OPERA models for predicting physicochemical properties and environmental fate endpoints

Kamel Mansouri^{1,2,3*}, Chris M. Grulke¹, Richard S. Judson¹ and Antony J. Williams¹

29 physiochemical properties in addition to the count of -CF₂- used to create prediction model

Monitoring suspect PFAS

Simultaneous Quantitation and Screening



Filter compounds that are Verified, Needs Review, Not Detected

Filter Targets and Suspects

Set outliers flags according to SANTE guidelines

MassHunter Quantitative 'Batch-at-a-glance' view


New Screener Tool

Review Isotope pattern match

Review Fragments

Monitoring suspect PFAS

Simultaneous Quantitation and Screening



Monitoring suspect PFAS

Screening summary PDF report

Screening Summary Report




Sample name:		16AGW06	Good		4	Warning		17	Error	354
Status	Screening Summary Report	Formula	R.T.	R.T. Diff.	Match Score	Target Ion	Mass Accuracy	# of Qualified Ions	Final Conc.	
!	(Heptafluoropropyl)trimethylsilane	C6H9F7Si	2.694	2.692		241.0289	3.95 PPM	2		
+	PFBA	C4 H F7 O2	2.079	0.041		212.9792	0.63 PPM	2	472.7851	
+	PFPeA	C5 H F9 O2	2.777	0.031		262.9760	0.47 PPM	2	448.9793	
!	4:2 FTS	C6 H5 F9 O3 S	3.162	0.048		326.9743	-0.22 PPM	1	880.0983	
+	PFHxA	C6 H F11 O2	3.363	0.036		312.9728	-0.28 PPM	2	475.3056	
!	PFBS	C4 H F9 O3 S	3.469	0.031		298.9430	-0.39 PPM	1	359.3393	
!	3H-Perfluorobutanoic acid	C4H2F6O2	3.530	0.499		194.9886	-1.07 PPM	1		
!	Perfluorooctanesulfonate	C8HF17O3S	5.933	1.754		498.9302	-1.26 PPM	2		
+	6:2 FTS	C8 H5 F13 O3 S	4.266	0.076		426.9679	-0.59 PPM	2	911.3406	
!	2H-Perfluoro(2-methylpentane)	C6HF13	3.956	0.505		318.9798	-0.97 PPM	2		
!	Perfluoro(2-ethoxyethane)sulfonic acid	C4HF9O4S	3.785	0.778		314.9379	0.46 PPM	2		
!	Perfluoropentanesulfonic acid	C5HF11O3S	4.165	0.729		348.9398	-0.65 PPM	2		
!	1-Hydroperfluoroheptane	C7HF15	4.511	0.662		368.9766	-0.44 PPM	2		
!	PFNA	C9 H F17 O2	5.058	0.143		462.9632	-0.38 PPM	2	303.9080	
!	2,3,3,3-Tetrafluoro-2-(perfluoropentoxy)propan-1-ol	C8H3F15O2	4.526	0.718		414.9821	0.81 PPM	2		
!	1H-Perfluorohexane	C6HF13	3.956	1.326		318.9798	-0.97 PPM	2		
!	((Perfluorooctyl)ethyl)phosphonic acid	C10H6F17O3P	5.300	0.485		526.9710	4.37 PPM	1		
!	4-[3-(Perfluorobutyl)-1-propyloxy]benzyl alcohol	C14H13F9O2	6.167	0.221		383.0699	2.89 PPM	1		
!	(Perfluorooctyl)propanoyl chloride	C11H4ClF17O	5.927	0.137		508.9606	-2.22 PPM	1		
!	PFOS	C8 H F17 O3 S	5.933	0.167		498.9302	-1.23 PPM	2	63.0760	
!	FOSA	C8 H F17 O3 S	7.271	0.015				1	0.4484	


Flagging RT outlier

Flagging number of verified ions


Targeted Suspect Screening





Custom PFAS database with >150 compounds
MS/MS spectra and retention time data available for a subset of compounds





Targeted Suspect Screening Results



WWTP1 Final Effluent






Standard




Targeted Suspect Screening Results

Additional PFAS identified with LC-Q/TOF



Targeted Suspect Screening Results




- **Targeted:** Only captured ~6% of compounds
- **Suspects:** 10-16% of fluorinated compounds
- **Unknowns:** 66-81% of fluorinated compounds
- Large portion of fluorinated anionic compounds in sample unknowns

60-70% of fluorine containing molecules in these samples is 'unknown'

Complete Workflow


Known knowns
Expanded targeted
list (~50 compounds)



Continually refining
targeted method



Known unknowns
Expanded PFAS database



Unknown unknowns
Adding identified
compounds to database list

Measuring PFAS by LC-IM-Q/TOF

The ultimate tool for PFAS identification


“Rapid Assessment of Isomeric Diversity in PFAS by Ion Mobility Spectrometry-Mass Spectrometry (IMS-MS)”

James Dodds & Erin Baker , NC State




Agilent 6560 IM-Q/TOF

Ion path design



Agilent 6560 IM-Q/TOF


Basic operational principle of Ion Mobility for conventional DC uniform field IMS



$$v = K E \propto \frac{e E}{P \sqrt{T} \Omega}$$

Measuring PFAS by LC-IM-Q/TOF

An additional dimension of separation




Combining the Power of
Accurate Mass and Structural Characterization
With Current UHPLC (and GC) Separations

Measuring PFAS by LC-IM-Q/TOF


Mass Spectrometry data

PFOS Branched/Linear Isomers


(I) Linear (*n*) PFOS
~ECF,TLM




(II) Branched (*m*) PFOS
~Only ECF







(III) Impure PFOS



(IV) Linear (*n*) - PFOS




PFOS MS/MS [M-H]⁻



Measuring PFAS by LC-IM-Q/TOF

PFOS Conformations resolved by IM



Summary

- TQ→QTOF→IM/QTOF resolving PFAS reveals complexity of the problem
- Simultaneous target quantitation and suspect screening workflow
- Need for direct injection methods for broad screening
- Adding confidence to putative identifications using RT projection and prediction



Acknowledgements

Timothy Coggan and Bradley Clarke
RMIT University, VIC Australia

James Dodds and Erin Baker
North Carolina State University, NC USA