

Dr. Andrew Eaton
MWH Labs
Monrovia, California

CHROMIUM⁺⁶

Deja Vu All Over Again



NEMC
August 15, 2011

Presentation Outline

1. History – Science and Society...
2. Chemistry of Chromium Sources and Forms
3. Analytical Issues
4. National Occurrence
5. Treatment process occurrence “anomalies”
6. Conclusions

Cr⁶: A 10 year+ Regulatory Timeline





Chart concept courtesy of Nicole Blute (Arcadis)

Cr⁶⁺: A 15 year+ Analytical Timeline



Chromium Exists as Several Chemical Species

Most common oxidation states: 0, +3, +6

0: Elemental Chromium (Cr)

+3: Trivalent Chromium

Species: Cr^{+3} , Cr_2O_3


+6: Hexavalent Chromium

Species: CrO_4^{2-} , $Cr_2O_7^-$

Groundwaters – largely Cr^{+6}

Surface waters – more Cr^{+3}

Seidel (2004) AWWARF Study



Forms of Chromium and Toxicity

- Trivalent (Cr^{+3})

- Essential Nutrient



- Converts to Hexavalent in the presence of oxidants

- Hexavalent (Cr^{+6})

- Highly water soluble

- Natural AND industrial sources

- Carcinogenic by inhalation for sure; ingestion??


- OSHA PEL = 5 $\mu\text{g}/\text{M}^3$

- Converts to Trivalent in the stomach

- Converts to Trivalent in the presence of Ferrous Iron



Chromium Speciation in Water (eH-pH Diagram for Chromium)




In most natural waters Cr⁺⁶ is the most soluble and stable form!

As will be shown, analytical data on occurrence supports this... thermodynamics works!

Analytical Issues for the Measurement of Hexavalent Chromium

- Ion Chromatography Methods for Cr^{+6} are **EXTREMELY** sensitive
 - Quantitative data at **0.020** ug/L



- Cr^{+6} is very stable in most waters

- The holding time WAS 24 hours; now it is proposed as **2 weeks (EPA DRAFT method 218.7)**

ne	Peak Name	Height AU	Area AU*min	Rel.Area %	Amount
	CRVI-LOW	0.000	0.0000200	95.08	0.018

- For Cr^{+6} lab blanks are not a big issue
- The challenge is to prevent oxidation of Cr^{+3}

Preservation Issues Depend on the eH-pH Diagram and Preventing Redox Reactions


- Bring the pH up
 - 218.6 said >9 (or 9.2-9.7)
 - **218.7 says >8**
 - 10 years worth of data in our lab supports the fact that you don't need to have the pH above 9...
- Minimize the ability of free chlorine to oxidize Cr^{+3}
 - Add NH_4SO_4 to form chloramines
 - Work from EPA region 6 and EPA-Cincinnati demonstrates that this is effective
 - Without adding the ammonia, oxidation will occur
 - Also no impact from Fe (III) reduction

Buffering Options Have Varied Over the Years

- 218.6 (original) used strong $\text{NH}_4\text{SO}_4/\text{NH}_4\text{OH}$ buffer
 - Impacted column capacity
 - Caused some signal suppression
- CDPH proposed a borate buffer
 - Did not address chlorine issues because CDPH measurements were source waters
- Newer methods (e.g. 218.7) use weaker buffers
 - Dilute $\text{NH}_4\text{SO}_4/\text{NH}_4\text{OH}$ (**liquid**)
 - Sodium carbonate/bicarbonate + NH_4SO_4 (**solid**)

So: Analytical Issues for Cr⁺⁶- Some Real, Some “Imagined”



- Holding Time?
 - 24 hours (original EPA 218.6), CAUCMR
 - 5 days (EPA DW Guidance)
 - 28 days (40CFR136)
 - 14 days (EPA 218.7)
- pH adjustment?
 - Greater than 9 (drinking water)
 - 9.2 to 9.7 (wastewater)
 - >8 (EPA 218.7)
- Field Filtration in 218.6 (wastewater vs drinking water)
 - Cross Contamination. Is it even necessary?



In Draft Method 218.7 there Are Still Analytical “Issues” with Cost Implications

- How frequently do you need to prepare the preservative?
- Can you send liquid preservative to the field? (e.g, pre-preserve samples?)
- Do you need to filter samples before injection on the IC?
- Do you need to chill samples during transport?

Comparing LCMRLs for Cr⁶ and Total Chromium. Cr⁶ method is 5-10X more sensitive




**MWH Total Chromium LCMRL =
120 ng/L
(DI water with digestion)**

**EPA 218.7 Chromium-6 LCMRL =
22 ng/L**

Total Chromium Measurements Are More Problematic

- Chromium has 2 major stable isotopes at mass 52 (83.8%) and mass 53 (9.5%)
- Traditional ICPMS is prone to interferences from ArC (mass 52) and ClO (mass 53)
- For UCMR3 EPA recommends digestion (regardless of turbidity) to minimize the ArC interference
- Collision cell technology is not yet permitted for drinking waters

Does Digestion Solve the Carbon Interference Problem?




Occurrence Data Impact Will Depend to a Large Extent on Toxicity Decisions


- Some existing public data gives a sense of what's present, but not the whole picture.
- What if you start looking at lower levels?
- What about Cr⁺⁶ vs Cr⁺³ occurrence?
- What happens in treatment plants?

Chromium⁺⁶ Occurrence

- EPA National Occurrence (6-year review data)
 - 18% of sources >10 ppb (total chromium-not ⁺⁶ only)
- California (CA-UCMR data for Cr⁺⁶)
 - 11% of sources >10 ppb (hexavalent)
 - 33% of sources >1 ppb
- MWH (Cr⁺⁶) (1500 non California samples, representing 2011 data)
 - 1% >10 ppb
 - 20% >1 ppb
 - 50% >0.1 ppb
 - 70% >0.05 ppb




Cr⁶⁺ Occurrence – the Lower We Look the More We Will Find....




How Much Does Cr⁺⁶ Vary Over Time?

Answer: Not Very Much



Paired Distribution System Samples (Cr-T and Cr⁶⁺) From GW or SW Show Trends



GW is Higher Than SW


GW is more likely to be predominantly Cr⁶⁺

Most of these samples are <1 ug/L Total Chromium (semi-quantitative range for Total Chromium)

Potential for Chromium Cross Contamination (or analytical errors)?

- Stainless Steel (10%+ Chromium)
 - Could leach Cr^{+3} under acidic conditions
- Buffer Chemicals used for pH adjustment
 - NaOH has traces of Cr^{+6}
 - KCO_3 has traces of Cr^{+6}
- ICP/MS Carbon Correction (Eaton, 2001)
 - ArC in ICPMS has same mass as Cr (52)
- Treatment process adding Cr^{+6}
 - Source water 0.15 ppb
 - Treated after coagulant – 1.5 ppb

A Few Interesting Tidbits of Cr⁺⁶ Occurrence Through a Conventional Treatment System



Note that this is UNUSUAL – we **have not seen** the increase in all plants we have studied.

Total Chromium results mimic the Cr⁺⁶, so this is likely Cr being added from treatment chemicals

What About RO Treatment?

- In general RO is very effective at removing chromium.....
- But again, we have on occasion seen cases where Cr^{+6} is detected in RO permeate at levels in excess of 0.2 ug/L.
- This seems more likely to be leaching from pipes than actual RO breakthrough, but the jury is still out.

Conclusions- Analytical Methods Assessment

- Hexavalent chromium by ion chromatography (e.g. 218.6 or 218.7) is a rugged and sensitive analytical method.
- Total chromium by ICPMS using drinking water approved methods is rugged down to 0.2 ppb, IF samples are digested. Otherwise you may get positive bias (up to low ppb levels).
- Collision cell technology should also address the bias issues, but is not approved for DW yet.

Conclusions – Precautions and Occurrence

- The analytical issues relate more to sample collection than analysis per se.
- There is ample evidence that Cr^{+6} is widespread at low levels. Health effects and treatment costs will trump occurrence as drivers.

Any Questions?

Contact Information

Dr. Andrew Eaton
Technical Director
MWH Laboratories
750 Royal Oaks Drive
Monrovia, CA 91016

Andrew.d.eaton@us.mwhglobal.com