

Metal speciation by IC-ICP-MS: Arsenic and Chromium case studies

Tatyana S. Pinyayev¹, Robert A. Wilson², Karen Herbin-Davis³,
Michael J. Kohan³, John T. Creed¹ and David J. Thomas³

¹ *United States Environmental Protection Agency, National Exposure Research Laboratory,
Microbiological and Chemical Exposure Assessment Research Division, Cincinnati, OH, 45268*

² *Student Services Contractor, NERL, MCEARD, Cincinnati, OH, 45268*

³ *United States Environmental Protection Agency, National Health and Environmental Effects Research
Laboratory, Experimental Toxicology Division,
Research Triangle Park, NC, 27711*

*Mention of trade names or commercial products does not
constitute endorsement or recommendation for use*

Arsenic


Chromium

As^V / As^{III} Type 1 carcinogen	Cr^{VI} / Cr^{III} Type 1 carcinogen AND essential nutrient
World Health Organization:	
10 ppb <u>guideline</u> for drinking water	50 ppb <u>guideline</u> for total chromium in drinking water
US EPA:	
<u>enforceable</u> maximum contaminant level (MCL) of 10 ppb - ultimate goal (MCLG) – 0 ppb	<u>enforceable</u> maximum contaminant level (MCL) of 100 ppb -Office of Environmental Health Hazard Assessment recently proposed 20 ppt public health goal (PHG)
Research objective:	
<i>To study metabolism of arsenic oxides by intestinal microbiota and to understand the exposure and toxicity implications of these transformations</i>	<i>To develop a reliable analytical method for determination of Cr^{III} and Cr^{VI} in drinking water with low ppt detection limit</i>
Analytical challenges:	
<ul style="list-style-type: none"> - Speciate multiple (up to 12) unknown metabolites from a complex sample matrix in one chromatographic run - Synthesize and characterize STDs (no SRMs) 	<ul style="list-style-type: none"> - Assure preservation of Cr^{VI} during sample handling - Designing a separation system that has good long term stability and minimizes polyatomic interference due to mobile phase and matrix constituents: ⁴⁰Ar¹²C, ³⁴S¹⁸O, ¹⁶O³⁷Cl


Part 1: Arsenic analysis by IC-ICP-MS

- *Multiple metabolites are isolated from a complex matrix (2-phase)*
- *Ion-chromatographic separation with elemental detection by ICP-MS*
- *Distribution of metabolites within sample phases*


Arsenic species (metabolic pathway)



Pre-systemic metabolism of MMA in GI tract by anaerobic bacteria



- Cecum contains H_2S -producing bacteria
- Potential thiolation site




As, ppb	Hours at 37°C					
	0	3	6	12	24	48
0	3	—	—	—	3	3
20	3	—	—	—	3	3
200	3	3	3	3	3	3
1000	3	3	3	3	3	3
2000	6	3	3	3	3	6


Snap-freeze @ -70°C → ship

Analyze supernatant by IC-ICP-MS


IC-ICP-MS mass chromatogram of MMA metabolites in **supernatant**



Preparation of cecum samples for IC-ICP-MS analysis







IC-ICP-MS mass chromatogram of MMA metabolites in **extract of cecum solids**




= unique to extract
from cecum solids

Mass-balance distribution of MMA metabolites in supernatant and cecum solids (at 48 hours)



Part 1: conclusions



Part 2: Chromium analysis by IC-ICP-MS

- *Analytical capabilities of ICP-MS for Cr speciation studies*
 - *Polyatomic signal interference (mobile phase, matrix)*
 - *Separation/Sensitivity*
- *Preservation of native Cr^{III} and Cr^{VI} during shipping and analysis*

Potential polyatomic interferences associated with detection of Cr by ICP-MS at m/z 52 and m/z 53


- Chromium isotopes: 50 (4.3%), **52** (83.8%), **53** (9.5%) and 54 (2.4%)
- Any diatomic species at these masses can contribute to background signal and chromatographic peaks

Origin of the diatomic background interference:	Possible interferences at m/z 52:	Possible interferences at m/z 53:	Polyatomic induced chromatographic peaks from matrix
plasma gas	$^{36}\text{Ar}^{16}\text{O}$	$^{36}\text{Ar}^{17}\text{O}$	
	$^{38}\text{Ar}^{14}\text{N}$	$^{38}\text{Ar}^{15}\text{N}$	
mobile phase (sulfates)	$^{36}\text{S}^{16}\text{O}$; $^{34}\text{S}^{18}\text{O}$	$^{36}\text{S}^{17}\text{O}$	Sulfate (SO_4^{2-})
mobile phase (carboxylates)	$^{40}\text{Ar}^{12}\text{C}$	$^{40}\text{Ar}^{13}\text{C}$	Carbonate (CO_3^{2-})
		$^{37}\text{Cl}^{16}\text{O}$	Chloride (Cl^-)

Effect of mobile phase composition on polyatomic background signal

Considered mobile phases:

1. Ammonium oxalate (AO)
 2. Ammonium sulfate (SO₄)
 3. Ammonium nitrate (NO₃)
- 1) Oxalate contributes to high background signal; also reacts with chromium (complexation)
 - 2) Sulfate also reacts with chromium to produce unidentified chromium compound (elutes under 4 minutes); provides less background at m/z 52 compared to oxalate
 - 3) Ammonium nitrate allows separation of Cr^{III} and Cr^V under 14 minutes and produces least background signal




*All chromium compounds were spiked at 1 ppb level


* Naturally occurring matrix ions contributing to diatomics such as $^{40}\text{Ar}^{12}\text{C}$, $^{34}\text{S}^{16}\text{O}$, $^{16}\text{O}^{37}\text{Cl}$ are base-line resolved at these conditions (not shown)

** Dionex AS7 ion-exchange column, 1mL/min flow rate, ambient temperature

Effect of collision gas flow rate on Cr^{VI} detection limit (3σ)



Sample preservation issue: redox interconversion of chromium species




1. Cr^{VI} is only predominant in oxidative waters
2. Disinfectants added to portable water are oxidizing in nature, favoring formation of Cr^{VI} from Cr^{III}
3. Analytical challenge is to develop a method that preserves native Cr^{III} and Cr^{VI} from the time of sample collection until detection

**Pourbaix diagram for chromium adapted from Kotaś, J.; Stasicka, Z
Environmental Pollution (2000) 107 (3): 263–283*

Stability of Cr^{III} in drinking water as a function of time


1. Water samples were collected in different states and transported to lab
2. Cr^{III} standard spike was added to waters in lab and held at room temperature
3. Aliquots of each water were taken at time intervals, heated with EDTA and a formed Cr^{III}-EDTA complex was detected
4. Cr^{III} recoveries were 0% for all waters (except Tennessee, 50%) after 8 days



** 1 ppb standard spikes, 0.8 mL/min He*

Stability of Cr^{III} in EDTA-spiked drinking water at room temperature

1. Water samples were collected in different states and transported to lab
2. EDTA followed by Cr^{III} standard spike was added to waters, held at room temperature
3. Aliquots of each water were taken at time intervals, heated for 1 hour and a formed Cr^{III}-EDTA complex was detected



** 6 ppb Cr^{III} standard spike, 0.8 mL/min He*

Conclusions

- Ion-exchange chromatography is suitable for separation of Cr^{III} and Cr^{VI} in drinking water
- Preservation of Cr^{III} spike in water samples is challenging unless immediately heated after addition of EDTA

Acknowledgement



Jack
Creed



David
Thomas



Carol
Brockhoff



Madhavi
Mantha




Patty
Creed



Tony
Wilson



Questions?



"I think you should be more explicit here in step two."