

Introduction

EPA Method 8270¹ is commonly used for the analysis of a wide variety of semivolatile organic compounds (SVOCs) in environmental matrices. Due to the large number of compounds the analytical run can be as long as 45 minutes. Using a GC tandem mass spectrometer² (GC/MS/MS or GC-Triple Quad) with electron ionization (EI) and dynamic multiple reaction monitoring (dMRM) a 10-minute analysis time was developed to increase efficiency in the laboratory. The analyte list has six internal standards, eight surrogates and 72 target compounds including polycyclic aromatic hydrocarbons³ (PAHs) and other acid, base, neutral compounds. Results of the fast 8270 method using GC/MS/MS were compared to traditional GC/MS on real extracts from environmental samples prepared by liquid/liquid extraction⁴ with dichloromethane (DCM) at pH ≤ 2.

Equipment

Instrument: Agilent 7890/7010A Triple Quadrupole GC/MS

Software: MassHunter B.09

Liner: Agilent 5190-2295 Ultra Inert, wool

Column: Agilent DB5MS UI 20 m x 0.18 mm x 0.18 μm

Sample Preparation: Corning AOS™ liquid-liquid extractors

Reagents and Standards

Calibration standards: Supelco

DCM: Foxpure D454-4

Alternate source: Restek

Reagent Water: Milli-Q

Internal standards: Phenova

GC/MS/MS Operating Conditions

GC Conditions

Inlet temperature: 280°C

Oven temperature: 40°C

Ramp: 35°C/min to 320°C

Hold for 2 min

Injection volume: 1 μL

Split ratio: 100:1

Carrier gas: Helium

Flow: 1.0 mL/min

MS/MS Conditions

Ion source: EI

Ionization voltage: 70 eV

Interface temperature: 280°C

Ion source temperature: 325°C

dMRMs: 176

(Examples below)

Triple Quadrupole Acquisition Method - MS Parameters Report

Compound Table:


Compound Name	ISTD ?	Precursor ion	MS1 res	Product ion	MS2 res	RT	Left Delta RT	Right Delta RT	Dwell (ms)	CE (V)
2-Chlorophenol-d4		132	Wide	95.9	Wide	2.40	0.50	0.50	19.7	10
2-Chlorophenol-d4		132	Wide	67.8	Wide	2.40	0.50	0.50	19.7	10
Bis(2-chloroethyl)ether		95.1	Wide	65	Wide	2.43	0.10	0.10	10	5
Bis(2-chloroethyl)ether		93.1	Wide	63	Wide	2.43	0.10	0.10	10	0
2-Chlorophenol		128	Wide	64	Wide	2.46	0.10	0.15	9.2	15
2-Chlorophenol		128	Wide	63	Wide	2.46	0.10	0.15	9.2	30
1,3-Dichlorobenzene		146	Wide	111	Wide	2.55	0.20	0.42	8.2	15
1,3-Dichlorobenzene		146	Wide	75	Wide	2.55	0.20	0.42	8.2	30
1,4-Dichlorobenzene		146	Wide	111	Wide	2.60	0.42	0.42	10.3	15
1,4-Dichlorobenzene		146	Wide	75	Wide	2.60	0.42	0.42	10.3	30
1,4-Dichlorobenzene-d4	X	150	Wide	115	Wide	2.60	0.10	0.10	9.2	15
1,4-Dichlorobenzene-d4	X	150	Wide	78	Wide	2.60	0.10	0.10	9.2	30

Calibration



Calibration points: 9

Calibration range: 1-100 ug/ml in DCM

Curve fits: Average of Response factors, %RSD≤20 or Quadratic, r²≥0.99



GC/MS/MS versus GC/MS - Creosote sample



Initial Demonstration of Proficiency

Initial Demonstration of Proficiency, examples. Full list available upon request.

Compound, Spiked 50 ug/L	Quant dMRM	RT	IDP-1	IDP-2	IDP-3	IDP-4	AVG	STDEV	%RSD	AVG %REC
Phenol	94.0 -> 66.1	2.4	53	55	55	54	1	2	109	
1,2-Dichlorobenzene-d4	152.0 -> 78.0	2.7	41	38	41	35	39	3	8	78
2-Methylphenol	107.0 -> 77.0	2.7	53	55	56	54	55	1	2	109
Naphthalene	128.1 -> 102.1	3.4	48	48	48	48	48	0	1	96
Hexachlorocyclopentadiene	237.0 -> 119.0	4.0	47	48	49	48	48	1	2	95
2-Nitroaniline	138.0 -> 65.0	4.3	54	53	55	53	54	1	2	107
Acenaphthene	153.1 -> 77.0	4.6	49	49	50	49	49	0	0	99
4-Nitrophenol	138.9 -> 109.0	4.7	54	55	56	56	55	1	2	111
2,4-Dinitrotoluene	165.0 -> 63.0	4.7	54	54	55	54	54	0	1	108
Pentachlorophenol	265.9 -> 167.0	5.5	51	52	53	53	52	1	2	105
Phenanthrene	178.1 -> 152.1	5.6	51	51	51	51	51	0	0	102
Benzo (b) fluoranthene	252.1 -> 250.1	8.2	54	55	56	57	56	1	2	111
Benzo (a) pyrene	252.1 -> 250.1	8.4	57	57	57	57	57	0	1	114
Benzo (g,h,i) perylene	276.1 -> 274.1	9.5	56	57	57	58	57	1	1	114

Conclusions

The Fast 8270 method using GC/MS/MS has several advantages:

- Shortened run time increases efficiency in the laboratory.
- Tandem mass spectrometry reduces matrix interferences.
- Dynamic MRM eliminates time segments.

Concerns include reproducibility over time and lower limit of quantitation (LLOQ) due to the small amount of extract on column.

Acknowledgements

Special thanks to the following people for their support:

Diane Gregg Reese, EPA Senior Science Advisor, 8270E guidance

David Hiser, EPA Chemist, data review/LIMS assistance

Rachel Rathbone, EPA Chemist, sample preparation

Lisa Wool, EPA Environmental Scientist, procurement assistance

Tanya Slaughter, EPA Librarian, literature searching

Matthew Reynolds, EPA Facility Manager, printing

Bruce Quimby, Agilent Technologies, GC advice

Anastasia Andrianova, Agilent Technologies, MS/MS advice

Ron Tackett, Agilent Technologies, MassHunter support

References

¹U.S. EPA. Method 8270E Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS). SW-846 VI (2018)

²Yost & Enke. *J. Am. Chem. Soc.*, 100(7), (1978) 2274-2275

³Anderson et al. *J Chrom A*, 1419 (2015) 89-98

⁴U.S. EPA. Method 3520C Continuous Liquid-Liquid Extraction, Rev. 3; SW-846, (1996)

Disclaimer

The views expressed in this poster are those of the author and do not necessarily represent the views or the policies of the U.S. Environmental Protection Agency. Any mention of trade names, manufacturers, or products does not imply an endorsement by the United States Government or the U.S. EPA. The U.S. EPA and its employees do not endorse any commercial products, services, or enterprises.