





Environmental Incremental Sampling Methodology Update

Mark L. Bruce Ph.D.

August 8, 2019

Copyright © 2019, Eurofins TestAmerica. All rights reserved.

www.TestAmericalnc.com

ISM Purpose



Produce representative field samples and analytical subsamples

 Parameters / analytes in the right proportions to represent site conditions within the decision unit



eurofins Environment Testing TestAmerica

Environment Testing TestAmerica

🛟 eurofins

ISM Field Process (most common)

Systematic random increment collection

Collect surface plugs











ISM Lab Process (most common)



Air dry, disaggregate, sieve, mill, 2D slabcake subsampling







eurofins Environment Testing TestAmerica

New Understanding of Site Conditions



Presumed small scale heterogeneity was minimal



Except Pb on shooting range and other obvious heterogeneous sites

Beginning to realize that small scale heterogeneity is common



eurofins Environment Testing TestAmerica



What if we could "see" the contaminant distribution? Do we need all the details or just a good average for a decision unit?

ISM is designed to manage the small and moderate scale heterogeneity and produce the good average.





PCB infused tar balls

- How many "nuggets" in field sample?
- How many "nuggets" in lab subsample?



Collocated PCB Samples



Hawaii



Environment Testing TestAmerica **New insight – minimum # increments**



Original default – 30 increments

- Except for known highly heterogeneous sites
 - 50 100 increments

New default – 50 increments for unknown distributions

- Hawaii Technical Guidance Manual 2017
- Can use 30 increments for known low heterogeneity sites

Tools and techniques for increments



When push coring tools don't push well

- Cordless drill, Speedbor bit and collection plate
 - Not for dry sand





Tools and techniques for increments



When simple push coring tools don't push well

Slide hammers





Environment Testing TestAmerica **Tools and techniques for increments**



When push coring tools don't push well

Electric hammer and chisel bit





ITRC ISM update - selecting lab process options

- Air drying yes / no
- Disaggregation
 - Tool selection
- Sieving yes / no,
 - If yes, size selection
- Milling yes / no,
 - If yes, tool selection
- Subsampling
 - Tool selection



Concerns about analyte loss during air drying

- 2003 matrix spike study showed >80% loss of weakly absorbed low boiling SVOCs
- Led to 2D slabcake subsampling on as-received soil samples
 - Fast
 - Low cost
 - Reduces analyte loss
 - Better than "stir and scoop" subsampling
 - Not as reproducible as dry, disaggregate, sieve, 2D slabcake
- Triplicate comparison for PAH samples



Faster, cleaner air drying - Air drying towers







Environment Testing TestAmerica



Disaggregation and sieving





Old reliable mills

Puck mill & ball mill







Ball Milling PCB samples (7 replicates)

Aroclor 1242







Mini ball mill for Wet Chem subsamples











Disposable scoops for 2D slabcake subsampling



Environment Testing TestAmerica





Large scale (10 g) metals digestion / 250 mL digestion cup



ENVIRONMENTAL EXPRESS°

a Cole-Parmer company

🛟 eurofins **Environment Testing TestAmerica**



Process = air dry, disaggregate, sieve (#10), 2D slabcake subsample



Environment Testing TestAmerica



Surface soil that was exposed to air in the field



eurofins Environment Testing



No loss from surface soils during air drying



eurofins Environment Testing TestAmerica

Lab Process Quality Assurance



ITRC ISM update

- Precision
 - Lab replicates
- Bias
 - Surrogate spiking
 - Solid spike material
 - Contamination control
 - Blank materials





Field techniques & tools for sampling at more difficult sites

- More increments
- More sampling tools
- Lab tools for higher quality and more productive processing
 - More selection guidance
 - Tools update
 - Quality Assurance guidance

Acknowledgements



Roger Brewer - Hawaii Department of Health



Joe Boyd - Environmental Express

ENVIRONMENTAL EXPRESS[®] a Cole-Parmer company

Ed Corl - Navy



Jay Clausen - USACE







Mark Bruce

- Corporate Technical Director
- Eurofins TestAmerica
- **330-966-7267**
- Mark.Bruce@TestAmericalnc.com
- www.EurofinsUS.com | www.TestAmericaInc.com